ST.330.008

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

INDICATIONS FOR USE

This product is a surgical device in the form of temporary screw used as an aid in orthodontic movement procedures.

DESCRIPTION

The Self-Drilling Anchorage Neodent Implant is made of titanium alloy, as per standard ASTM F136.

MACROGEOMETRY: • Self-drilling conical geometry; • Head with two different transmucosal heights or belt (low = 0 mm and medium = 1 mm); • Hex for installation wrenches and channel with central hole for running wire sutures.

Diameter of Implant (mm)	Belt (in mm)	Lengths (mm)
1.3	Low	7, 9, and 11
	Medium	5, 7, 9, and 11
1.6	Low	7, 9, and 11
	Medium	5, 7, 9, and 11

APPLICATIONS

It is suitable for installation on mandible or maxilla, as an aid in anchoring orthodontic movement procedures.

CONTRAINDICATIONS

This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: titanium alloy Ti6Al4V-ELI. In the presence of acute inflammatory or infectious processes in live tissue, unsuitable bone volume or quality, serious medical problems such as; bone metabolism disorders, blood clotting disorders, unsuitable healing capacity, insufficient oral hygiene, patient uncooperative and not motivated, abuse of drugs or alcohol, psychosis, prolonged functional disorders which resist any treatment with medications, xerostomia, weakened immunological system, diseases which require the use of steroids, endocrinological diseases, pregnancy.

HANDLING


Select the length and diameter of the implant according to the bone and/or intracanal area available for its installation. The choice of the belt of the implant must be done in accordance with the thickness of the soft tissue. The surgical instrumentation is recommended according to bone density (bone classification of Lekholm and Zarb, 1985) and the selected implant, and must be performed according to the following table:

	Self-Drilling Anchorage Implant 1.3	Self-Drilling Anchorage Implant 1.6
Bone density I and II	Drilling with the 1.1 mm Bit.	Drilling with the 1.3 mm Bit.

	Self-Drilling Anchorage Implant 1.3	Self-Drilling Anchorage Implant 1.6
Bone density III and IV	It is necessary to use the punch key before installation of the implant. It is not necessary to use drills.	

When the use of drills is indicated, do the drilling of the bone bed using bits in good cutting condition with rotation between 80 and 100 rpm, under plenty of irrigation. The depth of the drilling must be in accordance with the planning of the final position of the implant. Laser marking on the bit represent the implant depths. Neodent Bits should be used according to their own Instructions for Use. For more information, consult the Neodent catalog.

INSTALLATION OF THE DENTAL IMPLANT

WARNING: The figures are merely illustrative. They do not represent the actual dimensions and features of the product.

- The cardboard and blister packagings should be opened, manually, without the use of sterile gloves.
- Break the seal of the cardboard packaging and remove the blister.
- Remove the protection from the blister pack.
- Deposit the sterile flask over the surgical field. **NOTE:** The clear tube and implant should be handled with a sterile surgical glove, in a surgical environment.
- Hold the bottle using the non-dominant hand and take the lid off.
- For installation using a surgical motor, hold the implant with the Orthodontic Anchorage connection for contra-angle handpiece, keeping the connection stable and slightly rotating bottle, searching for the perfect fit between the connection and the implant.
- Take the implant to the site of the installation.
- In the surgical motor, use maximum torque 10 N.cm and rotation of 30 rpm.
- For installation with the Ratchet Wrench, use the adapter Ratchet Connections for Contra-angle Handpiece. The insertion may be done manually, with the aid of the Plug Wrench.
- The maximum installation torque suggested is 10 N.cm.
- Ensure that the head of the implant is at supragingival level.

TRACEABILITY LABEL

This product is accompanied by three labels that allow its traceability and should be attached to the following documents: • Medical record; • Collection tax document; • Document to be

delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

PRESENTATION AND STERILIZATION

This product is indicated for single use and is supplied sterilized by gamma radiation method, uniformly packed in packaging that offers triple protection: cardboard, blister-type and clear tube.

PRE-OPERATORY AND PROSTHETIC PLANNING

The diameter, length, position and quantity of implants must be selected for each clinical case, taking into consideration the anatomy, quantity of bone and space available. When necessary execute the diagnostic wax-up of the patient.

PRECAUTIONS

- When the use of drills is indicated, the surgical drilling diameter must always be less than the diameter of the implant planned.
- Installation in keratinized mucosa region is recommended, whenever possible.
- After installation of the Implant for Orthodontic Anchorage, it is recommended to take periapical x-rays to make sure that the positioning of the implant is as planned.
- Note the conditions of the intra-oral tissue, the bone quantity and intraradicular space of the bed receiving the Implant, by means of radiographic and/or tomography examinations. Nonperformance of the pre-surgical assessment may compromise the success of the procedure.
- As for the systemic aspect, consider the general health of the patient. In particular, one must be careful in cases of patients who have allergies to drugs, local or systemic factors that may interfere with the healing process of the bone tissues or soft tissues. For example, the bone already exposed to radiation in the head and neck area, diabetes mellitus, anticoagulation drugs, hemorrhagic diathesis, bruxism, parafunctional habits, anatomically unfavorable bone situation, tobacco abuse, uncontrolled periodontitis, treatable jaws pathologies and oral mucosa abnormalities.
- Treatment with bisphosphonates results in potential risk of peri-implant osteonecrosis.
- Inadequate surgical planning can compromise the performance of the implant and/or orthodontic movement resulting in failures of the system.
- The material to be used during the procedure must be sterile.
- This product is of single use and cannot be re-sterilized.
- Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if the packaging is broken.
- The sterilization is guaranteed only if the blister packaging is not damaged.
- This product must be used immediately after opening the packaging, at the moment of surgery. If it is not used, discard it.

- Do not use the product with the validity expired.
- The maximum installation torque suggested is 10 N.cm. Insertion torque greater than that recommended can make the system inoperative.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Before each procedure check the conditions of the **Neodent** surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the **Neodent** product sequence. The use of instruments of other manufacturers does not ensure the perfect function of the **Neodent** Implant System and exempts any product warranty.
- It is the professional's responsibility to use the **Neodent** products according to the instructions for use.

ADVERSE EFFECTS

The installation of dental implants, as well as any other surgical procedure, may cause a slight discomfort and localized edema. More persistent symptoms can occur such as: chronic pain related to the dental implant, permanent paresthesia, dysesthesia, loss of maxillary/mandibular bone

edge, systemic or localized infection, oroantral or oronasal fistula, adjacent teeth affected unfavorably, irreversible damage to adjacent teeth, fracturing of the implant, jaw, bone or prosthesis, injury of the nerves, exfoliation, hyperplasia.

Any failure in the procedure during treatment may be caused by: inadequate osteotomy, infections, diseases or systemic problems, low quality or volume of the remaining bone, absence or failure of irrigation, use of nonspecific instruments and/or instruments without cutting power, poor oral hygiene, occlusal trauma and lack of specific training.

MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The **Neodent** Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the **Neodent** Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a

professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored, in its original packaging, in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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