ST.330.022

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

Italiano  Análogo

English  Analog

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Español  Analogo

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
	Requisito Brasileiro / Brazilian Requirement / Requisito Brasileño / Requisito Brasileiro / Brasilianische Anforderungen / Exigences brésiliennes

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Não estéril / Non-sterile / No estéril / Non sterile / Nicht steril / Non stérile
	Não reutilizar / Do not reuse / No lo reutilice / Non riutilizzare / Nicht wiederverwendbar / Ne pas réutiliser / Nicht wiederverwendbar / Ne pas réutiliser
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Consérvelo seco / Conservare in luogo asciutto / Trocken lagern / Garder au sec
	Mantener afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conserver à l'abri de la lumière solaire

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This device is intended for a specialized procedure, which should be performed by professionals trained in Prosthesis over implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions.

DESCRIPTION

It is a prosthetic component for laboratory use, made of titanium alloy, in accordance with ASTM standard F136, or of brass. It has a geometry with channels that favor retention on the plaster model.

It is available in the prosthetic interfaces shown in the table below:

Prosthetic interface	Available dimensions (mm)	Material
HE	Ø 3.3, Ø 4.1, Ø4.3 and Ø 5.0	Brass / Titanium
II Plus	Ø 4.3 and Ø 5.0	Brass
Transepithelial Pillar	Single	Titanium
Mini Conical Pillar	Ø4.1 and Ø 5.0	Titanium
Conical Pillar	Ø4.1 and Ø 5.0	Titanium
GT / CM Pillar	Single	Titanium
Universal Trunnion	3.3X4, 3.3X6, 4.5X4 and 4.5X6	Titanium
WS / Zygomatic	Single	Titanium
Micro Pillar	Single	Titanium
CM	Single	Titanium
Facility	Single	Titanium
Smart	Narrow, Regular, and Wide	Titanium

APPLICATIONS

It is indicated for use in the laboratory phase while obtaining the plaster model for making the prosthesis on implant, reproducing the dimensions and positioning of the implant or prosthetic intermediate.

CONTRAINDICATIONS

Wear of any type on this product is contraindicated.

HANDLING

Perform the impression of the implant or prosthetic component by using the corresponding transfer and following appropriate techniques. After obtaining the mold, screw/plug the Analog into the transfer used. When applicable, to tighten the screw, use a maximum torque of 10 N.cm. Pour rubbery material to simulate the soft tissue of the patient. Pour the plaster on the mold, covering the whole Analog. After the plaster setting phase, remove (separate) the mold from the plaster. Unscrew/unplug the transfer from the analog.

PRESENTATION AND STERILIZATION

This product is indicated for single use and supplied non-sterile, packed individually into a package that offers dual protection.

PRECAUTIONS

- Make sure the Analog is compatible with the Transfer selected.
- After making the prosthesis, evaluate its fit on the Analog.
- Inappropriate planning can jeopardize the performance of the implant/prosthesis unit leading to failures of the system, as loss or fracturing of the implant, loosening or fracturing of the prosthetic screws.
- This product is of single use and cannot be reused.
- Reuse of this product may cause:
 - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
 - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if the packaging is broken.
- The sterilization is guaranteed only if the blister packaging is not damaged.
- Do not use the product with the validity expired.
- When making the prosthesis, consult the torque indicated for the prosthetic component to be used. Excessive torque may cause undesirable results, making the system inoperable.

- Before each procedure, make sure the pieces are properly seated.
- Before each procedure check the conditions of the **Neodent** surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the **Neodent** product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the **Neodent** Implant System and exempts any product warranty.
- It is the professional's responsibility to use the **Neodent** products according to the instructions for use, and to determine whether it suits the individual situation of each patient.

ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The **Neodent** Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the **Neodent** Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in its original packaging, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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