ST.330.024

Português  Attachments

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English  Attachments

Français  Attachments

Español  Attachments

Simbologia Symbology Simbología Simbologia Symbolik Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
	Requisito Brasileiro / Brazilian Requirement / Requisito Brasileño / Requisito Brasiliano / Brasilianische Anforderungen / Exigences brésiliennes
	Esterilizado por óxido de etileno / Sterile - ethylene oxide / Esterilizado con óxido de etileno / Sterile - ossido di etilene / Sterilisiert mit Ethylenoxid / Stérilisé à l'oxyde d'éthylène

Simbologia Symbology Simbología Simbologia Symbolik Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Não reutilizar / Do not reuse / No lo reutilice / Non riutilizzare / Nicht wiederverwendbar / Ne pas réutiliser / Nicht wiederverwendbar / Ne pas réutiliser
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Consérvelo seco / Conservare in luogo asciutto / Trocken lagern / Garder au sec
	Manter afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conservar à l'abri de la lumière solaire
	Não utilizar se a embalagem estiver danificada / Do not use if package is damaged / No lo utilice si el envoltorio está dañado / Non utilizzare se la confezione originale è danneggiata / Bei beschädigter Verpackung Produkt nicht verwenden / Ne pas utiliser si l'emballage a été abimé

Simbologia Symbology Simbología Simbologia Symbolik Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Marcação CE para comercialização na Comunidade Europeia / CE Mark for European Community market / Marca CE para comercialización en la Comunidad Europea / Marchio CE per la commercializzazione in Europa / CE Kennzeichen zum Vertrieb in der Europäischen Gemeinschaft / Marquage CE pour une commercialisation en Europe
	Notificação exigida pelo FDA para comercialização nos Estados Unidos / Notification required by FDA for United States market / Notificación exigida por el FDA para comercialización en los Estados Unidos / Notifica necessaria per la commercializzazione negli Stati Uniti / Vom FDA geforderte Kennzeichnung für den Vertrieb in den USA / Notification exigée par FDA pour une commercialisation aux États-Unis
	Representante na Comunidade Europeia / Representative in the European Community / Representante en la Comunidad Europea / Rappresentante nella Comunità europea / Vertreter in der Europäischen Gemeinschaft / Représentant dans la Communauté Européenne
	Fabricante / Manufacturer / Fabricante / Produttore / Hersteller / Fabricant

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

The Mini Ball and Equator Attachments are manufactured in titanium alloy according to standard ASTM F136. They receive the application of Titanium Nitride coating to increase the wear resistance. The Attachment is an intermediary Prosthetic Component between the Implant and Prosthesis (crown). It presents variations in prosthetic interfaces: Morse Taper, External and Internal Hex, Facility and different transmucosal heights, complying with the variations in gingival thickness. The Mini Ball Attachment features fitting for retaining the O'ring and hex for installation with appropriate Driver. The Equator Attachment features fitting for retaining the O'ring and an inner hex for installation with appropriate Driver.

INDICATIONS FOR USE

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading.

APPLICATIONS

This product is indicated as the interocclusal space available, transmucosal height and 3D position of the Implant. Suitable for fitting on dentures with Implants placed in the jaw or maxilla. Can be used for overdentures or partial dentures. Can be used in immediate or delayed rehabilitation procedure.

CONTRAINDICATIONS

Contraindicated as support for fixed prostheses and unsatisfactory Implant position.

HANDLING

Using the Prosthetic Abutment Driver, capture the Mini Ball Attachment directly in the package and screw the product to the Implant, applying the corresponding torque, when applicable, according to the table. For the Equator Attachment, capture directly in the package using the 1.2mm Torque Connection and make installation. For the Facility Implant, capture using the 1.2mm Torque

Connection, is recommended three hits on the Attachment for its insertion.

Prosthetic Interface	Recommended Torque
CM	32 Ncm
Ti (External Hex)	32 Ncm
II Plus (Internal Hex)	20 Ncm

TRACEABILITY LABEL

The product is identified through numerical codes arranged on the label (REF and LOT). The product is fitted with these codes, allowing their characteristics to be identified since when it is a raw material. This product is fitted with three labels which enable it to be tracked and should be posted to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor).

PRESENTATION AND STERILIZATION

It is supplied unitarily, conditioned in a blister type package (surgical grade paper and film). It is delivered sterile and it of single use. Sterilization method: ethylene oxide.

PRECAUTIONS

- Regarding the systemic aspect, consider the patient's general health conditions according to the applicable literature. Regarding the local aspect, observe the conditions of the intraoral tissues.
- This product is to be used only once and may not be re-sterilized.
- Reuse of this product may cause:
 - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
 - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if the package is damaged.
- This product should be used immediately after the opening of the package, at the time of the surgery. In case it is not used, discard it.
- A higher insertion torque than the recommended and use of improper Driver may cause damage to the material and render the system inoperable.
- Best results are obtained with the use of the list of Neodent products. The use of Prosthetic Components and/ or instruments from other manufacturers does not guarantee the perfect operation of Neodent's Dental Implant System and waives any warranty of the product.
- Ensure you do not breathe this product.
- Be careful during the waxing process for the interocclusal space is not insufficient or excessive.
- Unsuitable planning may compromise

the performance of the Implant/Prosthesis Set, resulting in system failures such as loss or fracture of the Implant, loosening or fracture of the Prosthetic Screws.

- During installation, be sure to align it to the axis of Implant insertion, avoiding crashes and damage the thread.
- Make sure the Attachment is fully seated in the Implant. For this it is recommended periapical radiographs with the parallelism technique.
- For Morse Taper Implant, is recommended the use of Try-in components.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and the use of instrument are improper.

MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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