# ST.330.025

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

### DESCRIPTION

The Neodent Space Planning Instrument has two ends for measuring, one with 7mm and another with 9mm. Each end of the Space Planning Instrument has a small tip to ensure stable handling while using the product. The end of 7mm is indicated for planning installation of Implants with distance relation between canines, premolars and central incisors. In the center of this Space Planning Instrument, there is a laser marking with the measure of 3.5mm. The end of 9mm is indicated for planning installation of Implants with distance relation between premolar/molar. In the center of this Space Planning Instrument, there is a laser marking with the measure of 4.5mm.

### APPLICATIONS

Is an instrument suitable to assist in prosthetic/surgical planning, used to determine the mesiodistal position of the Implants, by measuring the distance between Implant/Implant, Implant/tooth or Implant/nerve. Can be used on study model and/or directly in mouth. This product is manufactured in surgical stainless steel.

### CONTRAINDICATIONS

This product does not show contraindications, provided that it be properly used for the shown purposes.

### HANDLING

Select the Space Planning Instrument according to the corresponding height to the prosthetic crown that will be placed on the Implant. To install an Implant close to a tooth, touch the Space Planning Instrument on the distal/mesial face of the tooth crown and the central marking of the Space Planning Instrument will determine the perforation position. To determine the position of an Implant in relation to another Implant, starting from the first Implant, put one end of the Space Planning Instrument in the center of the first Implant and in another end mark the perforation. To determine the position of an Implant in relation to a nerve, place the end of the Space

Planning Instrument on the mentonian foramen edge and start the perforation to the Implant on the laser marking or in the another end of the Space Planning Instrument.

### SANITATION

This product should be correctly sanitized after each use. Do not leave the product in contact with humidity more than the time necessary for the sanitation. Proceed as follows: 1<sup>st</sup> step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer). 2<sup>nd</sup> step: Wash it in an ultrasonic cleaner for approximately 10 to 15 minutes. 3<sup>rd</sup> step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is allowed. 4<sup>th</sup> step: Dry it with a clean, dry cloth or with compressed air. 5<sup>th</sup> step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent - 1<sup>st</sup> step - and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence. 6<sup>th</sup> step: Select the packing in accordance with the sterilization process. Preferably, use surgical grade paper self-sealing packing with laminated film. ATTENTION: Do not use anticrust products, avoid using enzymatic liquid in concentrations above 10% and dry parts that still have cleaning solution wastes, because these procedures cause oxidation. The use of steel brushes is also inadvisable.

### PRESENTATION AND STERILIZATION

This product is unitarily provided, conditioned in blister (surgical grade film and paper type) packing. It is reusable, delivered non-sterile and should, therefore, be sterilized before use. Sterilize the product on the day before or on the day of the procedure. We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard BS EN ISO 17665-1:2006 (Sterilization of health care products. Moist heat. Requirements for the development, validation and routine control of a sterilization process for medical devices). Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Suitable parameters: 121 °C, 30 min, at 1 atm. Sterilization validity: 7 to 15 days, provided that the product is conditioned in clean, dry and protected from the sunlight environment.

### PRECAUTIONS

- Do not use the product if the package has been violated.
- This product should be used sterile.
- Better results are achieved with the utilization of the Neodent products sequence. The utilization of instrumentations and/or Prosthetic Components from other manufacturers does not assure the perfect function of the Neodent Implant System and voids any product guarantee.
- Improper planning can jeopardize the performance of the Implant/prosthesis set, causing failures to the system, such as loss or break of the Implant, loosening or fracture of the prosthetic screws.

### ADVERSE EFFECTS

Adverse effects will only occur if the choice and the use of instrument are improper.

### POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

### STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40 °C and protected from direct sunlight.

### DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

### DATE OF EXPIRATION

Written on the label.

### LIFE CYCLE

This product should be disposed of after it loses its functionality.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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