# ST.330.037

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English  Bone Profile Drill

Español  Fresa Perfil Óseo

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
	Não estéril / Non-sterile / No estéril / Non sterile / Nicht steril / Non stérile

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Consérvelo seco / Conservare in luogo asciutto / Trocken lagern / Garder au sec
	Mantener afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conservar à l'abri de la lumière solaire
	Não utilizar se a embalagem estiver danificada / Do not use if package is damaged / No lo utilice si el envoltorio está dañado / Non utilizzare se la confezione originale è danneggiata / Bei beschädigter Verpackung Produkt nicht verwenden / Ne pas utiliser si l'emballage a été abîmé

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Marcação CE para comercialização na Comunidade Europeia / CE Mark for European Community market / Marca CE para comercialización en la Comunidad Europea / Marchio CE per la commercializzazione in Europa / CE Kennzeichen zum Vertrieb in der Europäischen Gemeinschaft / Marquage CE pour une commercialisation en Europe
	Notificação exigida pelo FDA para comercialização nos Estados Unidos / Notification required by FDA for United States market / Notificación exigida por el FDA para comercialización en los Estados Unidos / Notifica necessaria per la commercializzazione negli Stati Uniti / Vom FDA geforderte Kennzeichnung für den Vertrieb in den USA / Notification exigée par FDA pour une commercialisation aux États-Unis
	Representante na Comunidade Europeia / Representative in the European Community / Representante en la Comunidad Europea / Rappresentante nella Comunità europea / Vertreter in der Europäischen Gemeinschaft / Représentant dans la Communauté Européenne
	Fabricante / Manufacturer / Fabricante / Produttore / Hersteller / Fabricant

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

### DESCRIPTION

The Bone Profile Drill is manufactured in surgical stainless steel. It has fitting for Contra-Angle in accordance to standard ISO 1797-1 - Dental rotary instruments • Shank. In the other end there is a cavity that fits to the Guide of the Bone Profile Drill with biting profile. The Guide of the Bone Profile Drill fitting for the 1.2 mm hex driver for screwing and also, thread directly adapting to the Implants.

### APPLICATIONS

The Bone Profile Drill is a surgical instrument with the purpose of preparation and regularization of the bone tissue around the Implant platform. The Bone Profile Drill is a surgical instrument with the purpose of preparation and regularization of the bone tissue around the Implant platform. The Bone Profile Drill Guide has the same function of protecting the Implant prosthetics interface and guide the Bone Profile Drill.

### CONTRAINDICATIONS

This product does not present contraindications, once they are properly used for the indicated purposes.

### HANDLING

After removing the Implant cover screw, install the Bone Profile Drill Guide corresponding to the Implant with the aid of the 1.2mm Manual Screwdriver. Fit the Drill in the Contra-Angle or in the Manual Implant Driver - Contra-Angle. Choose 200 to 300 RPM rotation. Couple the Drill over the guide and advance until it approaches the maximum to the Implant, leaning over the guide. During drilling, the pressure should not be excessive and, abundant irrigation is necessary.

### SANITATION

The surgical instruments should be correctly sanitized after each use. Do not leave the product in contact with humidity more than the time necessary for the sanitation. Proceed as follows: 1st step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer) 2nd step: Wash in ultrasonic cleaner for approximately 10 to 15 minutes. 3rd

step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is recommended. 4th step: Dry with a clean and dry cloth or with compressed air. 5th step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent - 1st step - and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence. 6th step: If applicable, place the product in an appropriate case duly sanitized and dry. 7th step: Select the packing in accordance to the sterilization process. Use preferably surgical grade paper self-sealing packing with laminated film.

### PRESENTATION AND STERILIZATION

This product is unitarily provided, conditioned in blister (surgical grade film and paper type) packing. It is reusable, delivered non-sterile and should, therefore, be sterilized before use. Sterilize the product on the day before or on the day of the procedure. We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard ISO 17665-1:2006 • Sterilization of health care products. Moist heat. Requirements for the development, validation and routine control of a sterilization process for medical devices. Suitable parameters: 121°C, 30 min, at 1 atm. Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Sterilization validity: 7 to 15 days, provided that the product is conditioned in clean, dry and protected from the sunlight environment.

### PRECAUTIONS

• This product should be sterile. • The non-replacement of the Drills, as recommended by the manufacturer, can generate improper bone heating, compromising the osseointegration success. • Do not regrind. • Do not use solutions with ammonia, acids or hydrogen peroxide, because they can damage the Drill. • Do not use anticrust products, avoid using enzymatic liquid in concentrations above 10% and dry parts that still have cleaning solution wastes, because these procedures cause oxidation. • The use of steel brushes is inadvisable. We recommend the use of nylon brushes. • The improper planning can jeopardize the performance of the Implant/prosthesis set causing failures to the system, such as loss or break of the Implant, loosening

or fracture of the prosthetic screws. • Better results are obtained with the use of the Neodent products sequence. • The use of instruments and/or prosthetic abutments of other manufacturers do not guarantee the perfect function of the Neodent Implants System and exempts any guarantee of the product.

### ADVERSE EFFECTS

Adverse effects will only occur if the choice and the use of instrument are improper.

### POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

### STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

### DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

### DATE OF EXPIRATION

Written on the label.

### LIFE CYCLE

This product has a life cycle of up to 20 drillings for bone qualities I and II, and up to 30 drillings for bone qualities III and IV, since its use conditions are respected, including adequate irrigation and rotation. Apart of the life cycle period mentioned, the professional should check the sharpening conditions of the Drill. Drills with improper cutting conditions should be disposed. Note: Classification of the bone quality according to Lekholm and Zarb (1985).

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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