# ST.330.039

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
English  Trefine Bur

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Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
	Não estéril / Non-sterile / No estéril / Non sterile / Nicht steril / Non stérile

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Consérvelo seco / Conservare in luogo asciutto / Trodden lagern / Garder au sec
	Mantem afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conservar à l'abri de la lumière solaire
	Não utilizar se a embalagem estiver danificada / Do not use if package is damaged / No lo utilice si el envoltorio está dañado / Non utilizzare se la confezione originale è danneggiata / Bei beschädigter Verpackung Produkt nicht verwenden / Ne pas utiliser si l'emballage a été abimé

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Marcação CE para comercialização na Comunidade Europeia / CE Mark for European Community market / Marca CE para comercialización en la Comunidad Europea / Marchio CE per la commercializzazione in Europa / CE Kennzeichen zum Vertrieb in der Europäischen Gemeinschaft / Marquage CE pour une commercialisation en Europe
	Notificação exigida pelo FDA para comercialização nos Estados Unidos / Notification required by FDA for United States market / Notificación exigida por el FDA para comercialización en los Estados Unidos / Notifica necessaria per la commercializzazione negli Stati Uniti / Vom FDA geforderte Kennzeichnung für den Vertrieb in den USA / Notification exigée par FDA pour une commercialisation aux États-Unis
	Representante na Comunidade Europeia / Representative in the European Community / Representante en la Comunidad Europea / Rappresentante nella Comunità europea / Vertreter in der Europäischen Gemeinschaft / Représentant dans la Communauté Européenne
	Fabricante / Manufacturer / Fabricante / Produttore / Hersteller / Fabricant

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

### DESCRIPTION

The Trepine Bur is manufactured in surgical stainless steel. Its active end has a high cutting capacity. Its inner section has a cavity which allows the collection of the implant to be removed and also of the bone tissue block. Along the Bur, there are channels which allow the block to be removed after the procedure. In the other end there is fitting for Contra-Angle in accordance to standard ISO 1797-1 - Dental rotary instruments • Shank. Laser markings determine the drilling depth in accordance with the surgical planning.

### APPLICATIONS

The Trepine Bur is a surgical instrument used to remove dental Implants and bone tissue blocks for grafting.

### CONTRAINDICATIONS

This product does not present contraindications, once they are properly used for the indicated purposes.

### HANDLING

Fit the Bur in the Contra-Angle. Select rotation from 500 to 800 RPM. Perform the drilling with complete insertion and removal movements of the Bur. During drilling, the pressure should not be excessive and, an abundant irrigation, manual or combined, using the motor pump is necessary.

- For Implant removal: Select the Bur according to the diameter of the Implant to be removed. Insert the Bur so that it involves the implant, piercing the bone up to the planned depth.
- For block removal: Select the Bur diameter according to the diameter of the desired bone tissue block. Drill in accordance with the desired amount of bone. Note: Do not interrupt the motor rotation with the Bur within the surgical cavity, because this can make difficult the removal or cause the Bur break.

### SANITATION

The surgical instruments should be correctly sanitized after each use. Do not leave the product in contact with humidity more than the time necessary for the sanitation. Proceed as follows: 1st step: Fully immerse the part in

enzymatic detergent (diluted according to the manufacturer). 2nd step: Wash it in an ultrasonic cleaner for approximately 10 to 15 minutes. 3rd step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is recommended. 4th step: Dry it with a clean and dry cloth or with compressed air. 5th step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent - 1st step - and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence. 6th step: If applicable, place the product in an appropriate case duly sanitized and dry. 7th step: Select the packing in accordance to the sterilization process. Preferably, use surgical grade paper self-sealing packing with laminated film.

### PRESENTATION AND STERILIZATION

This product is unitarily provided, conditioned in blister (surgical grade film and paper type) packing. It is reusable, delivered non-sterile and should, therefore, be sterilized before use. Sterilize the product on the day before or on the day of the procedure. We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard ISO 17665-1:2006 • Sterilization of health care products. Moist heat. Requirements for the development, validation and routine control of a sterilization process for medical devices. Suitable parameters: 121 °C, 30 min, at 1 atm. Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Sterilization validity: 7 to 15 days, provided the product is conditioned in clean, dry and protected from the sunlight environment.

### PRECAUTIONS

- Do not use the product if the packing is damaged.
- This product should be used sterile.
- The non-replacement of the Burs, as recommended by the manufacturer, can generate improper bone heating, compromising the osseointegration success.
- Do not regrind.
- Do not use solutions with ammonia, acids or hydrogen peroxide, because they can damage the Bur.
- Do not use anticrust products, avoid using enzymatic liquid in concentrations above 10% and dry parts that still have cleaning solution wastes, because these procedures cause oxidation.
- The use of steel brushes is inadvisable. We recommend the use of nylon brushes.
- Improper planning can jeopardize the performance of the

implant/prosthesis set causing failures to the system, such as loss or break of the Implant, loosening or fracture of the prosthetic screws.

- Better results are obtained with the use of the Neodent products sequence.
- The use of instruments and/or prosthetic abutments of other manufacturers do not guarantee the perfect function of the Neodent Implants System and exempts any guarantee of the product.

### ADVERSE EFFECTS

Adverse effects will only occur if the choice and the use of instrument are improper.

### POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

### STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

### DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

### DATE OF EXPIRATION

Written on the label.

### LIFE CYCLE

This product has a life cycle of up to 20 drillings for bone qualities I and II, and up to 30 drillings for bone qualities III and IV, since its use conditions are respected, including adequate irrigation and rotation. Regardless of the life cycle period mentioned, the professional should always check the sharpening conditions of the Bur. Burs with improper cutting conditions should be disposed of. Note: Classification of bone quality according to Lekholm & Zarb (1985).

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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