# ST.330.044

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

### DESCRIPTION

This product has cylindrical shape and anatomic handle, favoring its grip. At its tip, there is a fast coupling for the fitting of the Phillips Connection for Manual ScrewDriver or the Punch for Bone Grafting and Orthodontic Anchorage.

### APPLICATIONS

This product is used for manual fixing of Screws used in cases of bone grafting in block with the Philips Connection for Manual ScrewDriver. It is also used to demarcate the insertion site of the Self-Drilling Graft Screw or the Self-Drilling Anchorage Implant with the Punch for Bone Grafting and Orthodontic Anchorage.

### CONTRAINDICATIONS

This product has no contraindications, providing its intended purposes and indications are correctly followed.

### HANDLING

Hold the Manual Driver and pull the fast coupling extremity backwards with the tip of your fingers. With your other hand, fit the Phillips Connection or the Punch for Bone Grafting and Orthodontic Anchorage on the Manual Driver. Loosen the extremity and make sure that the Connection / Punch is attached to the Manual Driver. Hold the Driver in the palm of your hand and turn the extremity with the fingers of the same hand.

### SANITATION

This product is reusable and should be correctly cleaned after each use. Follow these steps for cleaning: 1<sup>st</sup> step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer). 2<sup>nd</sup> step: Wash it in an ultrasonic

cleaner for approximately 10 to 15 minutes. 3<sup>rd</sup> step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is allowed. 4<sup>th</sup> step: Dry it with a clean, dry cloth or with compressed air. 5<sup>th</sup> step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent - 1<sup>st</sup> step - and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence. 6<sup>th</sup> step: Select the packing in accordance with the sterilization process. Preferably, use surgical grade paper self-sealing packing with laminated film. ATTENTION: Do not use antitrust products, avoid drying parts still containing cleaning solution wastes because such procedures favor the oxidation. The use of steel brushes is not recommended.

### PRESENTATION AND STERILIZATION

This product is reusable and supplied unitarily, non-sterile, conditioned in envelope (film type) packing and, therefore, should be sterilized before use. Sterilize the product on the day before or on the day of the procedure. We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard BS EN ISO 17665-1:2006 (Sterilization of health care products. Moist heat. Requirements for the development, validation and routine control of a sterilization process for medical devices). Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Suitable parameters: 121°C, 30 min, at 1 atm. Sterilization validity: 7 to 15 days, provided that the product is conditioned in clean, dry and protected from the sunlight environment.

### PRECAUTIONS

- Do not use the product if the package has been violated.
- Do not sterilize the product in its original package.
- Better results are achieved with the utilization of the Neodent products

sequence. The utilization of instrumentations and/or prosthetic components from other manufacturers does not assure the perfect function of the Neodent Implant System and voids any product guarantee.

### ADVERSE EFFECTS

No adverse effects will occur, since the indications of this product be respected.

### POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

### STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

### DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

### DATE OF EXPIRATION

Written on the label.

### LIFE CYCLE

This product should be disposed of after it loses its functionality.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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