ST.330.050

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English  Healing Abutment

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

It is a prosthetic component for temporary use, made of titanium alloy, in accordance with ASTM standard F136. On one end, it has a hex coupling (0.9 or 1.2 mm) and/or slot for its installation, and on the other, it has an interface compatible with the various implant prosthetic platforms, according to the table below.

Prosthetic Interface		Diameter (mm)	Transmucosal height (mm)
Morse Taper		3.3, 4.5	0.8, 1.5, 2.5, 3.5, 4.5, 5.5, 6.5
External Hex	Divergent	3.3, 4.1, 4.3, 5.0	2.0, 3.0, 4.0, 5.0
	SlimFIT	4.1	2.0, 3.0, 4.0, 5.0, 6.0
		5.0	2.0, 3.0, 4.0, 5.0
	Parallel	3.3, 4.3	2.0, 3.0, 4.0, 5.0
5.0		2.0, 3.0, 4.0	
Internal Hex		N/A	3.0, 4.0, 5.0
Zygomatic and WS		N/A	0.8, 1.5, 2.5, 3.5
Facility		N/A	1.5, 2.5, 3.5, 4.5

INDICATIONS FOR USE

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading.

APPLICATIONS

It is indicated for maintenance of the soft tissues during the osseointegration phase of Neodent implant to be rehabilitated using the late loading technique.

It can be used on the insertion of the Implant or re-opening surgery (second surgical phase). It is indicated as interocclusal space available, existing transmucosal height and three-dimensional position of the Implant.

The Facility 1.5 mm healing abutment can be used as a cover screw for the Facility implant.

CONTRAINDICATIONS

This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: Titanium alloy Ti6Al4V.

On implants with installation torque below 10 N.cm, use of the healing abutment on the same surgical occasion is contraindicated, except for use of the Facility healing abutment as a cover screw.

This product is not indicated for customization.

HANDLING

After opening the package, place the Healing Abutment on a sterile surface and pick it up with the aid of the Manual Screwdriver 0.9 or 1.2. Install it on the implant, with maximum torque of 10 N.cm. For the Facility Implant, we recommend one tap with the Facility Mallet on the Healing Abutment for its installation. When installed during the same surgical procedure as the installation of the implant, suturing must be performed after the mounting of the healing abutment.

Note: the Divergent 3.3 External Hex and Parallel healing abutments with a height of 2.0 are the only ones which use the 0.9 Driver for their installation.

PRESENTATION AND STERILIZATION

This product is indicated for single use and provided sterilized by the ethylene oxide method, packed individually into a package that offers dual protection.

PRECAUTIONS

- It is recommended that the device be removed within 6 months after its installation.
- During the installation, make sure to align it to the insertion axis of the implant, preventing jams and damage to the thread.
- This product should not be used as support for removable dentures.
- Note the conditions of the intra-oral tissue, the bone quality and quantity of the bed receiving the Implant, by means of radiographic and/or tomography examinations. Nonperformance of the pre-surgical assessment may compromise the success of the procedure.
- As for the systemic aspect, consider the general health of the patient. In particular, one must be careful in cases of patients who have allergies to drugs, local or systemic factors that may interfere with the healing process of the bone tissues or soft tissues, or the process of osseointegration. For example, the bone already exposed to radiation in the head and neck area, diabetes mellitus, anticoagulation drugs, hemorrhagic diathesis, bruxism, parafunctional habits, anatomically unfavorable bone situation, tobacco abuse, uncontrolled periodontitis, treatable jaws pathologies and oral mucosa abnormalities.
- Treatment with bisphosphonates results in potential risk of peri-implant osteonecrosis.
- Surgical planning can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of components and/or prosthetic screws.
- The material to be used during the procedure must be sterile.
- This product is of single use and cannot be re-sterilized.
- Reuse of this product may cause:
 - adverse biological effects of residual products,

microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.

- Do not use the product if the packaging is broken.
- The sterilization is guaranteed only if the blister packaging is not damaged.
- This product must be used immediately after opening the packaging, at the moment of surgery. If it is not used, discard it.
- Do not use the product with the validity expired.
- Insertion torque greater than that recommended can make the system inoperative.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Some of the possible causes of unsuccessful bone integration and loss of the prosthesis during the treatment are the following: unsuitable osteotomy, infections, deficient oral hygiene, occlusal trauma, systemic problems or diseases, low remaining bone quantity or quality, lack or failure of irrigation, use of instruments which are not specific and/or without cutting power and lack of specific training.
- Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the Neodent product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the professional's responsibility to use the Neodent products according to the instructions for use.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and the use of instrument are improper.

MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in its original packaging, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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