ST.330.057

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English  Protection Cylinder

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Simbologia Symbology Simbologia Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
	Requisito Brasileiro / Brazilian Requirement / Requisito Brasileño / Requisito Brasiliano / Brasilianische Anforderungen / Exigences brésiliennes
	Esterilizado por óxido de etileno / Sterile - ethylene oxide / Esterilizado con óxido de etileno / Sterile - ossido di etilene / Sterilisiert mit Ethylenoxid / Stérilisé à l'oxyde d'éthylène

Simbologia Symbology Simbologia Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Não reutilizar / Do not reuse / No lo reutilice / Non riutilizzare / Nicht wiederverwendbar / Ne pas réutiliser / Nicht wiederverwendbar / Ne pas réutiliser
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Consérvelo seco / Conservare in luogo asciutto / Trocken lagern / Garder au sec
	Manter afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conservar à l'abri de la lumière solaire
	Não utilizar se a embalagem estiver danificada / Do not use if package is damaged / No lo utilice si el envoltorio está dañado / Non utilizzare se la confezione originale è danneggiata / Bei beschädigter Verpackung Produkt nicht verwenden / Ne pas utiliser si l'emballage a été abîmé

Simbologia Symbology Simbologia Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Marcação CE para comercialização na Comunidade Europeia / CE Mark for European Community market / Marca CE para comercialización en la Comunidad Europea / Marchio CE per la commercializzazione in Europa / CE Kennzeichen zum Vertrieb in der Europäischen Gemeinschaft / Marquage CE pour une commercialisation en Europe
	Notificação exigida pelo FDA para comercialização nos Estados Unidos / Notification required by FDA for United States market / Notificación exigida por el FDA para comercialización en los Estados Unidos / Notifica necessaria per la commercializzazione negli Stati Uniti / Vom FDA geforderte Kennzeichnung für den Vertrieb in den USA / Notification exigée par FDA pour une commercialisation aux États-Unis
	Representante na Comunidade Europeia / Representative in the European Community / Representante en la Comunidad Europea / Rappresentante nella Comunità europea / Vertreter in der Europäischen Gemeinschaft / Représentant dans la Communauté Européenne
	Fabricante / Manufacturer / Fabricante / Produttore / Hersteller / Fabricant

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

Neodent titanium Prosthetic Components are used as intermediaries between the implant and the prosthesis, planned according to each case, respecting the prosthetic interfaces, aiming for aesthetics and function.

DESCRIPTION

The Protection Cylinder is a prosthetic component with a tapered shape, made of titanium alloy according to ASTM standard F136, installed on prosthetic components to protect them during the fabrication of the prosthesis and/or healing of the peri-implant tissues. It has a coupled screw for fastening it onto the corresponding prosthetic component.

Specifically for the GT implants, the Protection Cylinder should be installed directly on the implant.

It is available as shown in the table below:

Prosthetic interface	Available dimensions
Transepithelial Pillar	Single
Mini Conical Pillar	Ø4.1 and Ø 5.0
Conical Pillar	Ø4.1, Ø4.1 (Compatible) and Ø 5.0
GT / CM Pillar	Single
Micro Pillar	Single

INDICATIONS FOR USE

For conventional implants

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading. Multiple teeth applications can be rigidly splinted.

For Zygomatic implants

Zygomatic Implants are indicated for surgical installation in the zygoma region, in cases of severe jaw resorption, in order to restore patient esthetics and chewing function. Zygomatic Implants are recommended for the posterior (pre-molar/molar) region, one implant on each side, with at least two standard dental implants in the anterior region to support a fixed restoration. Zygomatic Implants may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading.

APPLICATIONS

It is a device for short-term use, indicated for protecting the Implant/Prosthetic component assembly during fabrication of the prosthesis and/or healing of the peri-implant tissues. Indicated after installation of intermediary prosthetic components: Conical Pillar, Mini Conical Pillar, CM Pillar, Micro Conical Pillar or Transepithelial Pillar, or GT Implant.

CONTRAINDICATIONS

This product is contraindicated for patients who present signs of allergy or hypersensitivity to the chemical components of the material: titanium alloy Ti6Al4V-ELI.

This product is not indicated for customization.

HANDLING

After opening the package, pick up the protection cylinder with the 1.2 hex wrench and install it onto the corresponding prosthetic component/implant, screwing it clockwise with the recommended torque of 10 N.cm. During the installation, ensure the correct fit between the protection cylinder and the prosthetic component/implant.

TRACEABILITY LABEL

This product is accompanied by three labels that allow its traceability and should be attached to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

PRESENTATION AND STERILIZATION

This product is indicated for single use and provided sterilized by the ethylene oxide method, packed individually into a package that offers dual protection.

PRECAUTIONS

- During the installation, make sure to align it to the insertion axis, avoiding jamming and damage to the threads.
- Inappropriate planning can jeopardize the performance of the implant/prosthesis unit leading to failures of the system, as loss or fracturing of the implant, loosening or fracturing of the prosthetic screws.
- This product is of single use and cannot be re-sterilized.
- Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if the packaging is broken.
- The sterilization is guaranteed only if the blister packaging is not damaged.
- This product must be used immediately after opening the packaging, at the moment of surgery. If it is not used, discard it.
- Do not use the product with the validity expired.

- See the torque to be given on the prosthetic component to be used. The excess or lack of torque can lead to undesirable result. Insertion torque greater than that recommended can make the system inoperative.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the Neodent product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the professional's responsibility to use the Neodent products according to the instructions for use, and to determine whether it suits the individual situation of each patient.

ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are respected.

MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored, in its original packaging, in a clean and dry location, in a maximum temperature of 40 °C and protected from direct sunlight.

DISPOSAL OF MATERIAL

After use, all the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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