# ST.330.100


Português  Medidor de Altura

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English  Height Measurer

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Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
	Não estéril / Non-sterile / No estéril / Non sterile / Nicht steril / Non stérile

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Consérvelo seco / Conservare in luogo asciutto / Trocken lagern / Garder au sec
	Manter afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conserver à l'abri de la lumière solaire
	Não utilizar se a embalagem estiver danificada / Do not use if package is damaged / No lo utilice si el envoltorio está dañado / Non utilizzare se la confezione originale è danneggiata / Bei beschädigter Verpackung Produkt nicht verwenden / Ne pas utiliser si l'emballage a été abimé

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Marcação CE para comercialização na Comunidade Europeia / CE Mark for European Community market / Marca CE para comercialización en la Comunidad Europea / Marchio CE per la commercializzazione in Europa / CE Kennzeichen zum Vertrieb in der Europäischen Gemeinschaft / Marquage CE pour une commercialisation en Europe
	Notificação exigida pelo FDA para comercialização nos Estados Unidos / Notification required by FDA for United States market / Notificación exigida por el FDA para comercialización en los Estados Unidos / Notifica necessaria per la commercializzazione negli Stati Uniti / Vom FDA geforderte Kennzeichnung für den Vertrieb in den USA / Notification exigée par FDA pour une commercialisation aux États-Unis
	Representante na Comunidade Europeia / Representative in the European Community / Representante en la Comunidad Europea / Rappresentante nella Comunità europea / Vertreter in der Europäischen Gemeinschaft / Représentant dans la Communauté Européenne
	Fabricante / Manufacturer / Fabricante / Produttore / Hersteller / Fabricant

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

## DESCRIPTION

The Neodent Height Measurer is manufactured in titanium alloy according to standard ASTM F136, has a pin shape with conical tip and a hole in its upper extremity to crossing safety floss (dental floss). It's available in three models: -CM Height Measurer - It has markings along the body, corresponding to the height of CM Prosthetic Components that will be used (heights: 0.8, 1.5, 2.5, 3.5, 4.5, 5.5 and 6.5mm). -CM Zygomatic Height Measurer - It has markings along the body, corresponding to the height of CM Zygomatic Prosthetic Components that will be used (heights: 3.0, 4.0, 5.0, 6.0mm). -Facility Height Measurer - It has recesses along the body, corresponding to the height of Facility Prosthetic Components that will be used (heights: 1.5, 2.5, 3.5, 4.5 e 5.5mm). This product can also be used as an X-Ray Positioner.

## APPLICATIONS

This product is an instrument used in dental surgery to perform the selection of mucosa height (gingiva) of the Prosthetic Components which will be installed over the CM, Zygomatic and Facility Implants.

## CONTRAINDICATIONS

This product has no contraindications, providing its intended purposes and indications are correctly followed.

## HANDLING

After install the Implant and suture the gingiva, introduce the Height Measurer by its conical area in the internal cone of Implant and verify the mucosal height, selecting thus the respective height of the components that will be used.

## SANITATION

It should be correctly sanitized after each use. Do not leave the product in contact with humidity more than the time necessary for the

sanitation. Proceed as follows: 1<sup>st</sup> step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer). 2<sup>nd</sup> step: Wash in ultrasonic cleaner for approximately 10 to 15 minutes. 3<sup>rd</sup> step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is allowed. 4<sup>th</sup> step: Dry it with a clean, dry cloth or with compressed air. 5<sup>th</sup> step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent - 1<sup>st</sup> step - and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence. 6<sup>th</sup> step: If applicable, place the product in an appropriate case, duly sanitized and dry. 7<sup>th</sup> step: Select the packing in accordance with the sterilization process. Preferably, use surgical grade paper self-sealing packing with laminated film. ATTENTION: Do not use antitrust products, avoid using enzymatic liquid in concentrations above 10% and dry parts that still have cleaning solution wastes, because these procedures cause oxidation. Use of steel brushes is also inadvisable.

## PRESENTATION AND STERILIZATION

This product is unitarily provided, conditioned in blister (surgical grade film and paper type) packing. It is reusable, delivered non-sterile and should, therefore, be sterilized before use. Sterilize the product on the day before or on the day of the procedure. We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard BS EN ISO 17665-1:2006 (Sterilization of health care products. Moist heat. Requirements for the development, validation and routine control of a sterilization process for medical devices). Suitable parameters: 121°C, 30 min, at 1 atm. Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Sterilization validity: 7 to 15 days, provided the product is conditioned in clean, dry and protected from sunlight environment.

## PRECAUTIONS

- Do not use the product if the package has been violated.
- This product should be

sterile. • Better results are achieved with the utilization of the Neodent products sequence. The utilization of instrumentations and/or Prosthetic Components from other manufacturers does not assure the perfect function of the Neodent Implant System and voids any product guarantee. • Improper planning can jeopardize the performance of the Implant/Prosthesis set, causing failures to the system, such as loss or break of the Implant, loosening or fracture of the Prosthetic Screws.

## ADVERSE EFFECTS

Adverse effects will only occur if the choice and the use of instrument are improper.

## POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

## STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

## DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

## DATE OF EXPIRATION

Written on the label.

## LIFE CYCLE

This product should be disposed of after losing its functionality.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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