# ST.330.102

Português 

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions.

Neodent titanium Prosthetic Components are used as intermediaries between the implant placed intraorally and the prosthesis, planned according to each case, respecting the prosthetic interfaces, aiming for aesthetics and function.

## DESCRIPTION

Intermediary prosthetic component to be installed between the Implant and the prosthesis (crown), made of titanium alloy according to standard ASTM F136. They are provided in a rotational form on the coupling with the prosthesis, in different prosthetic interfaces and in different transmucosal heights to match the variations in mucosal thickness. On the upper portion, it has an internal thread for subsequent mounting of the prosthetic screw, and a hex for its installation with an appropriate connection. On the External Hex interface, platform 3.3 and 4.3, the prosthetic component is accompanied by a specific screw for installation on the implant with an appropriate prosthetic connection. They are available according to the table below:

Description / Prosthetic Interface	Prosthetic platform	Transmucosal height (mm)
Mini Abutment	External Hex	3.3
		4.3
	External Hex Angled 17°	3.3, 4.1, 4.3, and 5.0
	External Hex Angled 30°	4.1, 4.3, and 5.0
	Slim Fit External Hex	4.1
		5.0
	II (Internal Hex)	4.3
		5.0
	II Plus	Single
	CM	Single
Micro Pillar	CM and CM Exact Angled 17° and 30°	Single
	WS	Single
	CM	Single
	Facility	Single

## INDICATIONS FOR USE

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading.

## APPLICATIONS

It is indicated for single-tooth prostheses screw-mounted onto implants, according to the available interocclusal space, transmucosal height, and three-dimensional position of the implant. It can be used in immediate or conventional rehabilitation procedures on the maxilla or mandible.

The CM Micro Conical Abutment is indicated for implants close to one another.

## CONTRAINDICATIONS

The Facility Micro Conical Abutment is contraindicated for final multiple-teeth applications in the region of upper incisors, canines, pre-molars, and molars.

This product is not indicated for customization.

This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: Titanium alloy Ti6Al4V.

Contraindicated for insufficient interocclusal space and unsatisfactory three-dimensional position of the implant.

## HANDLING

For use of the Mini Conical Abutment, CM Micro Conical Abutment, and Facility Micro Conical Abutment in two-stage procedures, a prior preparation can be done on the soft tissues with the use of a healing abutment.

Select the prosthetic component according to the planning and install it on the implant according to the recommended torque and connection (table below). Ensure the proper fit and install the compatible Protective Cylinder.

For the impression procedure, fit the corresponding Transfer onto the component, ensure the proper fit, and take the impression using appropriate materials. Prepare the plaster model. Prepare the prosthesis, using the Cylinders of the Mini / Micro Conical Abutment according to appropriate laboratory techniques. Tests should be carried out on passivity and the fit of the structure of the prosthesis.

For installation of the prosthesis, remove the Protective Cylinder and install it with torque of 10 N.cm on the Prosthetic Screw with the 1.2 Torque Connection. Protect the access for the screw with teflon and resin compound).

Component	Torque connection (mm)	Recommended Torque (N.cm)
CM and Exact Angled Mini Conical Abutment	0.9	15
Mini Conical Abutment (ø3.3)	Prosthetic Abutments	20
Mini Conical Abutment II Plus		
Mini Conical Abutment II 4.3		
Slim Fit Mini Conical Abutment (ø4.1)		
Mini Conical Abutment (ø4.3)		
CM Mini Conical Abutment	1.2	20
WS Mini Conical Abutment		
Angled Mini Conical Abutment (ø4.1)		
Mini Conical Abutment II 5.0		
Slim Fit Mini Conical Abutment (ø5.0)	1.6	32
Facility Micro Conical Abutment	Facility Abutment Placement Aid	3 taps*

\*NOTE: In the case of the Facility Micro Conical Abutment, its attachment to the Implant must be done using axial impact (tap) on the component, along the axis of the Implant.

## TRACEABILITY LABEL

This product is accompanied by three labels that allow its traceability and should be attached to the following documents:

- Medical record;
- Collection tax document;
- Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

## PRESENTATION AND STERILIZATION

This product is indicated for single use and provided sterilized by the ethylene oxide method, packed individually into a package.

## PRECAUTIONS

- It is clinically advisable to maintain the minimum height of 4 mm above the transmucosal height.
- In the case of the Morse Taper interface, it is recommended to use the CM Prosthetic Selection Kit. In the case of the WS interface, it is recommended to use the CM prosthetic selection kit, with trunnions compatible with the WS prosthetic interface. In the case of the Facility interface, the components for prosthetic selection available in the Facility Kit are recommended. CM Prosthetic Components are not compatible with the Facility and WS prosthetic interfaces.
- Surgical planning and/or inadequate prosthesis can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of components and/or prosthetic screws.
- The selection of the material of the structure of the prosthesis must take local aspects of the patient into account.
- This product is of single use and cannot be re-sterilized.
- Reuse of this product may cause:
  - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
  - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if the packaging is broken.
- The sterilization is guaranteed only if the packaging is not damaged.

- This product must be used immediately after opening the packaging, at the moment of the procedure. If it is not used, discard it.
- Do not use the product after expiry date.
- For immediate load application, check the torque indication of the implant installed.
- See the torque to be given on the prosthetic component to be used. The excess or lack of torque can lead to undesirable result.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- During installation, make sure to align it with the insertion axis of the implant. Make sure it is perfectly seated on the Implant. To do so, it is recommended to use periapical X-rays with the parallel technique.
- Make sure to use a prosthetic screw that is compatible with the prosthetic interface and with the prosthetic component.
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/prosthesis assembly.
- Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the Neodent product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the professional's responsibility to use the Neodent products according to the instructions for use.

## ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

## MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

## POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the procedure and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

## STORAGE CONDITIONS

This product should be stored, in its original packaging, in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

## DISPOSAL OF MATERIAL

After use, all the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

## DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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