ST.330.108

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
English  Graft Screw

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Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
	Requisito Brasileiro / Brazilian Requirement / Requisito Brasileño / Requisito Brasileiro / Brasilianische Anforderungen / Exigences brésiliennes
	Esterilizado por óxido de etileno / Sterile - ethylene oxide / Esterilizado con óxido de etileno / Sterile - ossido di etilene / Sterilisiert mit Ethylenoxid / Stérilisé à l'oxyde d'éthylène

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Não reutilizar / Do not reuse / No lo reutilice / Non riutilizzare / Nicht wiederverwendbar / Ne pas réutiliser / Nicht wiederverwendbar / Ne pas réutiliser
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Consérvelo seco / Conservare in luogo asciutto / Trocken lagern / Garder au sec
	Manter afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conservar à l'abri de la lumière solaire
	Não utilizar se a embalagem estiver danificada / Do not use if package is damaged / No lo utilice si el envoltorio está dañado / Non utilizzare se la confezione originale è danneggiata / Bei beschädigter Verpackung Produkt nicht verwenden / Ne pas utiliser si l'emballage a été abîmé

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Marcação CE para comercialização na Comunidade Europeia / CE Mark for European Community market / Marca CE para comercialización en la Comunidad Europea / Marchio CE per la commercializzazione in Europa / CE Kennzeichen zum Vertrieb in der Europäischen Gemeinschaft / Marquage CE pour une commercialisation en Europe
	Notificação exigida pelo FDA para comercialização nos Estados Unidos / Notification required by FDA for United States market / Notificación exigida por el FDA para comercialización en los Estados Unidos / Notifica necessaria per la commercializzazione negli Stati Uniti / Vom FDA geforderte Kennzeichnung für den Vertrieb in den USA / Notification exigée par FDA pour une commercialisation aux États-Unis
	Representante na Comunidade Europeia / Representative in the European Community / Representante en la Comunidad Europea / Rappresentante nella Comunità europea / Vertreter in der Europäischen Gemeinschaft / Représentant dans la Communauté Européenne
	Fabricante / Manufacturer / Fabricante / Produttore / Hersteller / Fabricant

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

This product has cylindrical format, self-drilling geometry and conical head with philips mortise. This product is manufactured in titanium alloy in accordance with ASTM F136 standard.

INDICATIONS FOR USE

This product is an implantable device used for fixation of bone blocks.

APPLICATIONS

This product serves for fixing bone blocks.

CONTRAINDICATIONS

The surgical technique involved in the use indication of this product it is not recommended in the presence of inflammatory or infectious processes of the intraoral tissues.

HANDLING

Select the length and the Graft Screw diameter according to the thickness of the bone block and the receiving region, in order to stabilize it. To install the Screw in bone densities grades III and IV, it is not necessary perforate the receiving region with Drills. Do the perforation in the block with a Drill 1.6, and the Screw is fixed on the receiving region. To install the Screw in bone densities grades I and II, do the perforation with a Drill 1.6 in the block and, depending on the bone cortical thickness of the receiving region, use the Drill 1.1 or 1.3 before placing the Graft Screw. The Screw with diameter 2.0 should only be used in cases where the Screw 1.5 it was not sufficient to stabilize the block. In this case, it is not necessary to do a previous perforation. The Screw with expanded head is indicated for donor blocks originated from the iliac crest (with inferior bone quality). Note: The surgical perforation diameter in the bed receptor should always be less than the Graft Screw diameter selected, so as to guarantee the initial stability of the fixation. In case of gingival graft, it is recommended to use the Graft Screw (wide) without pre-drilling of the tissue being fixed and receiving area. Should avoid areas of ischemia in the graft after tightening the Screw.

TRACEABILITY LABEL

The product is identified through numerical codes arranged on the label (REF and LOT). The product is fitted with these codes, allowing their characteristics to be identified since when it is a raw material. This product is fitted with three labels which enable it to be tracked and should be posted to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor).

PRESENTATION AND STERILIZATION

This product is protected by a blister packaging (film and medical grade paper) and supplied by unit. This product is supplied sterile and intended for single use. Sterilization method: ethylene oxide.

PRECAUTIONS

- Graft Screws may interfere in MRIs and CT Scans (image distortion). It is important that the professional warns their patients about this.
- Regarding the systemic aspect, consider the patient's general health conditions according to the applicable literature. Regarding the local aspect, observe the conditions of the intraoral tissues.
- This product is to be used only once and may not be re-sterilized.
- Reuse of this product may cause:
 - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
 - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
 - Do not use the product if the package is damaged.
 - This product must be used immediately after opening the packaging, at the moment of surgery. If it is not used, discard it.
 - A higher insertion torque than the recommended and use of improper Driver may cause damage to the material and render the system inoperable.
 - Best results are obtained with the use of the list of Neodent products. The use of Prosthetic Components and/or instruments from other manufacturers does not guarantee the perfect operation of Neodent's Dental Implant System and waives any warranty of the product.
 - Ensure you do not breathe this product.
 - The improper planning may endanger the Implant/prosthesis set performance resulting in system faults, such as Implant loss or fracturing, loosening or fracturing of the Prosthetic Screws.

ADVERSE EFFECTS

There may be mild discomfort and located edema.

MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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