# ST.330.120

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
English  Direction Indicator

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Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
	Não estéril / Non-sterile / No estéril / Non sterile / Nicht steril / Non stérile

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Consérvelo seco / Conservare in luogo asciutto / Trodden lagern / Garder au sec
	Manter afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conserver à l'abri de la lumière solaire
	Não utilizar se a embalagem estiver danificada / Do not use if package is damaged / No lo utilice si el envoltorio está dañado / Non utilizzare se la confezione originale è danneggiata / Bei beschädigter Verpackung Produkt nicht verwenden / Ne pas utiliser si l'emballage a été abîmé

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Marcação CE para comercialização na Comunidade Europeia / CE Mark for European Community market / Marca CE para comercialización en la Comunidad Europea / Marchio CE per la commercializzazione in Europa / CE Kennzeichen zum Vertrieb in der Europäischen Gemeinschaft / Marquage CE pour une commercialisation en Europe
	Notificação exigida pelo FDA para comercialização nos Estados Unidos / Notification required by FDA for United States market / Notificación exigida por el FDA para comercialización en los Estados Unidos / Notifica necessaria per la commercializzazione negli Stati Uniti / Vom FDA geforderte Kennzeichnung für den Vertrieb in den USA / Notification exigée par FDA pour une commercialisation aux États-Unis
	Representante na Comunidade Europeia / Representative in the European Community / Representante en la Comunidad Europea / Rappresentante nella Comunità europea / Vertreter in der Europäischen Gemeinschaft / Représentant dans la Communauté Européenne
	Fabricante / Manufacturer / Fabricante / Produttore / Hersteller / Fabricant

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

### DESCRIPTION

The Neodent Direction Indicators present cylindrical shape with different tip diameters separated by a ring with a larger diameter. This ring has a small hole for flossing and safety. These products are manufactured in titanium alloy according to ASTM F136.

### APPLICATIONS

The Neodent Direction Indicator is a surgical instrument used to verify the perforation angulation. The Neodent Neopronto Direction Indicator is a surgical instrument used to verify perforation position in relation to antagonist arch.

### CONTRAINDICATIONS

This product does not show contraindications, provided that it is properly used for the aimed purposes.

### HANDLING

Before using the Neodent Direction Indicator, it is recommended tie a floss on the central orifice of the product to avoid displacement and risk of swallowing by the patient. After drilling with Drill 2.0mm, the Direction Indicator is inserted into the perforation by the smaller tip to verify the perforation angulation. After drilling with Drill 3.0mm, it is introduced the Direction Indicator into the larger tip, verifying again the perforation angulation. When used with CM Implants, the middle portion also indicates the diameter of the cervical portion of the Implant. The Neodent Neopronto Direction Indicator is used with the Neopronto Surgical Guide to verify the parallelism between Implants and occlusion. After fitting the Surgical Guide, the Direction Indicator is inserted into the perforation by the smaller tip.

### SANITATION

The surgical instrumentation should be properly sanitized after its utilization. Do not let the product contact moisture beyond the time required for the sanitation. For such, proceed as follows: Step 1: Submerge the parts fully into

enzymatic detergent (diluted according to the manufacturer). Step 2: Wash in ultrasonic washer for approximately 10 to 15 minutes. Step 3: Rinse with plenty of distilled water until the solution wastes are fully removed. The utilization of nylon brushes is recommended. Step 4: Dry with a clean and dry cloth or with compressed air. Step 5: Perform a visual inspection, watching if there are faults in the cleaning process. If any wastes are still left, the part should be submerged again in detergent - step 1 - and if necessary, the cleaning should be made with the help of a nylon brush. Repeat the rinsing and drying sequence. Step 6: If applicable, place the product in a suitable case duly sanitized and dried. Step 7 - Select the package according to the sterilization process. Use preferably a self-sealing package if surgical grade paper with laminated film. ATTENTION: Do not use descaling products and avoid drying parts still containing cleaning solution wastes because such procedures favor oxidation. Also, the use of steel brushes is not recommended.

### STERILIZATION

This product is unitarily provided, conditioned in blister (surgical grade film and paper type) packing. It is reusable, delivered non-sterile and should, therefore, be sterilized before use. Sterilize the product on the day before or on the day of the procedure. We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard ISO 17665-1:2006 - Sterilization of health care products. Moist heat. Requirements for the development, validation and routine control of a sterilization process for medical devices. Suitable parameters: 121°C, 30 min, at 1 atm. Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Sterilization validity: 7 to 15 days, provided that the product is conditioned in a clean, dry environment away from the sunlight

### PRESENTATION AND STERILIZATION

It is supplied unitarily, conditioned in a blister type package (surgical grade paper and film) and delivered non-sterile.

### PRECAUTIONS

- Do not use the product if the package has been violated.
- This product should be used in its sterile condition.
- Better results are achieved with the utilization of the Neodent products sequence. The utilization of instrumentations and/or prosthetic components from other manufacturers does not assure the perfect function of the Neodent Implant System and voids any product guarantee.
- The improper planning may endanger the Implant/prosthesis set performance resulting in faults in the system, such as Implant loss or fracture, loosening or fracturing of the Prosthetic Screws.

### ADVERSE EFFECTS

Adverse effects will only occur if the instrumentation selection and use are improper.

### POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

### STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

### DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

### DATE OF EXPIRATION

Written on the label.

### LIFE CYCLE

This product should be disposed of when it loses its functionality.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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