# ST.330.121

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions.

Neodent titanium Prosthetic Components are used as intermediaries between the implant placed intraorally and the prosthesis, planned according to each case, respecting the prosthetic interfaces, aiming for aesthetics and function.

## DESCRIPTION

It is an intermediary prosthetic component to be installed between the Implant and the Prosthesis (crown), made of titanium alloy according to standard ASTM F136. They are provided in a form that is anti-rotational in relation to the prosthesis. It is a single piece with conical geometry and a platform to seating of the CM Abutment Cylinder on the upper portion. They are available according to the table below:

Description / Prosthetic Interface	Transmucosal height (mm)
CM Abutment	0.8, 1.5, 2.5, 3.5, 4.5, 5.5, 6.5
WS Abutment	0.8, 1.5, 2.5, 3.5

## INDICATIONS FOR USE

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading.

## APPLICATIONS

It is indicated for single-tooth prostheses screw-mounted onto implants, according to the available interocclusal space, transmucosal height, and three-dimensional position of the implant. It can be used in immediate or conventional rehabilitation procedures on the maxilla or mandible.

## CONTRAINDICATIONS

This product is not indicated for customization. This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: Titanium alloy Ti6Al4V. Contraindicated for multiple and/or cemented prostheses, for insufficient interocclusal space, and unsatisfactory three-dimensional position of

the implant.

## HANDLING

For use of the CM Abutment/WS Abutment in two-stage procedures, prior preparation of the soft tissues can be done with use of a healing cap. Select the CM Abutment/WS Abutment according to prosthetic planning (it is recommended to use the Kit of Prosthetic Selection in the mouth) and install it onto the implant, according to the recommended torque of 32 N.cm and Torque Connection 1.6. Ensure the proper fit and install the CM Abutment Protective Cylinder.

For the impression procedure, fit the CM Abutment Transfer onto the CM Abutment/WS Abutment, ensure the proper fit, and take the impression using appropriate materials. Prepare the plaster model. Prepare the prosthesis, using the Cylinders of the CM Abutment according to appropriate laboratory techniques. Tests should be carried out on passivity and the fit of the structure of the prosthesis.

For installation of the final prosthesis, remove the CM Abutment Protective Cylinder and install it with torque of 20 N.cm on the Prosthetic Screw with the 1.2 Torque Connection. Protect the access for the screw with teflon and resin compound).

## TRACEABILITY LABEL

This product is accompanied by three labels that allow its traceability and should be attached to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

## PRESENTATION AND STERILIZATION

This product is indicated for single use and provided sterilized by the ethylene oxide method, packed individually into a package.

## PRECAUTIONS

- It is clinically recommended to maintain a minimum height of 4 mm of the component assembly and cylinder, as well as the transmucosal height, in creation of the prosthesis.
- In the case of the Morse Taper interface, it is recommended to use the CM Prosthetic Selection Kit. In the case of the WS interface, use the CM prosthetic selection kit, with trunnions compatible with the WS abutment/implant interface. CM Prosthetic Components are not compatible with the WS abutment/implant.
- Surgical planning and/or inadequate prosthesis can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening

or fracture of components and/or prosthetic screws.

- The selection of the material of the structure of the prosthesis must take local aspects of the patient into account.
- This product is of single use and cannot be re-sterilized.
- Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if the packaging is broken.
- The sterilization is guaranteed only if the packaging is not damaged.
- This product must be used immediately after opening the packaging, at the moment of the procedure. If it is not used, discard it.
- Do not use the product after expiry date.
- For immediate load application, check the torque indication of the implant installed.
- See the torque to be given on the prosthetic component to be used. The excess or lack of torque can lead to undesirable result.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- During installation, make sure to align it with the insertion axis of the implant. Make sure it is perfectly seated on the Implant. To do so, it is recommended to use periapical X-rays with the parallel technique.
- Make sure to use a prosthetic screw that is compatible with the prosthetic interface and with the prosthetic component.
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/prosthesis assembly.
- Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the Neodent product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the professional's responsibility to use the Neodent products according to the instructions for use.

**ADVERSE EFFECTS**

No adverse effects will occur, provided that the indications are complied with.

**MAGNETIC RESONANCE IMAGING (MRI)  
- SAFETY INFORMATION**

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

**POST-OPERATIVE PRECAUTIONS AND  
MAINTENANCE**

Instruct the patient as to the need of a professional medical monitoring after the procedure and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

**STORAGE CONDITIONS**

This product should be stored, in its original packaging, in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

**DISPOSAL OF MATERIAL**

After use, all the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

**DATE OF EXPIRATION**

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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