# ST.330.124

Português  Pilar Transeptelial

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions.

## DESCRIPTION

The Transepithelial Abutment is a prosthetic component manufactured in Titanium alloy in accordance to standard ASTM F136. It is available in the External Hex or Morse Taper interface, in the heights of 3, 4, 5 and 6 mm, in the rotational form.

## INDICATIONS FOR USE

Zygomatic Implants are indicated for surgical installation in the zygoma region, in cases of severe jaw resorption, in order to restore patient aesthetics and chewing function. Zygomatic Implants are recommended for the posterior (pre-molar/molar) region, one implant on each side, with at least two standard dental implants in the anterior region to support a fixed restoration. Zygomatic Implants may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading.

## APPLICATIONS

This product is indicated depending on the interocclusal space available, existing transmucosal height, and three-dimensional position of the Implant. It is indicated for multiple screw-mounted prostheses over Implants installed in the maxilla. It can be used in immediate or delayed rehabilitation procedures. The External Hex and Morse Taper models of the Transepithelial Abutment are used exclusively with Neodent Zygomatic Implants.

## CONTRAINDICATIONS

This product is not indicated for customization. This product is not recommended for rehabilitation on WS implants.

This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: Titanium alloy Ti6Al4V.

Contraindicated for single-tooth, partial, and/or cemented prostheses, insufficient interocclusal space, and unsatisfactory three-dimensional positioning of the implant.

## HANDLING

For use of the Transepithelial Abutment in two-stage procedures, prior preparation of the soft tissues with use of a healing cap is recommended. Select the Transepithelial Abutment according to the prosthetic planning and install it on the implant using the Connection for Prosthetic Abutments, according to the recommended torque (table below). Ensure the proper fit and install the Protective Cylinder of the Transepithelial Abutment

For the impression procedure, fit the Transepithelial Abutment Transfer onto the Transepithelial Abutment, ensure the proper fit, and take the impression using appropriate materials. Prepare the plaster model. Prepare the prosthesis, using the Cylinders of the Transepithelial Abutment according to appropriate laboratory techniques. Tests should be carried out on passivity and the fit of the

structure of the prostheses.

For installation of the definitive prosthesis, remove the Protective Cylinder of the Transepithelial Abutment and install it with torque of 10 N.cm on the Prosthetic Screw with the 1.2 Torque Connection. Protect the access for the screw with teflon and resin compound).

Description	Recommended torque (N.cm)
External Hex Transepithelial Abutment	20
Morse Taper Transepithelial Abutment	32

## TRACEABILITY LABEL

This product is accompanied by three labels that allow its traceability and should be attached to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

## PRESENTATION AND STERILIZATION

This product is indicated for single use and provided sterilized by the ethylene oxide method, packed individually into a package.

## PRECAUTIONS

• CM Prosthetic Components are not compatible with the CM Zygomatic prosthetic interface. Para seleção dos componentes protéticos para o Implante Zigomático CM indica-se o Kit Seleção Protética CM com munhões compatíveis com a interface protética WS. • Prosthetic components for WS Implants are not recommended for rehabilitations on Zygomatic Implants. • Surgical planning and/or inadequate prosthesis can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of components and/or prosthetic screws. • The selection of the material of the structure of the prosthesis must take local aspects of the patient into account. • This product is of single use and cannot be re-sterilized. • Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products. • Do not use the product if the packaging is broken. • The sterilization is guaranteed only if the packaging is not damaged. • This product must be used immediately after opening the packaging, at the moment of the procedure. If it is not used, discard it. • Do not use the product after expiry date. • For immediate load application, check the torque indication of the implant installed. • See the torque to be given on the prosthetic component to be used. The excess or lack of torque can lead to undesirable result. • Before each procedure, make sure the pieces are properly seated.

• Ensure that the parts are not swallowed or aspirated by the patient. • During installation, make sure to align it with the insertion axis of the implant. Make sure it is perfectly seated on the Implant. To do so, it is recommended to use periapical X-rays with the parallel technique. • Make sure to use a prosthetic screw that is compatible with the prosthetic interface and with the prosthetic component. • Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/prosthesis assembly. • Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear. • Always use the Neodent product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty. • It is the professional's responsibility to use the Neodent products according to the instructions for use.

## ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

## MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

## POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the procedure and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

## STORAGE CONDITIONS

This product should be stored, in its original packaging, in a clean and dry location, in a maximum temperature of 40° C and protected from direct sunlight.

## DISPOSAL OF MATERIAL

After use, all the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

## DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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