# ST.330.142

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

**INDICATIONS FOR USE**

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading. The Implant Titamax WS is indicated for application of load after osseointegration.

**DESCRIPTION**

The Neodent Titamax WS Medular Implant is made of commercially pure titanium (Grade 4).  
**MACROGEOMETRY:** •Cone Morse prosthetic interface exclusive to System WS; • Unique prosthetic interface, regardless of implant diameter; • Cylindrical format with pyramid profile threads; • Cervical diameter equal to diameter of implant body; • Comes with assembler.  
**NEOPOROS MICROSURFACE:** rough surface obtained by double treatment: abrasive blasting and acid subtraction.

Diameter of the Implant (mm)	Lengths (mm)
5 e 6	5 e 6

**APPLICATIONS**

The Titamax WS Medular Implant is recommended for surgical intraoral installation in bone with density of III or IV, according to Lekholm & Zarb’s jawbone quality classification (1985), in posterior regions of the maxilla or mandible, with bone height less than or equal to 7 mm. It can be used as a support for single or multiple prostheses in conventional loading protocol.

**CONTRAINDICATIONS**

They are contraindicated for rehabilitation by the immediate load technique.  
 The installation of this implant is not recommended in situations which require angulation of the prosthetic component, as there are no corresponding prosthetic components for the correction of possible undesirable angulation. This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: titanium.  
 In the presence of acute inflammatory or infectious processes in live tissue, unsuitable bone volume or quality, serious medical problems such as; bone metabolism disorders, blood clotting disorders, unsuitable healing capacity, insufficient oral hygiene, incomplete jawbone growth, patient uncooperative and not motivated, abuse of drugs or alcohol, psychosis, prolonged functional disorders which resist any treatment with medications, xerostomia, weakened immunological system, diseases which require the use of steroids, endocrinological diseases, pregnancy.


**HANDLING**

Perform drilling using drill in a good cutting conditions and rotation between 200 and 300 rpm, under abundant irrigation. Select the sequence of drills according to the intended implant. The insertion depth of the drills must comply with planning of the final position of the implant. For positioning below the bone of the coronal portion of the implant, it is recommended to add 1 to 2 mm in length to the implant during the surgical instrumentation. The latest CM Pilot Drills, in accordance with the diameter of the implant, have three marks corresponding to the planned position of the implant relative to the bone crest: level, 1 or 2mm below the alveolar crest. Neodent drills should be used in accordance with their specific instructions. For more information, consult the Neodent catalog.

	Initial Drill	Twist Drill 2.0	Pilot Drill 2/3	Twist Drill 3.0	Twist Drill 3.3	CM Pilot Drill 3 / 3.75	CM Pilot Drill 3.3 / 4.0	Twist Drill 3.8	CM Pilot Drill 3.8 / 4.3	Twist Drill 4.3	CM Pilot Drill 4.3 / 5.0	CM Pilot Drill 4.3 / 5.3	Twist Drill 5.3	CM Pilot Drill 5.3 / 6.0
5.0 mm	✓	✓	✓	✓		✓		✓	✓	✓	*			
6.0 mm	✓	✓	✓	✓		✓		✓	✓	✓		✓	✓	*

✓ Sequence recommended  
 \* Optional

**INSTALLATION OF THE DENTAL IMPLANT**

**WARNING:** The figures are merely illustrative. They do not represent the actual dimensions and features of the product.

- The cardboard and blister packagings should be opened, manually, without the use of sterile gloves.
- Break the seal of the cardboard packaging and remove the blister.
- Remove the protection from the blister pack.
- Deposit the sterile flask over the surgical field. NOTE: The clear tube and implant should be handled with a sterile surgical glove, in a surgical environment.
- Hold the bottle using the non-dominant hand and take the lid off.
- Take with the titanium pincers the assembler/implant unit, gripping only the assembler.
- Fit the assembler/implant unit in the hex connection for contra-angle, seeking the fit between the connection and the assembler.
- Transport the assembler/implant unit to the surgical alveolus.
- In the surgical motor, use maximum torque of 35 N.cm and rotation of 30rpm.
- Conclude the installation of the implant with the torquemeter ratchet with the hex connection coupled.
- With the digital wrench 1.2, remove the screw of the assembler and disconnect the assembler of the implant.
- The maximum installation torque suggested is 60 N.cm.

**PROSTHETIC SEQUENCE**

In the immediate postoperative period, install the WS coverage screw. During the selection of the prosthetic component, use the CM prosthetic selection kit, with sleeves compatible with the WS prosthetic interface. \*Note: The CM Prosthetic Components are not compatible with the WS prosthetic interface.

## TRACEABILITY LABEL

This product is accompanied by three labels that allow its traceability and should be attached to the following documents:

- Medical record;
- Collection tax document;
- Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

## PRESENTATION AND STERILIZATION

This product is indicated for single use and is supplied sterilized by gamma radiation method, uniformly packed in packaging that offers triple protection: cardboard, blister-type and clear tube.

## PRE-OPERATORY AND PROSTHETIC PLANNING

The model, diameter, length, position and quantity of implants must be selected for each clinical case, taking into consideration the anatomy, quality and quantity of bone and space available. When necessary execute the diagnostic wax-up of the patient. In situations in which there are relatively high loads, special care must be taken to ensure the suitable alignment of the implant(s), component(s) prosthetics(s) and prosthesis.

## PRECAUTIONS

- CM Prosthetic Components are not compatible with the WS prosthetic interface.
- For short implants, clinicians should closely monitor patients for any of the following conditions: periimplant bone loss, changes to implant's response to percussion, or radiographic changes in bone to implant contact along the implant's length. If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal. If the clinicians choose a short implant, then clinicians should consider a two-stage surgical approach, splinting a short implant to an additional implant, and placement of the widest possible fixture. Allow longer periods for osseointegration and avoid immediate loading.
- Note the conditions of the intra-oral tissue, the bone quality and quantity of the bed receiving the Implant, by means of radiographic and/or tomography examinations. Nonperformance of the pre-surgical assessment may compromise the success of the procedure.
- As for the systemic aspect, consider the general health of the patient. In particular, one must be careful in cases of patients who have allergies to drugs, local or systemic factors that may interfere with the healing process of the bone tissues or soft tissues, or the process of osseointegration. For example, the bone already exposed to radiation in the head and neck area, diabetes mellitus, anticoagulation drugs, hemorrhagic diathesis, bruxism, parafunctional habits,

anatomically unfavorable bone situation, tobacco abuse, uncontrolled periodontitis, treatable jaws pathologies and oral mucosa abnormalities.

- Treatment with bisphosphonates results in potential risk of peri-implant osteonecrosis.
- Surgical planning and/or inadequate prosthesis can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of components and/or prosthetic screws.
- The material to be used during the procedure must be sterile.
- This product is of single use and cannot be re-sterilized.
- Reuse of this product may cause:
  - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
  - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if the packaging is broken.
- The sterilization is guaranteed only if the blister packaging is not damaged.
- This product must be used immediately after opening the packaging, at the moment of surgery. If it is not used, discard it.
- Do not use the product with the validity expired.
- The maximum installation torque suggested is 60 N.cm. Insertion torque greater than that recommended can make the system inoperative.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/prosthesis assembly.
- Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the Neodent product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the professional's responsibility to use the Neodent products according to the instructions for use.

## ADVERSE EFFECTS

The installation of dental implants, as well as any other surgical procedure, may cause a slight discomfort and localized edema. More persistent

symptoms can occur such as: chronic pain related to the dental implant, permanent paresthesia, dysesthesia, loss of maxillary/mandibular bone edge, systemic of localized infection, oroantral or oronasal fistula, adjacent teeth affected unfavorably, irreversible damage to adjacent teeth, fracturing of the implant, jaw, bone or prosthesis, esthetic problems, injury of the nerves, exfoliation, hyperplasia.

Any failure in osseointegration and loss of the prosthesis during treatment may be caused by: inadequate osteotomy, infections, diseases or systemic problems, low quality or volume of the remaining bone, absence or failure of irrigation, use of nonspecific instruments and/or instruments without cutting power, poor oral hygiene, occlusal trauma, lack of prosthetic passivity and lack of specific training.

## MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

## POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

## STORAGE CONDITIONS

This product should be stored, in its original packaging, in a clean and dry location, in a maximum temperature of 40 °C and protected from direct sunlight.

## DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

## DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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