ST.330.153


Português  Implante Zigomático

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English  Zygomatic Implant

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Simbologia Symbology Simbologia Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
Qty	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
PROIBIDO REPROCESSAR	Requisito Brasileiro / Brazilian Requirement / Requisito Brasileño / Requisito Brasiliano / Brasilianische Anforderungen / Exigences brésiliennes
	Produto esterilizado por radiação gama / Product sterilized through gamma rays / Producto esterilizado por radiación gama / Prodotto sterilizzato attraverso raggi gamma / Sterilisiert durch Gammastrahlung / Produit stérilisé par rayonnement gamma

Simbologia Symbology Simbologia Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Não reutilizar / Do not reuse / No lo reutilice / Non riutilizzare / Nicht wiederverwendbar / Ne pas réutiliser / Nicht wiederverwendbar / Ne pas réutiliser
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Conservelo seco / Conservare in luogo asciutto / Trocken lagern / Garder au sec
	Manter afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conservar à l'abri de la lumière solaire
	Não utilizar se a embalagem estiver danificada / Do not use if package is damaged / No lo utilice si el envoltorio está dañado / Non utilizzare se la confezione originale è danneggiata / Bei beschädigter Verpackung Produkt nicht verwenden / Ne pas utiliser si l'emballage a été abîmé

Simbologia Symbology Simbologia Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Marcação CE para comercialização na Comunidade Europeia / CE Mark for European Community market / Marca CE para comercialización en la Comunidad Europea / Marchio CE per la commercializzazione in Europa / CE Kennzeichen zum Vertrieb in der Europäischen Gemeinschaft / Marquage CE pour une commercialisation en Europe
Rx only	Notificação exigida pelo FDA para comercialização nos Estados Unidos / Notification required by FDA for United States market / Notificación exigida por el FDA para comercialización en los Estados Unidos / Notifica necessaria per la commercializzazione negli Stati Uniti / Vom FDA geforderte Kennzeichnung für den Vertrieb in den USA / Notification exigée par FDA pour une commercialisation aux États-Unis
	Representante na Comunidade Europeia / Representative in the European Community / Representante en la Comunidad Europea / Rappresentante nella Comunità europea / Vertreter in der Europäischen Gemeinschaft / Représentant dans la Communauté Européenne
	Fabricante / Manufacturer / Fabricante / Produttore / Hersteller / Fabricant

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants, with knowledge on the technique for installing Zygomatic implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere. **IMPORTANT:** This technique requires specific accreditation. For more information, visit the site at www.neodent.com.br.

DESCRIPTION

The Neodent Zygomatic Implant is made of commercially pure Titanium (Grade 4).
 • Cylindrical shape, pyramidal threads, smooth surface, and self-threading. • It has angling of 45° from the prosthetic interface, for compensation of the angle between the zygoma and maxilla. The diameter of the portion anchored in the zygoma is 3.9mm, and that of the portion anchored in the residual alveolar process of the maxilla is 4.4mm. • An External Hex and Morse Cone is available on the prosthetic interfaces (exclusively of the Zygomatic System), in the lengths of 30mm, 35mm, 40mm, 45mm, 47.5, 50mm, or 52.5mm. • It comes with a mount (pre-mounted).

INDICATIONS FOR USE

Zygomatic Implants are indicated for surgical installation in the zygoma region, in cases of severe jaw resorption, in order to restore patient esthetics and chewing function. Zygomatic Implants are recommended for the posterior (pre-molar/molar) region, one implant on each side, with at least two standard dental implants in the anterior region to support a fixed restoration. Zygomatic Implants may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading.

APPLICATIONS

The Zygomatic Implant is recommended for surgical installation and must be inserted in the posterior region of the maxilla and the zygoma. It is recommended in cases of severe reabsorption of the maxilla, and in total tooth loss, when the pneumatization of the maxillary sinus and reduced bone height reduce the posterior region of the maxilla (situations in which the insertion of the conventional implant and reconstructive procedures, employing bone grafts, are contraindicated), and must be used in association with two, three, or four conventional implants. Note: For application of immediate load and primary stability, it must reach at least 32N.cm.

CONTRAINDICATIONS

When there is symptomatology related to the facial sinuses (maxillary sinusitis, polyps or any other sinus pathologies); Cases of unfavorable maxillary relations. This product is not recommended for cases of unitary rehabilitations. The installation of this product is not recommended before the end of the growth and epiphyseal occlusion. Not recommended in the presence of acute inflammatory or infectious processes in live tissue, poor bone volume or quality, remains of roots in the place, serious medical problems such as; bone metabolism disorders, blood clotting disorders, unsuitable healing capacity, insufficient oral hygiene, incomplete jawbone growth, patient uncooperative and not motivated, abuse of drugs or alcohol, psychosis, prolonged functional disorders which resist any treatment with medications, weakened immunological system,

diseases which require the use of steroids, allergy or hypersensitivity to titanium, endocrinological diseases, allergy to drugs, bone already exposed to radiation, diabetes mellitus, anticoagulation medications/ hemorrhagic diathesis, teeth grinding, parafunctional habits, anatomically unfavorable bone situation, tobacco abuse, uncontrolled periodontitis, treatable pathological maxillary diseases, alterations of the oral mucous and pregnancy.

HANDLING

After the 45° palatal incision on the soft tissue along the entire maxillary crest, the soft tissue is fully folded back from the maxillary crest up to the Zygomatic bone post, and the infra-orbital nerve is located. A window is made by drilling at the upper boundary between the zygoma and the maxillary sinus. This window will be useful during the procedure to cool the perforations, in order to prevent overheating. Preparation of the recipient bed: it should be performed with the necessary precautions to reduce the surgical trauma. Perform drilling using drill in a good cutting conditions and rotation between 800 and 1200 rpm. Drilling must be performed along with substantial irrigation. Drills with progressive diameters are used, finishing with the insertion of the self-threading Zygomatic implant at low speed. The length of the implant is carefully chosen by means of a specific probe. Neodent drills should be used according to their instructions. During handling of the drills for Zygomatic, it is suggested to use the Zygomatic Lip Protector. For more information, see the Neodent catalog, available on the website www.neodent.com.br.

INSTALLATION OF THE DENTAL IMPLANT

• The cardboard and blister packagings should be opened without the use of sterile gloves.
 • Break the seal of the cardboard packaging and remove the blister. • Break the blister's sterility seal. • Place the bottle containing the implant on a sterile surface. NOTE: Handling the bottle and the implant must be performed under sterile conditions. • Hold the bottle using the non-dominant hand and take the lid off. • Hold the implant with the zygomatic connection • contra-angle, keeping the connection stable and slightly rotating bottle, searching for the perfect fit between the connection and the implant. • Take the implant to the surgical cavity. • For installation using a motor, use the zygomatic connection • contra-angle. Use maximum torque 35N.cm and rotation of 30rpm in the surgical motor • For manual installation, use the zygomatic installation driver. • Insert the whole body of the implant and assess, at the end of the installation, the positioning of the emergence of the screw of the mount in relation to the occlusal edge. The screw must be in a position that is favorable for the installation of the planned prosthetic component. • The maximum installation torque suggested is 60N.cm. • With Manual Screwdriver 1.2, remove the screw from the mount and disconnect the mount from the implant and discard it. The choice between installing the cover screw, the healing abutment or the prosthetic component is at the professional's discretion.

PROSTHETIC SEQUENCE

Select Transepithelial Abutment, checking it's indication and limitations. CM Prosthetic Component are not compatible with the

Zygomatic CM prosthetic interface. For selection of the prosthetic components for the CM Zygomatic Implant, the CM Prosthetic Selection Kit, with abutments compatible with the WS prosthetic interface, is recommended.

TRACEABILITY LABEL

This product is accompanied by three labels that allow its traceability and should be attached to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

PRESENTATION AND STERILIZATION

This product is indicated for single use and is supplied sterilized by gamma radiation method, uniformly packed in packaging that offers triple protection: tertiary packaging (cardboard), secondary packaging (blister-type) and primary packaging (clear tube).

PRE-OPERATORY AND PROSTHETIC PLANNING

The model, diameter, length, position and quantity of implants must be selected individually, taking into consideration the anatomy, quality and quantity of bone and space available. When necessary execute the diagnostic wax-up of the patient. In situations in which there are relatively high loads, special care must be taken to ensure the suitable alignment of the implant /component and prosthesis.

PRECAUTIONS

• CM Prosthetic Components are not compatible with the Zygomatic CM prosthetic interface.
 • The Zygomatic Implants withstand functional load if they are rigidly connected to a minimum of two or more standard Implants. • In cases of severe reabsorption of the anterior region of the maxilla, in edentulous patients, associated with the use of bone grafting, the use of zygomatic implants reduces the dimensions of the graft.
 • Surgery for installation of Zygomatic Implants is not recommended without local and systematic evaluation of the patient. Besides that, it is important to have the assessment of a medical professional to carry out the surgery. • Three levels need to be investigated during the pre-operative planning: the maxillary level, the sinus, and the zygoma. The assessment of the computerized tomography allows for visualization of the health of the maxilla and the maxillary sinus. The density, length, and volume of the zygoma can be assessed and special guides for the insertion of the zygomatic implants can be prepared on prototype models to facilitate the orientation of the zygomatic implants during the surgery, with minimal errors of angling and positioning. • The surgical procedure of the Zygomatic Implant must involve a minimally traumatic surgery, preventing bone overheating in the zygoma and maxilla. • Although the surgery can be done under local anesthesia, for the comfort and safety of the patient, it is recommended that it be done with general anesthesia or neuroleptic analgesia, in a sterile environment. • Note the conditions of the intra-oral tissue, the bone quality and quantity of the bed receiving the Implant, by means of radiographic and/or tomography examinations. Failure to perform the pre-surgical assessment may compromise the success of the procedure.
 • Regarding the systemic aspect consider the general state of health of the patient, who

must undergo a careful clinical and radiological examination before surgery, assessing the psychological and physical status. Special attention should be paid in cases of patients who have local or systemic factors which may interfere with the healing processes of the bone or soft tissue, or in the process of bone integration (refer to CONTRAINDICATIONS). The Non-execution of the pre-surgical evaluation can lead to it being impossible to perceive pre-existing diseases.

- Inappropriate planning can jeopardize the performance of the implant/prosthesis unit leading to failures of the system, as loss or fracturing of the implant, loosening or fracturing of the prosthetic screws.
- The handling as well as all the material to be used during the procedure must be sterile.
- This product is of single use and cannot be re-sterilized.
- Reuse of this product may cause:
 - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
 - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if the packaging is broken.
- The sterilization is guaranteed only if the blister packaging is not damaged.
- This product must be used immediately after opening the packaging, at the moment of surgery. If it is not used, discard it.
- Do not use the product with the validity expired.
- For application of immediate load the primary stability must attain at least 32N.cm.
- The maximum installation torque suggested is 60 N.cm.
- See the torque to be given on the prosthetic component to be used. The excess or lack of torque can lead to undesirable result.
- Adjust the surgical motor used in the procedure (torque and RPM) as per specification of the implant to be used.
- Insertion torque greater than that recommended can make the system inoperative.
- The maximum angulation allowed for Neodent implants is up to 30°

- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Make the occlusal and interproximal adjustment after the installation of the prosth.
- The lack of adjustment or inappropriate occlusal adjustment can endanger the implant/prosthesis unit.
- Some of the possible causes of unsuccessful bone integration and loss of the prosthesis during the treatment are the following: unsuitable osteotomy, infections, deficient oral hygiene, occlusal trauma, systemic problems or diseases, low remaining bone quantity or quality, lack or failure of irrigation, use of instruments which are not specific and/or without cutting power and lack of specific training.
- Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the Neodent product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the professional's responsibility to use the Neodent products according to the instructions for use, and to determine whether it suits the individual situation of each patient.

ADVERSE EFFECTS

The installation of Zygomatic Implants is an invasive surgical procedure with associated risks, and can provoke discomfort and edema. The installation of dental implants, as well as any other surgical procedure, may cause a slight discomfort and localized edema. More persistent symptoms can occur such as: chronic pain related to the dental implant, permanent paresthesia, dysesthesia, loss of maxillary/mandibular bone edge, systemic or localized infection, oroantral or oronasal fistula, adjacent teeth affected unfavorably, irreversible damage to adjacent teeth, fracturing of the implant, jaw, bone or prosthesis, esthetic problems, injury of the nerves, exfoliation, hyperplasia.

MAGNETIC RESONANCE IMAGING (MRI)

• SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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