ST.330.176

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English  Universal Abutment

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions.

Neodent titanium Prosthetic Components are used as intermediaries between the implant placed intraorally and the prosthesis, planned according to each case, respecting the prosthetic interfaces, aiming for aesthetics and function.

DESCRIPTION

It is an intermediary prosthetic component to be installed between the Implant and the prosthesis (crown), made of titanium alloy according to standard ASTM F136. They are provided in a form that is anti-rotational in relation to the prosthesis, according to the table below:

| Interface | Diameter (mm) | Cementable Area (mm) | Trans-mucosal (mm) | |
|---------------------|---------------|----------------------|--|------------------------|
| CM | 3.3 and 4.5 | 4 and 6 | 0.8, 1.5, 2.5, 3.5, 4.5, 5.5, 6.5 | |
| CM Angled 17° | | | 1.5, 2.5, and 3.5 | |
| CM Angled 30° | | | 1.5, 2.5, and 3.5 | |
| CM thru screw | | | 0.8, 1.5, 2.5, 3.5, 3.5, 4.5, 5.5, 6.5 | |
| CM Exact | | | 0.8, 1.5, 2.5, 3.5, 4.5, 5.5 and 6.5 | |
| CM Exact Angled 17° | | | 1.5, 2.5, and 3.5 | |
| CM Exact Angled 30° | | | 1.5, 2.5, and 3.5 | |
| II Plus | | | 1, 2, and 3 | |
| WS | | | 4.5 | 0.8, 1.5, 2.5, and 3.5 |
| Munhão Universal SF | | | | |

INDICATIONS FOR USE

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading.

APPLICATIONS

It is indicated for single-tooth prostheses cemented onto implants, according to the available interocclusal space, transmucosal height, and three-dimensional position of the implant. It can be used in immediate or conventional rehabilitation procedures on the maxilla or mandible.

CONTRAINDICATIONS

The 30° CM Angled Universal Abutment in combination with implants of 3.5 mm diameter are contraindicated for single-tooth final rehabilitations in the molar region.

This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: Titanium alloy Ti6Al4V.

Contraindicated for insufficient interocclusal space and unsatisfactory three-dimensional position of the implant.

HANDLING

For use of the Universal Abutment in two-stage procedures, prior preparation of the soft tissues can be done with use of a healing cap.

Select the Universal Abutment according to the prosthetic planning (it is recommended to use the Kit of Prosthetic Selection in the mouth) and install it on the implant using the recommended torque and connection (table below). Make sure of the correct fit between the Universal Trunnion and the implant. Protect the access for the screw (with teflon and resin compound), if applicable. It is recommended to make up a provisional prosthesis (using the corresponding acrylic Cylinder of the Universal Trunnion) for gingival conditioning.

For the impression procedure, fit the Transfer of the corresponding Universal Trunnion, make sure of the fit with a “click” and do the impression with the closed technique using appropriate materials. Remove the mold with the Transfer of the Universal Trunnion and fit the Analog of the Universal Trunnion on the Transfer of the Universal Trunnion. Prepare the plaster model. Use the Universal Abutment Coping for creation of the structure of the prosthesis. Tests should be carried out on passivity and the fit of the structure of the prosthesis. Cement the final prosthesis on the Universal Abutment. Cementation is done in the mouth following the instructions from the cement manufacturer. In this case, it is worth remembering that it is necessary to take note of the technique to prevent excess cement and chemical burns on the tissues adjacent to the Implant. Protect the access of the screw (with teflon and resin compound) during the cementing process, if applicable.

| Component | Recommended Torque (N.cm) | Hex Screwdriver Torque Connection (mm) |
|----------------------------------|---------------------------|--|
| CM Exact Universal Abutment | 15 | 0.9 |
| CM Angled Universal Abutment | | |
| CM Thru-Screw Universal Abutment | | |
| Universal Abutment II Plus | 20 | 1.2 |
| CM and WS Universal Abutment | 32 | 1.2 |
| SF Universal Abutment | | |

TRACEABILITY LABEL

This product is accompanied by three labels that allow its traceability and should be attached to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

PRESENTATION AND STERILIZATION

This product is indicated for single use and provided sterilized by the ethylene oxide method, packed individually into a package.

PRECAUTIONS

- In the case of the Morse Taper interface, it is recommended to use the CM Prosthetic Selection Kit. In the case of the WS interface, use the CM prosthetic selection kit, with trunnions compatible with the WS abutment/implant interface. CM Prosthetic Components are not compatible with the WS abutment/implant.
- Surgical planning and/or inadequate prosthesis can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of components and/or prosthetic screws.
- The selection of the material of the structure of the prosthesis must take local aspects of the patient into account.
- This product is of single use and cannot be re-sterilized.
- Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if the packaging is broken.
- The sterilization is guaranteed only if the packaging is not damaged.
- This product must be used immediately after opening the packaging, at the moment of the procedure. If it is not used, discard it.
- Do not use the product after expiry date.
- For immediate load application, check the torque indication of the implant installed.
- See the torque to be given on the prosthetic component to be used. The excess or lack of torque can lead to undesirable result.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- During installation, make sure to align it with the insertion axis of the implant. Make sure it is perfectly seated on the Implant. To do so, it is recommended to use periapical X-rays with the parallel technique.
- Make sure to use a prosthetic screw that is compatible with the prosthetic interface and with the prosthetic component.
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/prosthesis assembly.
- Before each procedure check the conditions

of the **Neodent** surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.

- Always use the **Neodent** product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the **Neodent** Implant System and exempts any product warranty.
- It is the professional's responsibility to use the **Neodent** products according to the instructions for use.

ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The **Neodent** Implant System has not been

evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the **Neodent** Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the procedure and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored, in its original packaging, in a clean and dry location, in a

maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

After use, all the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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