ST.330.184

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
English  Sinus Lift Curette

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

This product is an instrument in the shape of a spoon or spade, with non-cutting edges on its ends.

APPLICATIONS

This product is an instrument used in dental surgery, indicated for maxillary sinus lifting.

CONTRAINDICATIONS

This product has no contraindications, providing their intended purposes are correctly followed.

HANDLING

This product is a digital grasping instrument, used according to its angulation, to separate the sinus mucosa from the bone that constitutes the maxillary sinus. Thus, an adequate space is created for local bone grafting.

SANITATION

The surgical instrumentation should be properly sanitized after its utilization. Do not let the product contact moisture beyond the time required for the sanitation. For such, proceed as follows: Step 1: Submerge the parts fully into enzymatic detergent (diluted according to the manufacturer). Step 2: Wash in ultrasonic washer for approximately 10 to 15 minutes. Step 3: Rinse with plenty of distilled water until the solution wastes are fully removed. The utilization of nylon brushes is recommended. Step 4: Dry with a clean and dry cloth or with compressed air. Step 5: Perform a visual inspection, watching if there are faults in the cleaning process. If any wastes are still left, the part should be submerged again in detergent - step 1 - and if necessary, the cleaning should be made with the help of a nylon brush.

Repeat the rinsing and drying sequence. Step 6: If applicable, place the product in a suitable case duly sanitized and dried. Step 7 - Select the package according to the sterilization process. Use preferably a self-sealing package if surgical grade paper with laminated film. ATTENTION: Do not use descaling products and avoid drying parts still containing cleaning solution wastes because such procedures favor oxidation. Also, the use of steel brushes is not recommended.

PRESENTATION AND STERILIZATION

This product is manufactured in stainless steel, sold non-sterile and can be reused in accordance with the guidelines and precautions contained in these usage instructions. This product is unitarily provided, conditioned in blister (surgical grade film and paper type) packing. It is reusable, delivered non-sterile and should, therefore, be sterilized before use. Sterilize the product on the day before or on the day of the procedure. We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard ISO 17665-1:2006 • Sterilization of health care products. Moist heat Requirements for the development, validation and routine control of a sterilization process for medical devices. Suitable parameters: 121 °C, 30 min, at 1 atm. Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Sterilization validity: 7 to 15 days, provided that the product is conditioned in a clean, dry environment away from the sunlight.

PRECAUTIONS

• Do not use the product if the package has been violated. • Better results are achieved with the utilization of the Neodent products sequence. The utilization of instrumentations and/or prosthetic components from other manufacturers does not assure the perfect function of the Neodent Implant System and voids any product guarantee. • Improper planning can jeopardize the performance of the Implant/Prosthesis set,

causing failures to the system, such as loss or break of the Implant, loosening or fracture of the Prosthetic Screws.

ADVERSE EFFECTS

This product is a surgical instrument used in surgical procedures. Adverse effects will only occur if the instruments are chosen incorrectly, which can cause damage to the surgical sequence of the implant to be installed.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40 °C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

LIFE CYCLE

This product should be disposed of when it loses its functionality.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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