# ST.330.189

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

## DESCRIPTION

**Coping:** prosthetic component made of castable polymer material for laboratory use in the phases of waxing and casting of infrastructure of prosthesis screwed or cemented on dental implants. It features external retentions. **UCLA:** prosthetic component made of castable polymer material for laboratory use in the phases of waxing and casting of infrastructure of prostheses screwed or cemented on dental implants. Available in rotational and anti-rotational versions.

## APPLICATIONS

**Cylinder:** Suitable for manufacturing infrastructures of multiple or single prostheses (table) screwed or cemented on the Prosthetic Component.

Prosthetic Component	Indication
CM Abutment	Single Prostheses
SF Conical Abutment	
Conical Abutment II Plus	
CM Mini Conical Abutment	Multiple Prostheses
CM Micro Conical Abutment	
SF Mini Conical Abutment	
Angled Mini Conical Abutment	
Mini Conical Abutment II Plus	
Angled Conical Abutment II Plus	Single or Multiple Prostheses

**UCLA:** Suitable for manufacturing infrastructures of multiple or single prostheses (table) screwed or cemented on the dental Implant.

## CONTRAINDICATIONS

This product does not show contraindications, provided that it be properly used for the shown purposes.

## HANDLING

After installation of Implant or intermediate in the mouth, transfer its position by means of molding to a working model, which will have its corresponding Analog. Place the Prosthetic Component on the analog and perform waxing. The obtained pattern shall be included in coating compatible with the alloy to be used in the casting process. Thus, the structure that will serve as the basis for acrylization or porcelain application is obtained. Recommended torque:

Prosthetic Component	Screwdriver	Recommended Torque
Coping	1.2	10 Ncm
UCLA External Hex		32 Ncm

The use of titanium Screw in the laboratory and Neotorque Screw for intraoral fixation is recommended.

## PRESENTATION AND STERILIZATION

This product is unitarily provided, conditioned in blister (surgical grade film and paper type) packing. It is delivered non-sterile and it of single use.

## PRECAUTIONS

- This product is of single use.
- Reuse of this product may cause:
  - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
  - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if the package has been violated.
- This product should not be used as a plastic temporary element.
- The metal choice should be done by a dentist in order to minimize the possible adverse events for each patient.
- An insertion torque higher than that recommended and improper Drivers may cause damages to the

material and make the system unusable.

- Better results are achieved with the utilization of the Neodent products sequence. The utilization of instrumentations and/or Prosthetic Components from other manufacturers does not assure the perfect function of the Neodent Implant System and voids any product guarantee.
- A special precaution should be taken during the closure so that the interocclusal space is not insufficient or excessive.
- Improper planning can compromise the performance of all the teeth or Implants/ Prosthesis resulting in system failures, such as the loss or breakage of a tooth/Implant, the loosening or fracture of a Prosthetic Screw.

## ADVERSE EFFECTS

Adverse effects will only occur if the choice and the use of instrument are improper.

## POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

## STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

## DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

## DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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