ST.330.190

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

Coping: Prosthetic Component made of titanium alloy suitable for manufacturing temporary prosthesis screwed to Abutments. It features external retention points for acrylic resin fixation. **UCLA:** Prosthetic Component made of titanium alloy suitable for manufacturing temporary prostheses cemented on dental Implants. It features a faceted side to promote prosthesis anti-rotational effect. Available in rotational and anti-rotational versions.

INDICATIONS FOR USE

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading.

APPLICATIONS

Coping: Suitable for manufacturing multiple or single temporary prosthesis (table). It can be used directly in the mouth or in the laboratory using plaster models and abutment analogs.

Prosthetic Component	Indication
CM Abutment	Single Prostheses
SF Conical Abutment	
Conical Abutment II Plus	
CM Mini Conical Abutment	Multiple Prostheses
CM Micro Conical Abutment	
SF Mini Conical Abutment	
Angled Mini Conical Abutment	
Mini Conical Abutment II Plus	
Angled Conical Abutment II Plus	Multiple Prostheses

UCLA: Suitable for manufacturing multiple or single temporary prosthesis screwed or cemented on dental Implants. It can be used directly in the mouth or in laboratory using plaster models and Implant Analogs.

CONTRAINDICATIONS

This product does not show contraindications, provided that it be properly used for the shown purposes.

HANDLING

This product can be used in two ways: 1) After installation of Implant or intermediate in the mouth, transfer its position by means of a molding to a working model, which will have

its corresponding Analog. Place the Prosthetic Component on the Analog and manufacture the temporary prosthesis. 2) The Prosthetic Component can be fixed directly in the mouth on the Implant or intermediate, and then manufacture the temporary prosthesis directly on it. In this case, care should be taken to avoid chemical burns to Implant adjacent tissues. An autopolymerizing opacifier may be used to mask the metal's gray appearance, preventing the temporary prosthesis from being gray due to translucence, which is very common when using acrylic resins. Recommended torque:

Prosthetic Component	Torque Connection	Recommended Torque
Coping	1.2	10 Ncm
UCLA External Hex		32 Ncm
UCLA Internal Hex		20 Ncm

The use of titanium Screw in the laboratory and Neotorque Screw for intraoral fixation is recommended.

TRACEABILITY LABEL

The product is identified through numerical codes arranged on the label (REF and LOT). The product is fitted with these codes, allowing their characteristics to be identified since when it is a raw material. This product is fitted with three labels which enable it to be tracked and should be posted to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor).

PRESENTATION AND STERILIZATION

This product is unitarily provided, conditioned in blister (surgical grade film and paper type) packing. It is delivered sterile and it of single use. Sterilization method: ethylene oxide.

PRECAUTIONS

• As to the systemic aspect, consider the patient's general health condition according to the applicable literature. As to the local aspect, please observe the intra-oral tissues conditions. • This product is of single use and it may not be re-sterilized. • Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products. • Do not use the product if the package has been violated. • An insertion torque higher than that recommended and improper Drivers may cause damages to the material and make the system unusable. • This product should be used immediately after opening the package, at the time of prosthesis's creation in a clinic or laboratory. If you do not use it,

discard it. • Better results are achieved with the utilization of the Neodent products sequence. The utilization of instrumentations and/or Prosthetic Components from other manufacturers does not assure the perfect function of the Neodent Implant System and voids any product guarantee.

• Make sure that the parts are not swallowed or inhaled by the person who is receiving the product (patient). • Improper planning can compromise the performance of all the teeth or Implants/Prosthesis resulting in system failures, such as the loss or breakage of a tooth/Implant, the loosening or fracture of a Prosthetic Screw. • During the installation, make sure to align it to the Implant insertion axle, avoiding jamming and damage to the threads. • Make sure that it is perfectly seated in the Implant. For that periapical radiographs are recommended with the parallelism technique.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and the use of instrument are improper.

MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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