
# ST.330.192

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

## DESCRIPTION

The Neodent Coping and UCLA Abutment have tapered shape, are manufactured in a Cobalt-Chromium alloy with high melting range and receive a cylindrical cover in castable polymer. They are used for definitive prosthesis confection. They are available in rotational (for bridge) and anti-rotational (for crown) versions. They are used for infrastructures overcasting of prosthesis in lab environment. Alloy characteristics in which the component is made • Melting range: 1287 to 1378 °C (2348,6 to 2512,4 °F) • Density: 8,28 g/cm<sup>3</sup> • Thermal expansion coefficient: 14,2 µm/°C (600 °C).

## INDICATIONS FOR USE

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading.

## APPLICATIONS

The Neodent cobalt-chromium Coping is used for waxing confection and overcasting of screw-retained prosthesis infrastructures installed over intermediaries. It is used in the laboratory phase for creating the prosthesis on plaster models and similar of the Implants. After casting, it becomes part of the permanent prosthesis. It is necessary to respect the diameters and tolerance indications of the components.

Prosthetic Component	Indication
CM Abutment	Single Protheses
SF Conical Abutment	
II Plus Conical Abutment	
CM Mini Abutment	Multiple Protheses
CM Micro Abutment	
SF Micro Conical Abutment	
Angled Mini Conical Abutment	
II Plus Mini Conical Abutment	
II Plus Angled Conical Abutment	

UCLA: Suitable for manufacturing multiple or single temporary prosthesis screwed or cemented on dental Implants. It is used in the laboratory phase for creating the prosthesis on plaster models and similar of the Implants. After casting, it becomes part of the permanent prosthesis.

## CONTRAINDICATIONS

This product does not show contraindications, provided that it be properly used for the shown purposes.

## HANDLING

The molding, the adaptation of the corresponding

intermediary Analogs and the modeling in plaster are made with transfers. The prosthetician must install this product over the corresponding intermediaries Analogs, join them with resin and perform the waxing. The obtained standard must be included in a compatible covering with the alloy that will be used in casting process. The chosen alloy must be an alloy with cobalt-chromium base and compatible with the esthetic material that is planned to the prosthesis. It is recommended to use phosphated covering to include the standard in wax. To obtain infrastructures for metaloceramic prostheses, must be chosen a compatible alloy, according to the norms ASTM F1537, ASTM F799, ISO 5832-12 and ISO 5832-4. The porcelain cannot be applied directly over the machined alloy of the prosthetic component, and, therefore, the waxing must ensure a recover of entire external area of the component over which the porcelain will be applied. The porcelain application (specific for alloys of this type) directly over an area not recovered by the metallic alloy of overcasting can cause crackings. Be careful during the waxing of the component margins to avoid invasion on the margins that will come into contact with the Implant. The waxing thickness must be, at least, of 0,5 mm, and can be reduced to 0,3 mm after overcasting. After casting, observe the corresponding area to the Prosthetic Screw seating. The original shape must be preserved, without the presence of metallic positive bubbles. During finishing and polishing procedures, the prosthetic interface of the component must be isolated by the respective Polishing Protector. During the process of alloy finishing, do not use corrosive materials, such as which contain iron particles. It is recommended to perform the alloy machine in low rotation and under cooling.

Alloy Composition CoCr	
Cobalt	62,5 a 68,65%
Chromium	26 a 30%
Carbon	0,2 a 0,3%
Molybdenium	0,2 a 0,3%
Nitrogen	0,15 a 0,2%

Recommended torque:

Prosthetic Component	Torque Connection	Recommended Torque
Coping		10 Ncm
UCLA External Hex	1.2mm	32 Ncm
UCLA Internal Hex		20 Ncm

The use of titanium Screw in the laboratory and Neotorque Screw for intraoral fixation is recommended.

## TRACEABILITY LABEL

The product is identified through numerical codes arranged on the label (REF and LOT). The product is fitted with these codes, allowing their characteristics to be identified since when it is a raw material. This product is fitted with three labels which enable it to be tracked and should be posted to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor).

## PRESENTATION AND STERILIZATION

The CoCr Coping is accompanied by a Neotorque Screw. The CoCr UCLA is supplied by unit. The Coping and UCLA are protected by a blister type packaging

(film and medical grade paper). They are delivered non-sterile and are for single use.

## PRECAUTIONS

CAUTION: For the overcasting, always use the cobalt-chromium alloy with melting range lower than showed above. Another metallic alloys do not ensure an adequated union. -The porcelain should never be applied directly over the adaptation belt of the prosthetic component. -Welding must be done at the range of 1120 to 1150 °C (2048 - 2102 °F). -This product is of single use and it may not be re-sterilized. • Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products. -Do not use the product if the package has been violated. -Better results are achieved with the utilization of the Neodent products sequence. The utilization of instrumentations and/or Prosthetic Components from other manufacturers does not assure the perfect function of the Neodent Implant System and voids any product guarantee. • Make sure that the parts are not swallowed or inhaled by the person who is receiving the product (patient). -Improper planning can compromise the performance of all the teeth or Implants/Prosthesis resulting in system failures, such as the loss or breakage of a tooth/Implant, the loosening or fracture of a Prosthetic Screw. • Make sure that it is perfectly seated in the Implant. For that periapical radiographs are recommended with the parallelism technique.

## ADVERSE EFFECTS

Adverse effects will only occur if the choice and the use of instrument are improper.

## POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

## STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40 °C and protected from direct sunlight.

## DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

## DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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