ST.330.201

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

Prosthetic Component suitable for manufacturing prostheses over the Universal Abutment. They are available in the following models: Provisional (Polymer) and Castable (Casting). The Universal Abutment Provisional Coping is polymer-based, single use and supplied sterile. The Universal Abutment Castable Coping is made from castable polymer material, single use and supplied non-sterile. However, when supplied as part of the Universal Abutment Kit, it's supplied sterile.

APPLICATIONS

The Provisional Coping is suitable for manufacturing temporary prostheses over the Universal Abutment. The Castable Coping is suitable for manufacturing metallic infrastructure of prostheses cemented over Universal Abutment.

CONTRAINDICATIONS

This product does not show contraindications, provided that it is properly used for the aimed purposes.

HANDLING

Provisional Coping: Place Coping on Universal Abutment Analog in the working model for indirect polymerization, or on the Abutment itself in the mouth for direct polymerization, in the temporary prosthesis manufacturing process. The indirect manufacturing technique may still go through the waxing stage before polymerization. However, in all three options, the material must undergo surface polishing before being cemented in the mouth with temporary

cement. Castable Coping: Place the Coping on the Universal Abutment Analog. Then, perform waxing and overcasting. The prosthesis metallic infrastructure is thus obtained.

PRESENTATION AND STERILIZATION

These products are packaged in blisters (film and surgical grade paper). The Universal Abutment Provisional Coping is single use and supplied sterile.

PRECAUTIONS

- As to the systemic aspect, consider the patient's general health condition according to the applicable literature. As to the local aspect, please observe the intra-oral tissues conditions.
- This product is of single use and it may not be re-sterilized.
- Reuse of this product may cause:
 - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
 - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
 - Do not use the product if the package is violated.
 - An insertion torque higher than that recommended and improper Drivers may cause damages to the material and make the system unusable.
 - Better results are achieved with the utilization of Neodent products. The utilization of instruments and/or Implants and/or Prosthetic Components of other systems do not assure the perfect function and voids any product guarantee.
 - Make sure that the parts are not swallowed or inhaled by the person who is receiving the product (patient).
 - A special precaution should be taken during the closure so that the interocclusal space is not insufficient or excessive.
 - The improper planning may endanger the Implant/Prosthesis set performance resulting in system faults,

such as implant loss or fracturing, loosening or fracturing of the Prosthetic Screws.

- During the installation, make sure to align it to the Implant insertion axle.
- Make sure that it is perfectly seated in the Implant. For that periapical radiographs are recommended with the parallelism technique.

ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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