ST.330.214

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English  Impression Coping

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Español  Transfer

Simbologia Symbology Simbología Simbologia Symbolik Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
	Requisito Brasileiro / Brazilian Requirement / Requisito Brasileño / Requisito Brasiliano / Brasilianische Anforderungen / Exigences brésiliennes
	Esterilizado por óxido de etileno / Sterile - ethylene oxide / Esterilizado con óxido de etileno / Sterile - ossido di etilene / Sterilisiert mit Ethylenoxid / Stérilisé à l'oxyde d'éthylène

Simbologia Symbology Simbología Simbologia Symbolik Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Não reutilizar / Do not reuse / No lo reutilice / Non riutilizzare / Nicht wiederverwendbar / Ne pas réutiliser / Nicht wiederverwendbar / Ne pas réutiliser
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Consérvelo seco / Conservare in luogo asciutto / Trocken lagern / Garder au sec
	Manter afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conservar à l'abri de la lumière solaire
	Não utilizar se a embalagem estiver danificada / Do not use if package is damaged / No lo utilice si el envoltorio está dañado / Non utilizzare se la confezione originale è danneggiata / Bei beschädigter Verpackung Produkt nicht verwenden / Ne pas utiliser si l'emballage a été abimé

Simbologia Symbology Simbología Simbologia Symbolik Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Marcação CE para comercialização na Comunidade Europeia / CE Mark for European Community market / Marca CE para comercialización en la Comunidad Europea / Marchio CE per la commercializzazione in Europa / CE Kennzeichen zum Vertrieb in der Europäischen Gemeinschaft / Marquage CE pour une commercialisation en Europe
	Notificação exigida pelo FDA para comercialização nos Estados Unidos / Notification required by FDA for United States market / Notificación exigida por el FDA para comercialización en los Estados Unidos / Notifica necessaria per la commercializzazione negli Stati Uniti / Vom FDA geforderte Kennzeichnung für den Vertrieb in den USA / Notification exigée par FDA pour une commercialisation aux États-Unis
	Representante na Comunidade Europeia / Representative in the European Community / Representante en la Comunidad Europea / Rappresentante nella Comunità europea / Vertreter in der Europäischen Gemeinschaft / Représentant dans la Communauté Européenne
	Fabricante / Manufacturer / Fabricante / Produttore / Hersteller / Fabricant

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

This product is a Prosthetic Component with the specific function of transfer (molding). The Impression Coping used over prosthetics components: Conical Abutment, Mini Conical Abutment, Micro Conical Abutment, Transepithelial Abutments, and Implants are manufactured in titanium alloy according to standard ASTM F136. It is indicated to the open tray (drag) or closed tray (replacement) techniques. It is presented as a single piece or two pieces, for Bridge or for Crown, in relation to the Prosthetic Component or Implant fitting. Universal Abutment Impression Coping: It is manufactured in polymer. It is indicated to the closed tray (replacement) technique. It is presented as a single piece for Crown, in relation to the Universal Abutment fitting. Smart Molding Impression Coping: It is manufactured in polymer. It is specific for Smart Implant line. It is indicated to the closed tray (replacement) technique. This component can be used as Provisional Coping for Smart Implant Mount.

APPLICATIONS

This product is a prosthetic component used to transfer the position and shape of the implant or prosthetic intermediary in accordance with each planned case, respecting diameters and tolerance indications of the components, aimed at aesthetics and function.

CONTRAINDICATIONS

This product has no contraindications, providing their intended purposes are correctly followed.

HANDLING

Fit the selected Neodent Impression Coping according to the planned molding technique, and install it on the Implant or Prosthetic Component, according to the planning. Open Tray Impression Coping: Fit the Impression Coping body over the prosthetic component or Implant selected. Fasten only the Screw using the 1.2mm Screwdriver until well stuck. Closed Tray Impression Coping:

Fasten the Impression Coping over the prosthetic component or Implant using the Impression Coping Driver. Make sure it is well seated on the component. Abutment Impression Coping: Fit the Impression Coping through a “click” on the Universal Abutment, then performing molding. Molding Impression Coping or Smart Provisional Coping: If you decide for molding, use the Impression Coping in the closed tray technique. In case of using the Impression Coping to make the provisional prosthesis, you should remove the triangle at the top of the Coping with a cutting tool. After this procedure you should use the Provisional Smart Cylinder which will be cemented over the Mount. The Smart Cylinder is connected to the Provisional Abutment (Mount) using a click system. If the provisional prosthesis is produced in lab, the specific Smart Implant Mount Analog must be used.

PRESENTATION AND STERILIZATION

This product is protected by a blister type packaging (film and medical grade paper) and supplied by unit, or if it is part of the Universal Abutment Kit, in a single package along with Provisional and Castable Copings of Universal Abutment. It is delivered sterile and It of single use. Sterilization method: ethylene oxide.

PRECAUTIONS

- This product is of single use.
- Reuse of this product may cause:
 - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
 - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
 - As to the systemic aspect, consider the patient’s general health condition according to the applicable literature. As to the local aspect, please observe the intra-oral tissues conditions.
- This product is of single use and it may not be re-sterilized.
- Do not use the product if the package is violated.
- This product must be used immediately after opening the packaging, at the moment of surgery. If it is not used, discard it.
- An insertion torque higher than that recommended and improper Drivers may cause damages to the material and make the system unusable.
- Better results

are achieved with the utilization of Neodent products. The utilization of instruments and/or Implants and/or Prosthetic Components of other systems do not assure the perfect function and voids any product guarantee.

- Make sure that the parts are not swallowed or inhaled by the person who is receiving the product (patient).
- The improper planning may endanger the Implant/prosthesis set performance resulting in system faults, such as Implant loss or fracturing, loosening or fracturing of the Prosthetic Screws.
- During the installation, make sure to align it to the Implant insertion axle, avoiding jamming and damage to the threads.
- Make sure that it is perfectly seated in the Implant. For that periapical radiographs are recommended with the parallelism technique.

ADVERSE EFFECTS

No adverse effects will occur, provided that the indications of this product are followed.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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