# ST.330.219

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Español  Tornillo Protésico

English  Prosthetic Screw

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Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
	Requisito Brasileiro / Brazilian Requirement / Requisito Brasileño / Requisito Brasiliano / Brasilianische Anforderungen / Exigences brésiliennes

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Esterilizado por óxido de etileno / Sterile - ethylene oxide / Esterilizado con óxido de etileno / Sterile - ossido di etilene / Sterilisiert mit Ethylenoxid / Stérilisé à l'oxyde d'éthylène
	Não reutilizar / Do not reuse / No lo reutilice / Non riutilizzare / Nicht wiederverwendbar / Ne pas réutiliser / Nicht wiederverwendbar / Ne pas réutiliser
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Consérvelo seco / Conservare in luogo asciutto / Trocken lagern / Garder au sec
	Manter afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conserver à l'abri de la lumière solaire

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

### DESCRIPTION

This product is manufactured in titanium alloy in accordance with ISO 5832-2 and ASTM F136. The Working Screw can also be manufactured in brass.

Prosthetic Screw	Description
Hex Screw* ▲	They have adequate geometry to be internally threaded to the Implant and with head seating in the inner portion of the corresponding intermediates (Abutments)
Angled Abutment Screw	
Conical Abutment Screw	
Mini Conical Abutment Screw	
Angled Mini Conical Abutment Screw	
Square Screw	
WS Universal Abutment Screw	
Transepithelial Abutment Screw	
Working Screw	They have adequate geometry to be internally threaded to the Analog and with head seating in the inner portion of the Coping in brass
CM Abutment Coping Screw and WS Abutment Screw*	They have adequate geometry to be internally threaded to the Abutment and with head seating in the inner portion of the corresponding Abutment Coping.
Conical Abutment Coping Screw* ▲	
Mini Conical Abutment Coping Screw and Micro Conical Abutment Screw* ▲	
Open-Tray Impression Coping Screw	They have adequate geometry to be internally threaded to the intermediate or Implant and with head seating in the inner portion of the Impression Copings.
Closed-Tray Impression Coping Screw	

\* Available with Neotorque cover: application of a thin film (covering), carbon-based which reduces the friction coefficient resulting an increase of the preload. ▲ Available in Torqceram model: recommended only for zirconia Prosthetic Components. The Screw head is purple and has conical seating. They present Neotorque cover.

### INDICATIONS FOR USE

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading.

### APPLICATIONS

Indicated as the following chart:

Prosthetic Screw	Indications
Hex Screw	Fix the intermediate to the Implant
Angled Abutment Screw	
Conical Abutment Screw	
Mini Conical Abutment Screw	
Angled Mini Conical Abutment Screw	
Square Screw	
WS Universal Abutment Screw	
Transepithelial Abutment Screw	
Working Screw	Fix the intermediate to the Implant Lab. use.
CM Abutment Coping Screw and WS Abutment Screw **	Fix the Prosthetic Coping to the intermediate
Conical Abutment Coping Screw	
Mini Conical Abutment Coping Screw and Micro Conical Abutment Screw	
Open-Tray Impression Coping Screw	Fix the transfer to the intermediate or Implant
Closed-Tray Impression Coping Screw	

\*\* They have adequate geometry to be internally threaded to GT Implant. Note: The Neodent Mini Conical Abutment Coping Screw 5.0 is used for CM Coping, CM Abutment Coping and WS Abutment.

### CONTRAINDICATIONS

The Neodent Prosthetic Screw does not present contraindications, providing their intended purposes and indications are correctly followed.

### HANDLING

It must be handled using adequate drivers and connections, as shown in table. The insertion torque values, when applied, must be respected:

Prosthetic Screw	Insertion Torque	Driver / Connection
Hex Screw	External Hex: 32 N.cm	Hex 1.2
	II Plus: 20 N.cm	
Angled Abutment Screw	External Hex: 32 N.cm	Hex 1.2
	II Plus: 20N.cm	
Conical Abutment Screw	20 N.cm	Ø3.3, 4.1 e 4.3: Connection for Prosthetic Abutments Ø5.0: Hex 1.6
Mini Conical Abutment Screw	20 N.cm	Ø3.3, 4.1 e 4.3: Connection for Prosthetic Abutments Ø5.0: Hex 1.6
Angled Mini Conical Abutment Screw	20 N.cm	Hex 1.2
Square Screw	32 N.cm	Square
WS Universal Abutment Screw	32 N.cm	Hex 1.2
Transepithelial Abutment Screw	20 N.cm	Connection for Prosthetic Abutments
Working Screw	-	Hex 1.2
CM Abutment Coping Screw and WS Abutment Screw	20 N.cm	Hex 1.2
Conical Abutment Coping Screw	10 N.cm	Hex 1.2
Mini Conical Abutment Coping Screw and Micro Conical Abutment Coping Screw	10 N.cm	Hex 1.2
Open-Tray Impression Coping Screw	-	Hex 1.2
Closed-Tray Impression Coping Screw	-	Impression Coping Driver

## TRACEABILITY LABEL

The product is identified through numerical codes arranged on the label (REF and LOT). The product is fitted with these codes, allowing their characteristics to be identified since when it is a raw material. This product is fitted with three labels which enable it to be tracked and should be posted to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor).

## PRESENTATION AND STERILIZATION

This product is packed by unit, protected by a blister packaging (film and medical grade paper). It is supplied sterile, by the Ethylene Oxide sterilization method. The Working Screws are non-sterile, they are indicated only for lab use.

## PRECAUTIONS

• Regarding the systemic aspect, consider the general condition of patient's health as the applicable literature. Related to the local aspect, observe the conditions of intra-oral tissues. • This product is intended for single use, except for the Working Screw. • Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products. • Do not use the product if the package is damaged. • Se produto estéril, este deve ser utilizado imediatamente após a abertura da embalagem, no momento da cirurgia. Se não utilizá-lo, descarte-o. • Insertion

torque higher than recommended and inadequate Connections can damage the material and make the system unusable. • Best results are obtained with the use of Neodent products. The use of instrumentals and/or Implants and/or Prosthetic Components of other systems do not ensure perfect function and disclaims any warranty of the product. • Ensure that the parts are not swallowed or inhaled by the person who is receiving the product (patient), when applicable. • Take particular care during the waxing to avoid excessive or insufficient interocclusal space. • The incorrect planning can compromise the performance of joint Implant / Prosthesis resulting in system failures, as a loss or fracture of the Implant, loosening or fracture of the Prosthetic Screws. • During installation, be sure to align it with the insertion axis of the Implant, preventing crashes and damages to the thread. • Be sure it is perfectly fixed to the implant. For this periapical radiographs with the parallelism technique is recommended, when applicable.

## ADVERSE EFFECTS

No adverse effects will occur, since the indications are respected.

## MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

## POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

## STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

## DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

## DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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