ST.330.222

Português  Posicionador Radiográfico

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
English  X-Ray Positioner

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

The Neodent X-Ray Positioner has as function evaluate and ensure through periapical radiograph the depth of osteotomy in installation procedures of Neodent Implants. This product is frequently used in flapless surgery, ie, after extraction and immediate implant installation. As Titanium has radiopaque properties, the channels which indicate the depth achieved (Implant length) are evident on radiograph. This product are manufactured in titanium alloy according to standard ASTM F136, supplied non-sterile and can be reused in accordance with the guidelines contained in these usage instructions.

APPLICATIONS

This product is an instrument used in dental surgery to verify depth in flapless osteotomy. It is suggested to perform periapical radiograph for evaluation.

CONTRAINDICATIONS

This product has no contraindications, providing its intended purposes and indications are correctly followed.

HANDLING

After perform the planned osteotomy of the Implant, the instrument should be inserted in the perforated cavity, performing the radiograph for evaluation.

SANITATION

This product should be correctly cleaned after each use. Do not leave it in contact with dampness for prolonged periods. For cleaning, follow these steps: Step 1: Submerge the parts fully into enzymatic detergent (diluted according to the manufacturer). Step 2: Wash in ultrasonic washer for approximately 10 to 15 minutes. Step

3: Rinse with plenty of distilled water until the solution wastes are fully removed. The utilization of nylon brushes is recommended. Step 4: Dry with a clean and dry cloth or with compressed air. Step 5: Perform a visual inspection, watching if there are faults in the cleaning process. If any wastes are still left, the part should be submerged again in detergent - step 1 - and if necessary, the cleaning should be made with the help of a nylon brush. Repeat the rinsing and drying sequence. Step 6: If applicable, place the product in a suitable case duly sanitized and dried. Step 7 - Select the package according to the sterilization process. Use preferably a self-sealing package if surgical grade paper with laminated film. ATTENTION: Do not use descaling products and avoid drying parts still containing cleaning solution wastes because such procedures favor oxidation. Also, the use of steel brushes is not recommended.

STERILIZATION

Sterilize the product on the day before or on the day of the procedure. We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard ISO 17665-1:2006 • Sterilization of health care products. Moist heat. Requirements for the development, validation and routine control of a sterilization process for medical devices. Suitable parameters: 121°C, 30 min, at 1 atm. Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Sterilization validity: 7 to 15 days, provided that the product is conditioned in a clean, dry environment away from the sunlight.

PRESENTATION AND STERILIZATION

This product is protected by a blister type packaging (film and medical grade paper) and supplied by unit. This product is unitarily provided, conditioned in blister (surgical grade film and paper type) packing. It is reusable, delivered non-sterile and should, therefore, be sterilized before use.

PRECAUTIONS

- Do not use the product if the package has been violated.
- This product should be used in its sterile condition.
- Better results are achieved with the utilization of the Neodent products sequence. The utilization of instrumentations and/or Prosthetic Components from other manufacturers does not assure the perfect function of the Neodent Implant System and voids any product guarantee.
- The improper planning may endanger the Implant/Prosthesis set performance resulting in faults in the system, such as Implant loss or fracture, loosening or fracturing of the Prosthetic Screws.

ADVERSE EFFECTS

Adverse effects will only occur if the instruments are chosen incorrectly.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

LIFE CYCLE

This product should be discarded when it loses its functionality.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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