ST.330.243

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English  Castable Attachment

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Español  Attachment Calcinable

Français  Attachment Calcinable

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
	Requisito Brasileiro / Brazilian Requirement / Requisito Brasileño / Requisito Brasileiro / Brasilianische Anforderungen / Exigences brésiliennes

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Não estéril / Non-sterile / No estéril / Non sterile / Nicht steril / Non stérile
	Não reutilizar / Do not reuse / No lo reutilice / Non riutilizzare / Nicht wiederverwendbar / Ne pas réutiliser / Nicht wiederverwendbar / Ne pas réutiliser
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Consérvelo seco / Conservare in luogo asciutto / Trocken lagern / Garder au sec
	Mantener afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conserver à l'abri de la lumière solaire

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

Prosthetic component used in the laboratory to obtain over denture retaining devices to be cemented to the roots of natural teeth or to be arranged in casting standards to obtain Bars, retained by teeth or Implants, to be also used in the retention of over dentures or partially removable prosthesis. Available in two types: Ball or Equator O-ring. In the upper portion it shows the chosen design and in the lower portion it has a conical design to be positioned in the wax standard. It is available in castable polymer.

APPLICATIONS

Indicated for removable prosthesis fitted to the teeth or installed in Implants installed in the maxilla or jaw. It can be used for partially removable prostheses or over dentures. It can be used in early or late rehabilitation procedures.

CONTRAINDICATIONS

It's contraindicated as supports for fixed prostheses or for direct use without a cast.

HANDLING

Hold the component and place it in a standard cast as indicated. In the case of multiple elements it will be necessary to use a prosthetic liner to achieve parallelism between the prosthetic components. The mixed standard should be included and merged with the alloy of choice. After divesting, finishing and polishing, it will be installed in the mouth. In case it has been used for obtaining cores for cementing in natural teeth, the cementing will take place and then the same retention components, of the female type withholding are fastened to the removable prosthesis. In the case of use for obtaining retaining bars in implant-retained over dentures, the Attachments are positioned according to the axis of insertion considering a future prosthetic

study and, after the bar is obtained by casting, it is screwed to the Implants, or intermediaries, and the retaining components are captured in the prosthesis.

PRESENTATION AND STERILIZATION

It's supplied by unit in a protective blister package made of medical grade film and paper. It is intended for single use and supplied non-sterile.

PRECAUTIONS

- Regarding the systemic aspect, consider the general health of the patient according to the applicable literature. As for the local aspect, observe the conditions of the intra-oral tissues.
- This product is single use and cannot be reused in more than one casting.
- The alloy chosen as well as the technique adopted for the casting procedure should enable obtaining of the prosthesis with the same dimensional characteristics as the castable component offered. Remember that prosthesis obtained (in the form of a core or bar) will have to adapt perfectly to the female type components, which will be used in the removable prosthesis.
- Reuse of this product may cause:
 - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
 - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
 - Do not use the product if the package has been violated. -Better results are achieved with the utilization of the Neodent products sequence. The utilization of instrumentations and/or prosthetic components from other manufacturers does not assure the perfect function of the Neodent Implant System and voids any product guarantee.
 - Ensure that the parts are not swallowed or breathed by the person receiving the product (patient).
 - Especial care should be taken when choosing the component so that the interocclusal space is neither insufficient nor excessive.
 - The number of components used

must be compatible with the performance you want from them and with the resistance setting of the retaining elements.

- Improper planning can compromise the performance of all the teeth or implants/prosthesis resulting in system failures, such as the loss or breakage of a tooth/implant, the loosening or fracture of a prosthetic screw.
- Make sure it is seated perfectly to the tooth prior to cementation or also the seating of the bar to the implants/intermediaries. For this we recommended using periapical x-rays with the use of a parallelism technique.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and the use of instrument are improper.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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