# ST.330.244

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Español  Implante Drive Ti

English  Drive Ti Implant

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Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
	Requisito Brasileiro / Brazilian Requirement / Requisito Brasileño / Requisito Brasiliano / Brasilianische Anforderungen / Exigences brésiliennes
	Produto esterilizado por radiação gama / Product sterilized through gamma rays / Producto esterilizado por radiación gama / Prodotto sterilizzato attraverso raggi gamma / Sterilisiert durch Gammastrahlung / Produit stérilisé par rayonnement gamma

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Não reutilizar / Do not reuse / No lo reutilice / Non riutilizzare / Nicht wiederverwendbar / Ne pas réutiliser / Nicht wiederverwendbar / Ne pas réutiliser
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Consérvelo seco / Conservare in luogo asciutto / Trocken lagern / Garder au sec
	Manter afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conservar à l'abri de la lumière solaire
	Não utilizar se a embalagem estiver danificada / Do not use if package is damaged / No lo utilice si el envoltorio está dañado / Non utilizzare se la confezione originale è danneggiata / Bei beschädigter Verpackung Produkt nicht verwenden / Ne pas utiliser si l'emballage a été abîmé

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Marcação CE para comercialização na Comunidade Europeia / CE Mark for European Community market / Marca CE para comercialización en la Comunidad Europea / Marchio CE per la commercializzazione in Europa / CE Kennzeichen zum Vertrieb in der Europäischen Gemeinschaft / Marquage CE pour une commercialisation en Europe
	Notificação exigida pelo FDA para comercialização nos Estados Unidos / Notification required by FDA for United States market / Notificación exigida por el FDA para comercialización en los Estados Unidos / Notifica necessaria per la commercializzazione negli Stati Uniti / Vom FDA geforderte Kennzeichnung für den Vertrieb in den USA / Notification exigée par FDA pour une commercialisation aux États-Unis
	Representante na Comunidade Europeia / Representative in the European Community / Representante en la Comunidad Europea / Rappresentante nella Comunità europea / Vertreter in der Europäischen Gemeinschaft / Représentant dans la Communauté Européenne
	Fabricante / Manufacturer / Fabricante / Produttore / Hersteller / Fabricant

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

## DESCRIPTION

The Drive Ti Implant is produced in commercially pure Titanium (Level 4). Macro-surface: • Outer Hexagon prosthetic interface with inner torque; •  $\varnothing$ 3.5, 4.3 or 5.0 Prosthetic Interface; • Double, progressive threads; • Cervical conicity to help reducing the compressive stress to the bone; □Cutting chambers in counterclockwise direction; • Blade-shaped apical thread; • High capacity of bone compression. Micro-surface: It is available in two surface technologies: NeoPoros or Acqua. • NeoPoros: porous surface obtained by double treatment: abrasive blasting and acid subtraction. • Acqua: hydrophilic surface obtained through a physical-chemical process aggregated to the NeoPoros surface topography. The result is an increased surface energy for titanium.

## INDICATIONS FOR USE

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading.

## APPLICATIONS

The Drive Ti Dental Implant is indicated for surgical installation to a bone which density is III or IV, according to Lekholm & Zarb's jawbone quality classification (1985). It may be used to support single or multiple prostheses in a conventional or immediate rehabilitation procedure. It may be installed immediately after the extraction of a natural tooth. Lekholm & Zarb's jawbone quality classification: • Bone Type I = Homogeneous compact bone. • Bone Type II = Compact bone with core of trabecular bone. • Bone Type III = Thin cortical bone layer surrounding dense trabecular bone of favorable strength. • Bone Type IV = Thin cortical bone layer surrounding low density trabecular bone. Note: For application of immediate load the primary stability must attain at least 32N.cm.

## CONTRAINDICATIONS

Not recommended in the presence of acute inflammatory or infectious processes in live tissue, poor bone volume or quality, remains of roots in the place, serious medical problems such as; bone metabolism disorders, blood clotting disorders, unsuitable healing capacity, insufficient oral hygiene, incomplete jawbone growth, patient uncooperative and not motivated, abuse of drugs or alcohol, psychosis, prolonged functional disorders which resist any treatment with medications, weakened immunological system, diseases which require the use of steroids, allergy or hypersensitivity to titanium, endocrinological diseases, allergy to drugs, bone already exposed to radiation, diabetes mellitus, anticoagulation medications/ hemorrhagic diathesis, teeth grinding, parafunctional habits, anatomically unfavorable bone situation, tobacco abuse, uncontrolled periodontitis, treatable pathological maxillary diseases, alterations of the oral mucous and pregnancy.

## HANDLING

Preparation of the recipient bed: this should be performed with the necessary precautions to reduce surgical trauma, including Drills in good cutting conditions, spinning between 500 and 800rpm, and under copious irrigation. Sequence of drill bits: Select the sequence of drill bits according to the diameter of the planned dental implant. For more information, see the Neodent catalog. Countersink drill bits should be used in adherence to the measure of the prosthetic interface.

## INSTALLATION OF THE DENTAL IMPLANT

- Remove the blister from the packaging.
- After opened the sterility seal, holds the primary packaging with the non-dominant hand and removes the cover. • For 4.1/4.3mm platform, capture the Implant using the 4.1/4.3 Ti Connection - Contra-Angle (blue), for 3.5mm platform capture the Implant using the 3.5 Ti Connection - Contra-Angle (purple) or for 5.0mm platform capture the Implant using the 5.0 Ti Connection - Contra-Angle (yellow) seeking the perfect fit between the Driver and the Implant. • Move the Implant to the bone bed. • In order to install it with the motor, use the Ti Connection - Contra-Angle according the diameter of the Implant chosen. For manual installation use the short or long 4.1/4.3, 3.5 or

5.0mm Ti Connection - Torque Wrench according the diameter of the Implant chosen. • In the surgical motor, use a 35N.cm maximum torque and rotation from 15 to 30rpm. • Preferably, complete the installation of the Implant with the Torque Wrench. • The maximum final torque recommended is 60N.cm. • The choice among installing the Cover Screw, Healing Abutment or Prosthetic Component of Ti ( external hex) line is up to the professional.

## PROSTHETIC SEQUENCE

Select the intermediaries between the Implant and the prosthesis, adhering to its indications and limitations, according to the applicable literature.

## TRACEABILITY LABEL

The product is identified through numerical codes arranged on the label (REF and LOT). The product is fitted with these codes, allowing their characteristics to be identified since when it is a raw material. This product is fitted with three labels which enable it to be tracked and should be posted to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor).

## PRESENTATION AND STERILIZATION

This product is supplied sterile (sterilization method: gamma rays) and is intended for single use, unitarily packaged in a container offering triple protection: tertiary packaging (cardboard), blister-type secondary packaging (surgical grade paper and film) and primary packaging (clear tube).

## PRE-OPERATORY AND PROSTHETIC PLANNING

The model, diameter, length, position and quantity of implants must be selected individually, taking into consideration the anatomy, quality and quantity of bone and space available. When necessary execute the diagnostic wax-up of the patient. In situations in which there are relatively high loads, special care must be taken to ensure the suitable alignment of the implant /component and prosthesis. The maximum angulation allowed for Neodent implants is up to 30° degrees.

## PRECAUTIONS

• Note the conditions of the intra-oral tissue, the bone quality and quantity of the bed receiving the Implant, by means of radiographic and/or tomography examinations. The Non-execution of the pre-surgical evaluation can lead to it being impossible to perceive pre-existing diseases. • Regarding the systemic aspect, consider the general state of health of the patient, who must undergo a careful clinical and radiological examination before surgery, assessing the psychological and physical status. Special attention should be paid in cases of patients who have local or systemic factors which may interfere with the healing processes of the bone or soft tissue, or in the process of bone integration (refer to CONTRAINDICATIONS). • The handling as well as all the material to be used during the procedure must be sterile. • This product is of single use and cannot be re-sterilized. • Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products. • Do not use the product if the packaging is broken. The sterilization is only guaranteed if the secondary packaging is not damaged. • This product must be used immediately after opening the packaging, at the moment of surgery. If it is not used, discard it. • Do not use the product with the validity expired. • For application of immediate load the primary stability must attain at least 32N.cm; if there is lower torque do not apply load. • See the torque to be given on the prosthetic component to be used. The excess or lack of torque can lead to undesirable result. • Implants of reduced diameter (e.g. 3.5mm) and angled components are not recommended in posterior regions of the mouth. • The maximum angulation allowed for

Neodent implants is up to 30° degrees. • Insertion torque greater than that recommended can make the system inoperative. • Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear. • Adjust the surgical motor used in the procedure (torque and RPM) as per specification of the implant to be used. • Always use the Neodent product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty. • Ensure against product aspiration. • Some of the possible causes of unsuccessful bone integration and loss of the prosthesis during the treatment are the following: unsuitable osteotomy, infections, deficient oral hygiene, occlusal trauma, systemic problems or diseases, low remaining bone quantity or quality, lack or failure of irrigation, use of instruments which are not specific and/or without cutting power and lack of specific training. • Inappropriate planning can jeopardize the performance of the Implant/prosthesis unit leading to failures of the system, as loss or fracturing of the Implant, loosening or fracturing of the Prosthetic Screws. • The lack of adjustment or inappropriate occlusal adjustment can endanger the implant/prosthesis unit. • It is the dentist's responsibility to use the Neodent products in compliance with the instruction for use, as well as establish if it is suited to the individual situation of each patient.

## ADVERSE EFFECTS

The installation of dental implants, as well as any other surgical procedure, may cause a slight discomfort and localized edema. More persistent symptoms can occur such as: chronic pain related to the dental implant, permanent paresthesia, dysesthesia, loss of maxillary/mandibular bone edge, systemic or localized infection, oroantral or oronasal fistula, adjacent teeth affected unfavorably, irreversible damage to adjacent

teeth, fracturing of the implant, jaw, bone or prosthesis, esthetic problems, injury of the nerves, exfoliation, hyperplasia.

## MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

## POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

## STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

## DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

## DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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