ST.330.245

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions. **Neodent** titanium surgical instruments are used in dental surgeries with presentation corresponding to their indication.

DESCRIPTION

Abutment Try-ins are prosthetic components made of titanium alloy, according to standard ASTM F136, which aid in selection of the prosthetic component to be installed on the implant. They are available in different colors to facilitate the identification of each family (CM, WS, and Facility) and also to distinguish between the two heights of the cementable portion (4 and 6 mm).

For the CM line, they are available for selection of Angled Abutments and Anatomical Abutments, and for the Facility line, they are available for selection of Anatomic Abutments, according to the table below:

Interface	Transmucosal height (mm)	Color	Height of the cementable portion (mm)
CM Universal Abutment	0.8, 1.5, 2.5, 3.5, 4.5, 5.5, 6.5	Yellow	4
		Dark blue	6
CM Universal Abutment 17° and 30°	1.5, 2.5, 3.5	Dark blue	6
WS Universal Abutment	0.8, 1.5, 2.5, 3.5	Green	6
CM Anatomical Abutment	1.5, 2.5, and 3.5	Purple	N.A.
CM Lateral Anatomical Abutment		Light blue	N.A.
Facility Anatomical Abutment	1.5, 2.5, and 3.5	Yellow	N.A.

INDICATIONS FOR USE

The Abutment Try-In is indicated in the procedure for selection of prosthetic component to be installed over implants. It may be used with singlestage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading.

APPLICATIONS

They are indicated for selecting the diameter and angle of the Prosthetic Component, the transmucosal height (mucosa), as for assessing the height of cementable area, thereby selecting the best prosthetic component for each case.

CONTRAINDICATIONS

This product is not indicated for customization.

HANDLING

Position the Abutment Try-In over the Implant of the CM, WS, or Facility line with the aid of the 1.2 mm Manual Screwdriver. During installation, apply pressure over the Abutment Try-In to make sure that the assembly is secure. After installation, perform a clinical evaluation and, when necessary, an X-ray evaluation, to ensure the seating of the Abutment Try-In on the Implant and confirm its dimensions. Remove the Abutment Try-In with the aid of the 1.2 Manual Screwdriver. Proceed with installation of the Abutment selected.

SANITATION

This product must be correctly sanitized after each use.

For sanitization, proceed as follows: 1st step: Immerse completely in enzymatic detergent (diluted according to the manufacturer). 2nd step: Wash in an ultrasound washer for approximately 10 to 15 minutes. 3rd step: Rinse with plenty of distilled water, until the solution residues are completely removed. It is recommended to use nylon brushes. 4th step: Dry with a clean, dry cloth or with compressed air. 5th step: Conduct a visual inspection, observing whether there are failures in the cleaning process. If there are still residues, the piece must be immersed in detergent again - 1st step - and, if necessary, the cleaning should be done with the aid of a nylon brush. Repeat the sequence of rinsing and drying. 6th step: Select the packaging according to the sterilization process. Preferable use self-sealing surgical grade paper packaging with laminated film. **ATTENTION:** Do not use descaling agents, and avoid drying pieces that still contain residues of cleaning solution, because these procedures promote oxidation. The use of steel brushes is not recommended.

PRESENTATION AND STERILIZATION

This product is reusable and supplied non-sterile, packed individually. This product must be correctly sanitized and sterilized before use. Sterilize the product on the day before or on the day of the procedure. We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard BS EN ISO 17665-1:2006 (Sterilization of health care products. Moist heat Requirements for the development, validation and routine control of a sterilization process for medical devices). Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Suitable parameters: 121 °C, 30 min, at 1 atm. Sterilization validity: 7 to 15 days, provided that the product is conditioned in clean, dry and protected from the sunlight environment.

PRECAUTIONS

- It is recommended to use periapical x-rays for evaluation of the relation between the contour of the bone ridge and the Abutment.
- Be aware in cases of patients that show signs of allergy or hypersensitivity to the chemical components of the material: titanium alloy Ti6Al4V-ELI.
- Surgical planning and/or inadequate prosthesis can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of components and/or prosthetic screws.
- Do not use the product if the packaging is broken.
- Do not use the product after expiry date.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- During installation, make sure to align it with the insertion axis of the implant. Make sure it is perfectly seated on the Implant. To do so, it is recommended to use periapical X-rays with the parallel technique.
- Before each procedure check the conditions of the **Neodent** surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the **Neodent** product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the **Neodent** Implant System and exempts any product warranty.
- It is the professional's responsibility to use the **Neodent** products according to the instructions for use.

ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored, in its original packaging, in a clean and dry location, in a maximum temperature of 40° C and protected from direct sunlight.

DISPOSAL OF MATERIAL

After use, all the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

LIFE CYCLE

This product should be disposed of when it loses its functionality.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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