# ST.330.246

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

### DESCRIPTION

The Zygomatic Labial Protector is manufactured in Stainless Steel, and it has ergonomically designed geometry for use during the surgery of Zygomatic Implants. It has a rod with holdbacks for better hardness during the use and in its end there is a tube through which the Drills for Zygomatic Implant are passed.

### APPLICATIONS

It is indicated for providing a higher protection of the patient's lip during the surgical drilling for Zygomatic Implants.

### CONTRAINDICATIONS

This product does not present contraindications, provided that it is used properly for the assigned purposes.

### HANDLING

The Zygomatic Labial Protector is used along with the Drills for Zygomatic Implants. The Drills are fastened to the Contra-Angle and next, they cross the Zygomatic Labial Protector tube. During all the drillings it continues involving the Drills so that the patient lip is not injured. Such procedure is repeated in all Drills models and diameters.

### SANITATION

This product should be correctly sanitized after each use. Proceed as follows: 1st step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer). 2nd step: Wash in ultrasonic cleaner for approximately 10 to 15 minutes. 3rd step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is recommended. 4th step: Dry with a clean and dry

cloth or with compressed air. 5th step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent - 1st step - and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence. 6th step: Select the packing in accordance to the sterilization process. Use preferably surgical grade paper self-sealing packing with laminated film. ATTENTION: Do not use anticrust products, avoid using enzymatic liquid in concentrations above 10% and dry parts that still have cleaning solution wastes, because these procedures cause oxidation. The use of steel brushes is also inadvisable.

### PRESENTATION AND STERILIZATION

This product is conditioned in a blister type package (surgical grade paper and film) sold unitarily. It is not recommended to sterilize the product in autoclave in its original package. It is reusable, delivered non-sterile and it should, therefore, be sterilized before its use. Sterilize the product on the eve or in the procedure day. Preferably, it is recommended to follow the steam autoclave sterilization method and also, the parameters and procedures established in the standard BS EN ISO 17665-1:2006 (Sterilization of health care products. Moist heat. Requirements for the development, validation and routine control of a sterilization process for medical devices). Recommended parameters: 121 °C, 30 min, at 1 atm. Alternatively, you may proceed according to the autoclave manufacturer's instruction. Sterilization validity: 7 to 15 days, provided that the product is conditioned in a clean, dry environment away from the sunlight.

### PRECAUTIONS

• Do not use the product if the package is violated. • This product should be used in its sterile condition. • Better results are achieved with the utilization of Neodent products. -The utilization of instruments and/or Implants and/or Prosthetic Abutments of other systems does

not assure the perfect function and voids any product guarantee. • The improper planning may endanger the implant/prosthesis set performance resulting in system faults, such as Implant loss or fracturing, loosening or fracturing of the Prosthetic Screws.

### ADVERSE EFFECTS

Adverse effects will only occur if the selection and use of the instrumentation are improper.

### POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

### STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40 °C and protected from direct sunlight.

### DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

### DATE OF EXPIRATION

Written on the label.

### LIFE CYCLE

This product should be disposed of when it loses its functionality.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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