# ST.330.263

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This device is intended for specialized lab procedure, which should only be executed by qualified professionals with specific training in CAD/CAM technique of scanning models. To avoid unsatisfactory results for the patient, do not use the product without mastering the appropriate techniques. Always perform the correct procedures in appropriate conditions and with suitable equipment.

## DESCRIPTION

The Scanbody is manufactured in high performance polymer and the Fixation Screw is manufactured in titanium alloy according to ASTM F136. The Scanbody is exclusively used during the scanning process of a plaster model, with the purpose to transfer to the CAD/CAM software the exact position and orientation of models Analogs. It's available with specific Screw, can be rotational (for bridge) or anti-rotational (for crown) in relation to the prosthetic component or Implant fitting.

## APPLICATIONS

This product is a prosthetic component used to transfer the position of Analogs or intermediaries components during the plaster model scanning through the CAD/CAM technique according each planned case, respecting diameters and components tolerance indications aiming aesthetics and functionality.

## CONTRAINDICATIONS

This product is not indicated for intraoral scanning.

## HANDLING

Fit the Scanbody over the Analog in the plaster model. Fasten the Screw using the 1.2mm Manual Screwdriver applying a torque (maximum torque of 10N.cm). Make sure the Scanbody is well fit in the Analog. There should be no vertical or rotational gap. Gaps indicate that the Scanbody must be repositioned. Follow the instructions delivered with the scanner to

scanning instructions. If you use a scanning spray, avoid spraying on Scanbody. Identification of the component by prosthetic interface:

Prosthetic interface	Identification number on the component
HE 3.3	1
HE 4.1/4.3	2
HE 5.0	3
II Plus	5
CM Exact	8
Micro Pillar	9
Conical Mini Pillar 4.1/4.3	10
Conical Mini Pilar 5.0	11
Conical Pillar 4.1/4.3	12
Conical Pillar 5.0	13
Pillar CM/GT/WS	14
WS	15

## FORMS OF PRESENTATION AND STERILIZATION

This product is conditioned in a blister type packaging (film and medical grade paper). It goes along with a titanium Screw. It is delivered non-sterile and for single use.

## PRECAUTIONS

- Before using the product, make sure that it clean and has the same prosthetic platform of the Analog to be scanned.
- Make sure the Scanbody used is compatible with Analog being scanned.
- This product is for single use and must not be reused and reprocessed.
- Check carefully if the Scanbody is well adapted to the Analog of the model, with no rotational or vertical gaps.
- Do not use the product if the packing is damaged.
- Overtightening and improper Drivers may cause damages to the material and make the system unusable. (see Handling)
- The use of Scanbody with a damaged platform can lead to unsatisfactory scanning results.
- During

the installation, make sure to align it to the Analog insertion axle, avoiding crashes and system damage

- The improper planning may endanger the Implant/Prosthesis set performance resulting in system faults, such as implant loss or fracturing, loosening or fracturing of the Prosthetic Screws.
- Always test the adaptation of the prosthesis made on the plaster model.

## ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

## POSTOPERATIVE CARE AND MAINTENANCE

Instruct the patient as to the need of performing a professional follow-up and complying with the instructions on precautions, hygiene and medication prescription. Those instructions are the professional's responsibility.

## STORAGE CONDITIONS

This product should be stored in a clean and dry location, at a maximum temperature of 40°C and protected from the direct sunlight radiation.

## DISPOSAL OF MATERIALS

All products and consumption materials used in the surgery for the installation of dental implants or prosthetic components may endanger the health of those who handle them, after their utilization. Before disposing them in the environment, the consultation and fulfillment of the legislation in force is recommended.

## SHELF-LIFE

Suitable on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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