ST.330.266

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Italiano  Kit

English  Cases

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Español  Estuches

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Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
	Não estéril / Non-sterile / No estéril / Non sterile / Nicht steril / Non stérile

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Consérvelo seco / Conservare in luogo asciutto / Trodden lagern / Garder au sec
	Mantener afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conserver à l'abri de la lumière solaire
	Não utilizar se a embalagem estiver danificada / Do not use if package is damaged / No lo utilice si el envoltorio está dañado / Non utilizzare se la confezione originale è danneggiata / Bei beschädigter Verpackung Produkt nicht verwenden / Ne pas utiliser si l'emballage a été abimé

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Marcação CE para comercialização na Comunidade Europeia / CE Mark for European Community market / Marca CE para comercialización en la Comunidad Europea / Marchio CE per la commercializzazione in Europa / CE Kennzeichen zum Vertrieb in der Europäischen Gemeinschaft / Marquage CE pour une commercialisation en Europe
	Notificação exigida pelo FDA para comercialização nos Estados Unidos / Notification required by FDA for United States market / Notificación exigida por el FDA para comercialización en los Estados Unidos / Notifica necessaria per la commercializzazione negli Stati Uniti / Vom FDA geforderte Kennzeichnung für den Vertrieb in den USA / Notification exigée par FDA pour une commercialisation aux États-Unis
	Representante na Comunidade Europeia / Representative in the European Community / Representante en la Comunidad Europea / Rappresentante nella Comunità europea / Vertreter in der Europäischen Gemeinschaft / Représentant dans la Communauté Européenne
	Fabricante / Manufacturer / Fabricante / Produttore / Hersteller / Fabricant

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

These are boxes for the storage of surgical and/or prosthetic instruments. They are made in polymer that can be sterilized in autoclave. They are provided with silicon supports for conditioning and safe fastening of each tool. They have markings guiding the tools use during the procedures. The Case may be equipped according to the work procedure.

APPLICATIONS

They are intended for safe conditioning of surgical and/or prosthetic instruments, as well as support for their sterilization.

CONTRAINDICATIONS

It is contraindicated for sterilization under dry heat (ovens).

HANDLING

Remove the instrumentation access cover pressing lightly the locking system between the base and the cover. To close the Case, position the cover on the base and press on it. For sterilization, the disinfected, clean and dry instruments should be conditioned in the Case, which should be packed in surgical grade crepe paper or in plastic bags suitable for autoclave. During the sterilization, the Case should not lean on the autoclave walls.

SANITATION

The carrying cases should be properly sanitized after its utilization. Do not let the product contact the moisture beyond the time required for the sanitation. For such, you should proceed as follows: Step 1: Disassemble the Case. Step 2: Submerge the part fully into enzymatic detergent (diluted according to the manufacturer). Step 3: Wash in ultrasonic washer for approximately 10 to 15 minutes. Step 4: Rinse with plenty of distilled water until the solution wastes are fully removed. The utilization of nylon brushes is

recommended. Step 5: Dry with a clean and dry cloth or with compressed air. Step 6: Perform a visual inspection, watching if there are faults in the cleaning process. If any wastes are still left, the part should be submerged again in detergent - step 1 - and if necessary, the cleaning should be made with the help of a nylon brush. Repeat the rinsing and drying sequence. Step 7: If applicable, place the product in a suitable Case duly sanitized and dried. Step 8: Select the package according to the sterilization process. Use preferably surgical grade paper self-sealing package with laminated film. ATTENTION: Do not use descaling products and avoid drying parts still containing cleaning solution wastes because such procedures favor oxidation. The use of steel brushes is not recommended.

PRESENTATION AND STERILIZATION

This product is supplied unitarily, conditioned in plastic package. It is reusable, delivered non-sterile and it should, therefore, be sterilized before its use. Sterilize the product on the eve or in the procedure day. Preferably, it is recommended to follow the steam autoclave sterilization method and also, the parameters and procedures established in the standard BS EN ISO 17665-1:2006 (Sterilization of health care products. Moist heat. Requirements for the development, validation and routine control of a sterilization process for medical devices). Recommended parameters: 121°C, 30 min, at 1 atm. Alternatively, you may proceed according to the autoclave manufacturer's instruction. Sterilization validity: 7 to 15 days, provided that the product is conditioned in a clean, dry environment away from the sunlight.

PRECAUTIONS

- This product should be used in the sterile condition.
- Do not sterilize under dry heat (oven).
- The Case should not lean on the Autoclave walls, to avoid the plastic from melting.
- Place the water measure for autoclave recommended by the manufacturer. The lack of water during the sterilization cycle may lead to the melting of the Case polymer as well as to the damage of the instruments.
- The contact of the Case with acrylic resin liquid results in damages of the Case polymer structure.
- The chemical

sterilization is not recommended.

- Better results are achieved with the utilization of Neodent products. The utilization of instruments and/or implants and/or prosthetic abutments of other systems do not assure the perfect function and voids any product guarantee.
- The improper planning may endanger the implant/prosthesis set performance resulting in system faults, such as implant loss or fracturing, loosening or fracturing of the prosthetic screws.

ADVERSE EFFECTS

Adverse effects will only occur if the selection and use of the product are improper.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

LIFE CYCLE

This product should be disposed of when it loses its functionality.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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