ST.330.268

Português 

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

Four exclusive Prosthetic Abutments are required for accomplishing the One Step Hybrid Technique: Brass Coping: a Prosthetic Abutment for laboratory use having a tapered geometry in the outer part for fitting the Castable Coping . In its inner part it has a slot compatible with the planned Prosthetic Abutment. Castable Coping: a cylindrical Prosthetic Abutment having retentions in its outer part for the waxing. In its inner part it has a slot compatible with the planned Prosthetic Abutment. Titanium Coping: a Prosthetic Abutment with outer retentions and treated surface (either porous or rugged) for retaining cement. It is the prosthetic intermediate of the One Step Hybrid technique that is cemented to the metallic infrastructure and screwed between the intermediate and the prosthesis. Working Screw - Long (laboratory room): with a 1.2 mm hexagon fit in one end and thread compatible with the prosthetic intermediate diameter in the other end.

INDICATIONS FOR USE

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading. Multiple tooth applications may be rigidly splinted.

APPLICATIONS

The One Step Hybrid Technique is indicated for multiple screwed prostheses when there is a need to obtain passivity in the metallic infrastructure fitting without using weld and independent from the metal dimensional change after the melting process.

CONTRAINDICATIONS

This product does not present contraindications provided that it is used properly for the stated purposes.

HANDLING

The Brass Coping placed on the Analog corresponding to the intermediate Prosthetic Abutment. The Castable Coping is placed on the Brass Coping, fastened by the Working Screw-Long one step hybrid. The Castable Coping is joined together by acrylic resin. Next the infrastructure waxing is performed. Once the waxing is completed, the Working Screw one step hybrid is removed and the infrastructure is submitted to the melting process with specific metallic alloy. The molten infrastructure is placed on the Brass Coping to check the fitting passivity. If necessary, internal wear is made in the region corresponding to the Castable Coping that is now molten in metal in order to reach the infrastructure passivity on the Brass Coping. Once the infrastructure passivity is achieved, internal retentions should be made in the molten Coping area, in order to promote mechanic retentions for the cement. A thin layer of primer (alloy primer) specific for metal should be applied in that area. The Brass Coping has a dimension slightly higher (0.10 mm) than the Titanium Coping that are replaced by the same and screwed on the respective Analogs with the matching Screws (short ones). The Titanium Coping are coated with an application of metal specific primer (alloy primer) in its outer part, and the screw input hole should be obliterated by wax in order to avoid the entry of resinous cement. A portion of dual resinous cement is applied to the outer part of the Titanium Coping, as well as in the inner part corresponding to the Coping in the infrastructure. It is cemented (pressed) on the Titanium Coping and one should immediately remove the cement surplus that leaks out from its orifice. After the cement setting, the infrastructure is unscrewed from the model and the surplus cement remaining in the Titanium Coping edge is removed. This cementation process may also be performed on the Coping's mouth provided that the requirements of the used cement are complied with. NOTE: For the cementation, we suggest the dual resinous cement Panavia F (Kuraray Co Ltd Tokyo-Japan) and the Alloy Primer (Kuraray Co Ltd Tokyo-Japan).

PRESENTATION AND STERILIZATION

They are conditioned in a blister type package (surgical grade paper and film), single use, non-sterile and delivered unitarily.

PRECAUTIONS

- Do not use the product if the package is violated.
- This product is of single use and it cannot be reused.
- Reuse of this product may cause:
 - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
 - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
 - The result shall only be achieved if the sequence of Neodent products is used. The utilization of instruments and/or implants and/or prosthetic abutments from other manufacturers does not assure the perfect function of the Neodent Implant System and voids any product guarantee.
 - The improper planning may endanger the Implant/prosthesis set performance resulting in system faults, such as Implant loss or fracturing, loosening or fracturing of the Prosthetic Screws and/or metallic infrastructure, or releasing of the titanium Coping from the metallic infrastructure.

ADVERSE EFFECTS

Adverse effects will only occur if the selection and use of the product are improper.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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