ST.330.269

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

The Provisional Coping/ SMART Provisional Coping are manufactured in polymer, indicated for cemented Provisional Prosthesis confection over intermediaries.

INDICATIONS FOR USE

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading. Multiple tooth applications may be rigidly splinted.

APPLICATIONS

Provisional Coping :is indicated to confection of Provisional Prosthesis over Universal Abutment. SMART Provisional Coping: : is indicated to confection of cemented Provisional Prosthesis over Provisional Abutment (SMART Implant Mount).

CONTRAINDICATIONS

This product does not show contraindications, provided that it be properly used for the shown purposes.

HANDLING

The Provisional Coping/Provisional SMART Coping is placed over the Analog in a work model for indirect acrylization, or over the own Prosthetic Component for direct acrylization. The technique of indirect confection still can pass by waxing

before acrylization. The Provisional Coping should be polished before cementation in mouth. Connect the Provisional SMART Cylinder over the Provisional SMART Abutment by either a snap connection/click system or by cement retention. The Provisional Coping/SMART Provisional Coping is intended to be used straight with no angulation.

PRESENTATION AND STERILIZATION

This product is supplied by unit in a protective blister package made of medical grade film and paper. Is intended for single use and supplied sterile. Sterilization method: Ethylene oxide.

PRECAUTIONS

- As to the systemic aspect, consider the patient's general health condition according to the applicable literature. As to the local aspect, please observe the intra-oral tissues conditions.
- This product is of single use and it may not be re-sterilized.
- Reuse of this product may cause:
 - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
 - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
 - Do not use the product if the package is violated.
 - An insertion torque higher than that recommended and improper Drivers may cause damages to the material and make the system unusable.
 - Better results are achieved with the utilization of Neodent products. The utilization of instruments and/or Implants and/or Prosthetic Components of other systems do not assure the perfect function and voids any product guarantee.
 - Make sure that the parts are not swallowed or inhaled by the person who is receiving the product (patient).
 - A special precaution should be taken during the closure so that the interocclusal space is not insufficient or excessive.
 - The improper

planning may endanger the Implant/Prosthesis set performance resulting in system faults, such as implant loss or fracturing, loosening or fracturing of the Prosthetic Screws.

- During the installation, make sure to align it to the Implant insertion axle.
- Make sure that it is perfectly seated in the Implant. For that periapical radiographs are recommended with the parallelism technique.

ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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