ST.330.272

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

The Distal Bar Coping and the Distal Bar are manufactured in titanium alloy according to the ASTM F136 standard. The Distal Bar Copings have external retentions for the fixation of acrylic resin. The Distal Bar has rectangular shape and is fixed in the Distal Bar Coping. It has holes along of its body facilitating the acrylic resin retention. The Distal Bar Copings are indicated for central regions. The distal Bar is indicated for posterior regions.

APPLICATIONS

This product is used for making temporary Prosthesis on Mini Conical Abutments or Conical Abutments. It can be used directly in the mouth or in the laboratory using plaster models and Abutments Analogs. It is necessary to respect the diameters and tolerance indications of the components.

CONTRAINDICATIONS

This product does not show contraindications, provided that it be properly used for the shown purposes.

HANDLING

This product can be used in two ways: 1) After the installing of the intermediary in the mouth over the Implant, transfer its position by means of a molding to a working model, which will have a replica or Analog of the chosen intermediary. On the Analog the titanium Coping will be positioned, and the latter will be made on the temporary Prosthesis. 2) The titanium Coping can be set directly in the mouth screwed on the intermediary Abutment, and then the fabrication of the temporary prosthesis will be done immediately on the Coping. In this case, remember that you must take care to avoid chemical burns to tissues adjacent to the Implant. For this, the Coping has a hollow at its base that allows the stabilization of a rubber dam to protect the adjacent soft tissues. In both cases, the final torque applied on the Prosthetic

Screw is 10 Ncm. Note: Observe the compatibility of the intraoral intermediary with the Analog and the selected Coping.

TRACEABILITY LABEL

The product is identified through numerical codes arranged on the label (REF and LOT). The product is fitted with these codes, allowing their characteristics to be identified since when it is a raw material. This product is fitted with three labels which enable it to be tracked and should be posted to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor).

PRESENTATION AND STERILIZATION

It is supplied unitarily, conditioned in a blister type package (surgical grade paper and film). It is delivered sterile and it is of single use. Sterilization method: ethylene oxide.

PRECAUTIONS

- Regarding the systemic aspect, consider the general health of the patient according to the applicable literature. As for the local aspect, observe the conditions of the intra-oral tissues.
- This product is of single use and it may not be re-sterilized. • Reuse of this product may cause:
 - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
 - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products. • Do not use the product if the package has been violated. • This product must be used immediately after opening the packaging, at the moment of surgery. If it is not used, discard it. • An insertion torque higher than that recommended and improper Drivers may cause damages to the material and make the system unusable. • This product should be used immediately after opening the package at the time of creation of the prosthesis in a clinic or laboratory. If you don't use it, discard it. • Better results are achieved with the utilization of the Neodent products sequence. The utilization of instrumentations and/or prosthetic components from other manufacturers does not assure the perfect function of the Neodent Implant

System and voids any product guarantee. • This Cylinder is not suitable for making permanent prosthesis. Not to be used in overcast situations.

- Ensure that the parts are not swallowed or breathed by the person receiving the product (patient).
- Improper planning can compromise the performance of all the teeth or Implants/ Prosthesis resulting in system failures, such as the loss or breakage of a tooth/Implant, the loosening or fracture of a Prosthetic Screw.
- During the installation, make sure to align it to the insertion axis, avoiding jamming and damage to the threads.
- Ensure that it is perfectly seated on the intermediary. For this we recommended using periapical x-rays with the use of a parallelism technique.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and the use of instrument are improper.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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