# ST.330.277

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

## DESCRIPTION

The Facility Bone Tap is made of stainless steel. At one end, it has a coupling for a Contra-Angle handpiece and a hex that allows it to fit with the Facility Bone Tap Ratchet Connection. The active portion of the Facility Bone Tap has cutting edges for forming threads and laser markings relating to the length of Facility Implants: the first three (from the active portion) refer to the lengths of the 10, 12, and 14 mm implants, respectively. The fourth marking has a length of 16 mm and is recommended for surgeries without flap opening or intraosseous installation.

The Facility Bone Tap Connection is made of stainless steel. At one end, it has a square socket for coupling with the Ratchet Wrench. At the other end, it has a geometry suitable for coupling with the Facility Bone Tap.

## APPLICATIONS

Facility Bone Tap is indicated for threads creation in surgical alveolus before the Facility Implants placement in type I or II bone bed, according to Lekholm & Zarb's jawbone quality classification (1985):

- Bone Type I = Homogeneous compact bone.
- Bone Type II = Compact bone with core of trabecular bone.
- Bone Type III = Thin cortical bone layer surrounding dense trabecular bone of favorable strength.
- Bone Type IV = Thin cortical bone layer surrounding low density trabecular bone.


The Facility Bone Tap Connection is indicated for manual placement using the Torque Wrench.

## CONTRAINDICATIONS

The Facility Bone Tap should not be used on bone types III and IV. The Facility Bone Tap should not be used during the installation of Implants other than Facility Implants.

## HANDLING

To use the Facility Bone Tap, first use the sequence of drill according to the planning and the implant selected. Start the procedure with the Facility Bone Tap connected to the Contra-Angle handpiece, set the motor to a speed of 20-30 rpm and torque of 35 N.cm. Insert the Facility Bone Tap in the surgical alveolus, maintaining the drilling axis, until it is stable. Couple the Facility Bone Tap Ratchet Connection to the Facility Bone Tap. Continue the procedure with the use of the Ratchet Wrench. Make clockwise insertion motions maintaining the drilling axis, with light pressure on the Ratchet Wrench, with maximum torque of 60 N.cm. Insert the Facility Bone Tap until reaching the marking for the implant chosen (see image below). It may be necessary to alternate the direction between insertion and removal a few times to form the thread.

ATTENTION: The figure is merely illustrative. It does not represent the real dimension and characteristics of the product.

To remove the Facility Bone Tap from the surgical alveolus, reverse the direction of the Ratchet Wrench to counter-clockwise and proceed with complete removal of the Bone Tap; this procedure ensures the perfect formation of the thread.

## SANITIZATION

This product should be correctly sanitized after each use. Proceed as follows:

- 1<sup>st</sup> step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer)
- 2<sup>nd</sup> step: Wash in ultrasonic cleaner for approximately 10 to 15 minutes.
- 3<sup>rd</sup> step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is recommended.
- 4<sup>th</sup> step: Dry with a clean and dry cloth or with compressed air.
- 5<sup>th</sup> step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent - 1<sup>st</sup> step - and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence.
- 6<sup>th</sup> step: Select the packing in accordance to the sterilization process. Use preferably surgical grade paper self-sealing packing with laminated film. ATTENTION: Do not use anticrust products, avoid using enzymatic liquid in concentrations above 10% and dry parts that still have cleaning solution wastes, because these procedures cause oxidation. The use of steel brushes is also inadvisable.

## PRESENTATION AND STERILIZATION

This product is reusable and provided non-sterile, and must therefore be sterilized before each use. It is packaged in single units in a package that offers double protection: secondary packaging (carton) and primary packaging (blister). It is not recommended to autoclave the product in its original packaging. Sterilize the product on the day before or on the day of the procedure. We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard BS EN ISO 17665-1:2006 (Sterilization of health care products. Moist heat. Requirements for the development, validation and routine

control of a sterilization process for medical devices). Suitable parameters: 121 °C, 30 min, at 1 atm. Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Sterilization validity: 7 to 15 days, provided the product is conditioned in clean, dry and protected from sunlight environment.

## PRECAUTIONS

- Perform the preparation of the surgical alveolus according to the instrument sequence of the Facility implant chosen and surgical planning.
- To use with the Contra-Angle handpiece, adjust the surgical motor used in the procedure to a maximum speed of 20-30 rpm.
- This product must be handled securely, avoiding contact with nearby soft tissues beyond the area of interest.
- Do not sterilize the product in its original packaging.
- Note the conditions of the intra-oral tissue, the bone quality and quantity of the bed receiving the Implant, by means of radiographic and/or tomography examinations. Failure to perform the pre-surgical assessment may compromise the success of the procedure.
- Regarding the systemic aspect consider the general state of health of the patient, who must undergo a careful clinical and radiological examination before surgery, assessing the psychological and physical status. Special attention should be paid in cases of patients who have local or systemic factors which may interfere with the healing processes of the bone or soft tissue, or in the process of bone integration. The Non-execution of the pre-surgical evaluation can lead to it being impossible to perceive ore-existing diseases.
- Inappropriate planning can jeopardize the performance of the implant/prosthesis unit leading to failures of the system, as loss or fracturing of the implant, loosening or fracturing of the prosthetic screws.
- The handling as well as all the material to be used during the procedure must be sterile.
- Do not use the product if the packaging is broken.
- Do not use the product with the validity expired.
- For immediate load application, check the indication of the implant installed.
- See the torque to be given on the prosthetic component to be used. The excess or lack of torque can lead to undesirable result.
- Insertion torque greater than that recommended can make the system inoperative.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Make the occlusal and interproximal adjustment after the installation of the prosthetic.
- The lack of adjustment or inappropriate occlusal adjustment can endanger the implant/prosthesis unit.
- Some of the possible causes of unsuccessful bone integration and loss of the prosthesis during the treatment are the following: unsuitable osteotomy, infections, deficient oral hygiene, occlusal trauma, systemic problems or diseases, low remaining bone quantity or quality, lack or failure of irrigation, use of instruments which are not specific and/or without cutting power and lack of specific training.
- Before each procedure check the conditions of the Neodent

surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear. • Always use the **Neodent** product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the **Neodent** Implant System and exempts any product warranty. • It is the professional's responsibility to use the **Neodent** products according to the instructions for use, and to determine whether it suits the individual situation of each patient.

#### **ADVERSE EFFECTS**

Adverse effects will only occur if the choice and the use of instrument are improper.

#### **POST-OPERATIVE PRECAUTIONS AND MAINTENANCE**

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

#### **STORAGE CONDITIONS**

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

#### **DISPOSAL OF MATERIAL**

Every product and consumable used during the surgery for the installation of dental implants

may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

#### **DATE OF EXPIRATION**

Written on the label.

#### **LIFE CYCLE**

This product should be disposed of when it loses its functionality.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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