# ST.330.280

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

### DESCRIPTION

The Universal Abutment Set is composed of Impression Coping, Provisional Coping (Polymer) and Castable Coping (Casting), both used to make prostheses on the Universal Abutment. Impression Coping: it is a prosthetic component made of polymer with the specific function of transfer (impression) on the Universal Abutment. It is indicated for the closed tray technique. The Provisional Coping: it is made of polymer. It is indicated for making provisional prostheses. The Castable Coping: it is made of castable polymer. It is indicated for making final metal-ceramic prostheses. The Universal Abutment Set is available in four measures, as per the table below, which must be used in accordance with the height and diameter of the Universal Abutment selected.

Universal Abutment Set	
Ø 3.3 X 4	Ø 3.3 X 6
Ø 4.5 X 4	Ø 4.5 X 6

### APPLICATIONS

The Impression Coping is a prosthetic component used to transfer the position of the Universal Abutment. The Provisional Coping is indicated for making provisional prostheses on the Universal Abutment. The Castable Coping is indicated for making a metallic infrastructure of cemented prostheses on the Universal Abutment.

### CONTRAINDICATIONS

This product does not have any contraindications, provided that it is used correctly for the purposes indicated.

### HANDLING

Impression Coping: Fit by means of a “click”

the Impression Coping selected on the Universal Abutment, to then execute the impression.

Provisional Coping: Position the Coping on the Universal Abutment Analog in the work model for indirect polymerization, or on the Abutment itself in the mount for direct polymerization, in the process of making the provisional prosthesis. After concluding the provisional prosthesis, it must be cemented in the mouth with provisional cement. Castable Coping: Position the Coping on the Universal Abutment Analog. Then, execute the polishing and casting. Thereby obtaining the metallic infrastructure of the prosthesis.

### PRESENTATION AND STERILIZATION

These products are packed in blister type packaging (surgical grade paper and film) supplied, in a single packaging. It is furnished sterile and single use. Sterilization method: ethylene oxide.

### PRECAUTIONS

- In the case of personalization of the Universal Abutment, this unit loses its purpose.
- Regarding the systemic aspect, consider the patient's general state of health as per the applicable literature. Regarding the local aspect, note the conditions of the intraoral tissues.
- This product is for single use and must not be re-sterilized (when applicable).
- Reuse of this product may cause:
  - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
  - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
  - Do not use the product if the packaging has been damaged.
  - Insertion torque exceeding that recommended and unsuitable Driver can cause damage to the material and render the system unusable.
  - The best results are obtained by using Neodent products. The use of instruments and/or Implants and/or prosthetic components of the systems does not ensure the perfect function and invalidates any product warranty.
  - Ensure that the parts are not swallowed or sucked in by the

person who is receiving the product (patient).

- Special care must be taken during the polishing so that the interocclusal space is not insufficient or excessive.
- Incorrect planning can jeopardize the performance of the Implant/Prosthesis unit leading to failures of the system, as loss or fracture of the Implant, loosening or fracture of the Prosthetic Screws.
- During the installation ensure that it is aligned with the insertion axis of the Implant, avoiding locking and damage to the system.
- Ensure that it is perfectly seated on the Implant. To do so, periapical X-rays are recommended with the parallelism technique.

### ADVERSE EFFECTS

Adverse effects do not occur, provided that the indications are respected.

### POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

### STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

### DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

### DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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