ST.330.281

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This device is intended for specialized procedure, which should only be executed by qualified dental professionals with specific training in CAD/CAM technique. For use in Brazil the use of the CAD/CAM Neodent Digital solution is recommended. For use in other countries the use of a CAD/CAM system using the Dental Wings (*Dental Wings Operating System - DWOS*) or CAD/CAM 3Shape system is recommended. To avoid unsatisfactory results for the patient, do not use the product without mastering the appropriate techniques. Always perform the correct procedures under conditions and with appropriate equipment.

DESCRIPTION

The Titanium Base is an intermediary prosthetic component between the implant and prosthesis manufactured in titanium alloy conforming to ASTM F136 standard. The Titanium Base has in one of its ends an interface for connection to the Implant. It can have different prosthetic interfaces on rotational (for bridge) and anti-rotational (for Crown/Coping) models. The height of the area to be cemented 4mm and must be used with no angulation.

Titanium Base (for Crown/Coping): It has cylindrical shape with outer retentions and three indexer guides to the Crown/Coping adjustment. It has different gingival heights (see table). The Titanium Base (for Crown/Coping) is available for Morse Taper, External Hex and Internal Hex Implants, prosthesis could be cemented or screwed.

Titanium Base (for bridge) - It has cylindrical shape with outer retentions. It has different gingival heights (see table). The Titanium Base (for bridge) is available for External Hex and Internal Hex Implants and it is indicated only for multiple cemented prostheses confection.

Prosthetic Interface	Gingival Height
Morse Taper	0,8; 1,5; 2,5; 3,5; 4,5
External Hex and Internal Hex	0,5; 2 ; 3 ;4

INDICATIONS FOR USE

Titanium Base Abutment is a titanium base placed onto Neodent dental implants to provide support for customized prosthetic restorations such as Copings, Crowns or Bridges.

It is indicated for cement-retained single or multi-unit restorations, or screw-retained single restorations.

APPLICATIONS

This product is indicated according to the available interocclusal space, existing gingival height and three-dimensional Implant position. It is indicated for cement-retained single or multi-unit restorations, or screw-retained single restorations in aesthetic areas installed in the jawbone or jaw.

It is recommended the use of CAD/CAM system using the Dental Wings (*Dental Wings Operating System - DWOS*) or CAD/CAM 3Shape system.

WARNING

The Titanium base abutment must not be modified, and must be used straight with no angulation.

Implants in combination with Titanium Bases with angled prosthetic structures are recommended in accordance with the table below:

Ø Implant (mm)	Region of the Mouth for installation	Zirconia	CoCr	IPS e.max
3.3	Mandibula 1-2, Maxilla 2	•	•	•
	Maxilla 1		•	
≥3.5	1-2	•	•	•
≥3.75	1-5	•	•	•
≥4.0	1-8	•	•	

Implants with a diameter equal to or less than 3.5 mm in combination with angled structures are recommended for the incisor region, with the exception of HE implants with 3.3 mm diameter in combination with angled structures made of zirconia or IPS e.max CAD, which are not recommended for the upper central incisor region.

Implants with 3.75 mm diameters, in combination with angled structures, are recommended only for the regions of incisors, canines, and pre-molars. Implants with diameters equal to or greater than 4.0 mm in combination with angled structures are recommended for all the regions of the mouth, with the exception of angled structures made of IPS e.max CAD, which are not recommended for the molar region.

CONTRAINDICATIONS

The Titanium Base abutment is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: Titanium alloy Ti6Al4V-ELI.

HANDLING

COPING/CROWN/BRIDGE CONFECTION BY CAD/CAM TECHNIQUE: Use the CAD software, according to the manufacturer's instructions, when making the Coping, Crown or Bridge. Always finish the Coping/Crown/Bridge before the Titanium Base fixation. When you use the digital platform select the Neodent Implant library - Titanium Base with an appropriate software platform, this way facilitating the interface drawing between the Coping/Crown/Bridge and Titanium Base. This library consists of a set of 3D files containing the milling array required for the Coping/Crown/Bridge geometry. The minimum thickness of the structure's wall varies according to the material chosen, as follows:

Material	Minimum thickness (mm)
Cobalt-Chromium	0,3
Zirconia	0,5
IPS e.max® CAD	0,6

COPING/CROWN/BRIDGE CONFECTION BY TRADITIONAL TECHNIQUE (E.G. PRESSING OR CASTING): Always finish the Coping/Crown/Bridge before the Titanium Base fixation. Place the Titanium Base on the analog and perform the waxing on the specific castable Coping for the selected Titanium Base. Perform the casting or pressing and finalize the Coping/Crown/Bridge. The Titanium must be used, in laboratory procedures, as follows: Place the Base component over the plaster model and position the Fixation Screw on its, applying a soft tight using the Screwdriver indicated (see chart). Protect access to the screw (with teflon and resinous compound) throughout the cementation process. The cement handling must follow the manufacturer instructions. Apply the cement in the outer portion of Titanium Base and press the fitting Coping/Crown/Bridge according to the three available indexer guides. The cemented Coping/Crown/Bridge must be pressed over the Titanium Base and the excess cement removed immediately. After the cement setting, the infrastructure is unscrewed from the Analog and the surplus cement remaining the Titanium Base edge is removed. NOTE: Is recommended to use Panavia - Kuraray cement. When using Lithium disilicate it is mandatory to use (IVOCLEAR Multilink) cement. After these procedures, must be handled using appropriate Drivers and adequate torque, as chart below:

Prosthetic Screw model that goes along the Titanium Base	Insertion Torque (N.cm)	Driver / Connection (mm)
Hex Screw	External Hex: 32	Hex 1.2
	II Plus: 20	
CM Screw	15	Hex 0.9

For more information about the preparation of the Coping/Crown/Bridge, please check:

- In Brazil: the guidebook "IMPLANT PREPARATION GUIDELINES", available on the Neodent website www.neodent.com.br/neodent-digital
- other countries: the preparation guidance from the specialized laboratory.

TRACEABILITY LABEL

This product is accompanied by three labels that allow its traceability and should be attached to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

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SANITATION

Before installing in mouth, you must sanitize and sterilize the final restoration. For Titanium bases in the Internal Hexagon and External Hexagon interfaces, disengage the screw before the sterilization. 1st step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer) 2nd step: Wash in ultrasonic cleaner for approximately 10 to 15 minutes. 3rd step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is recommended. 4th step: Dry with a clean and dry cloth or with compressed air. 5th step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent - 1st step - and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence.

PRESENTATION AND STERILIZATION

It is delivered sterile and It is of single use. Sterilization method: ethylene oxide. This product is unitarily provided, conditioned in blister (surgical grade film and paper type) packing. The following sterilization methods are recommended: moist heat (steam) autoclave, gravity-displacement or dynamic-air-removal (fractionated vacuum) cycle, unwrapped, 3 minute exposure at 132 °C (270 °F). The product must be unwrapped on an appropriate tray. Use the sterilized restoration immediately after sterilization, do not store.

PRECAUTIONS

- The Titanium Base must not be modified.
- For Titanium bases in the Internal Hexagon and External Hexagon interfaces, disengage the screw before the sterilization.
- Regarding the systemic aspect, consider the patient's general health conditions according to the applicable literature. Regarding the local aspect, observe the conditions of the intraoral tissues.
- This product is to be used only once, cannot be reused.
- Reuse of this product may cause:
 - adverse biological effects of residual products, microorganisms and / or

substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products. • Do not use the product if the package is damaged. • Use material and a technique suitable for the impression of the implant. • Make sure that the transfer used is compatible with the analog. • An insertion torque higher than the recommended and improper Drivers may cause damages to the material and make the system unusable. • Better results are achieved with the utilization of Neodent products. The utilization of instruments and/or Implants and/or Prosthetic Components of other systems do not assure the perfect function and voids any product guarantee. • Make sure that the parts are not swallowed or inhaled by the person who is receiving the product (patient). • The correct selection of the analog must be observed, according to the implant. • Local aspects of the patient must be taken into consideration in the selection of the prosthesis. • A special precaution should be taken during the closure so that the interocclusal space is not insufficient or excessive. • Make the occlusal and interproximal adjustment after the installation of the prosthetic. • Lack of adjustment or unsuitable occlusal adjustment may compromise the implant/prosthetic joint. • The improper planning may endanger the Implant/prosthesis set performance resulting in system faults, such as Implant loss or fracturing, loosening or fracturing of the Prosthetic Screws.

ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

Multilink® is a registered trademark of Ivoclar Vivadent AG, Liechtenstein.
Panavia™ is a trademark of Kuraray Co, LTD, JP.

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This device is intended for specialized procedure, which should only be executed by qualified dental professionals with specific training in CAD/CAM technique. The use of a CAD/CAM system using the Dental Wings (*Dental Wings Operating System - DWOS*) is recommended. To avoid unsatisfactory results for the patient, do not use the product without mastering the appropriate techniques. Always perform the correct procedures under conditions and with appropriate equipment.

DESCRIPTION

The Titanium Base is an intermediary prosthetic component between the implant and prosthesis manufactured in titanium alloy conforming to ASTM F136 standard. The Titanium Base has an appropriate interface for connection to the CM implant and is available in anti-rotational (for Crown/Coping) model. All Titanium Base abutments are provided with a prosthetic height of 4 mm and must be used straight with no customization or angulation. It has cylindrical shape with outer retentions and three indexer guides to the Coping adjustment. It is provided in prosthetic platform diameters of 3.5 mm or 4.5 mm, in gingival heights of 0.8, 1.5, 2.5, 3.5, and 4.5 mm. The Titanium Base for Morse Taper could be cemented or screwed.

INDICATIONS FOR USE

Titanium Base Abutment is a titanium base placed onto Neodent dental implants to provide support for customized prosthetic restorations. It is used with a Coping and Crown, or Crown alone, and is indicated for cement-retained or screw-retained single restorations.

APPLICATIONS

This product is indicated according to the available interocclusal space, existing gingival height and three-dimensional Implant position. It is indicated for cement-retained or screw-retained single restorations in aesthetic areas installed in the jawbone or jaw. It is recommended the use of CAD/CAM system using the Dental Wings (*Dental Wings Operating System - DWOS*).

WARNING

The Titanium base abutment must not be modified, and must be used straight with no angulation. Small diameter Implants (3.5 mm), in combination with angled structures, are recommended only for incisor regions.

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HANDLING

COPING/CROWN CONFECTION BY CAD/CAM TECHNIQUE: Use the CAD software, according to the manufacturer’s instructions, when making the Coping/Crown. Always finish the Coping or Crown before the Titanium Base fixation. When you use the digital platform select the Neodent Implant library - Titanium Base with an appropriate software platform, this way facilitating the interface drawing between the Coping/Crown and

Titanium Base. This library consists of a set of 3D files containing the milling array required for the Coping geometry. Note: For ceramic materials the structure of wall thickness cannot be less than 0,5 mm and the maximum angle to be no greater than 30°. The taper of the zirconia superstructure should not exceed 6°.

COPING/CROWN CONFECTION BY TRADITIONAL TECHNIQUE (E.G. PRESSING OR CASTING): Always finish the Coping/Crown before the Titanium Base fixation. Place the Titanium Base on the analog and perform the waxing on the specific castable Coping for the selected Titanium Base. Perform the casting or pressing and finalize the Coping/Crown.

The Titanium must be used, in laboratory procedures, as follows: Place the Base component over the plaster model and position the Fixation Screw on its, applying a soft tight using the Screwdriver indicated (see chart). Protect access to the screw (with teflon and resinous compound) throughout the cementation process. The cement handling must follow the manufacturer instructions. Apply the cement in the outer portion of Titanium Base and press the fitting Coping/Crown according to the three available indexer guides. The cemented Crown/Coping must be pressed over the Titanium Base and the excess cement removed immediately. After the cement setting, the infrastructure is unscrewed from the Analog and the surplus cement remaining the Titanium Base edge is removed. NOTE: Is recommended to use Panavia - Kuraray cement. When using Lithium disilicate it is mandatory to use (IVOCCLAR Multilink) cement. After these procedures, must be handled using appropriate Drivers and adequate torque, as chart below:

Prosthetic Screw model that goes along the Titanium Base	Insertion Torque (N.cm)	Driver / Connection (mm)
CM Screw	15	Hex 0.9

For more information on the preparation of the Crown/Coping follow the preparation guidelines of the specialized laboratory.

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SANITATION

Before installing in mouth, you must sanitize and sterilize the final restoration. 1st step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer) 2nd step: Wash in ultrasonic cleaner for approximately 10 to 15 minutes. 3rd step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is recommended. 4th step: Dry with a clean and dry cloth or with compressed air. 5th step: Perform a visual inspection, observing if there are failures in

the cleaning process. If there are still wastes, the part should be once again immersed in detergent - 1st step - and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence.

PRESENTATION AND STERILIZATION

It is delivered sterile and it is of single use. Sterilization method: ethylene oxide. This product is unitarily provided, conditioned in blister (surgical grade film and paper type) packing. The final restoration must be sterilized before installing in mouth. The following sterilization methods are recommended: moist heat (steam) autoclave, gravity-displacement or dynamic-air-removal (fractionated vacuum) cycle, unwrapped, 3 minute exposure at 132 °C (270 °F). The product must be unwrapped on an appropriate tray. Use the sterilized restoration immediately after sterilization, do not store.

PRECAUTIONS

- The Titanium Base must not be modified.
- Regarding the systemic aspect, consider the patient’s general health conditions according to the applicable literature. Regarding the local aspect, observe the conditions of the intraoral tissues.
- This product is to be used only once, cannot be reused.
- Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if the package is damaged.
- Use material and a technique suitable for the impression of the implant.
- Make sure that the transfer used is compatible with the analog.
- An insertion torque higher than the recommended and improper drivers may cause damages to the material and make the system unusable.
- Better results are achieved with the utilization of Neodent products. The utilization of instruments and/or Implants and/or prosthetic components of other systems do not assure the perfect function and voids any product guarantee.
- Make sure that the parts are not swallowed or inhaled by the person who is receiving the product (patient).
- The correct selection of the analog must be observed, according to the implant.
- Local aspects of the patient must be taken into consideration in the selection of the prosthesis.
- A special precaution should be taken during the closure so that the interocclusal space is not insufficient or excessive.
- Make the occlusal and interproximal adjustment after the installation of the prosthetic.
- Lack of adjustment or unsuitable occlusal adjustment may compromise the implant/prosthetic joint.
- The improper planning may endanger the Implant/prosthesis set performance resulting in system faults, such as Implant loss or fracturing, loosening or fracturing of the Prosthetic Screws.

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**MAGNETIC RESONANCE IMAGING (MRI)
- SAFETY INFORMATION**

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Instruct the patient as to the need of a professional medical monitoring after the

surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40° C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the

environment, it is recommended to take a look at the current legislation and adhere to it.

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