# ST.330.284

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This device, in its initial configuration, is intended exclusively for specialized laboratory procedures, which must be executed only by skilled professionals with specific training. For use in Brazil the use of the CAD/CAM Neodent Digital solution is recommended. For use in other countries the use of a CAD/CAM Dental Wings (Dental Wings Operating System - DWOS) system or CAD/CAM 3Shape system is recommended. To avoid unsatisfactory results for the patient, do not use the product without mastering the appropriate techniques. Always perform the correct procedures under conditions and with appropriate equipment.

**DESCRIPTION**

Preface is made of titanium alloy, according to ASTM F136, supplied in block, has a cylindrical shape and connection (interface) antirrotacional for connection with the implant. It is available for prosthetic interface Cone Morse interfaces and External Hexagon 4.1 / 4.3, in diameters 11.5mm and 15.8mm. After planning and customization, performed in a CAD/CAM software, it becomes an abutment to be used into Neodent implant and support the prosthesis (crown / coping). The CM Preface has an attached Screw. Preface with overall (milling) diameter of 15.8 mm is only used to mill an angled abutment.

**INDICATIONS FOR USE**

Preface Abutment is a titanium abutment to be used in fabricating a full custom abutment and placed onto Neodent dental implants to provide support for customized prosthetic restorations. Preface Abutments are indicated for screw-retained single restorations or cement-retained single or multi-unit restorations.

**APPLICATIONS**

This product is indicated according to the available interocclusal space, existing transmucosal height and three-dimensional Implant position. It is indicated for single restorations with screws or single or multiple cemented restorations on Implants installed in the jawbone or jaw. Preface must be designed and milled using CAD/CAM Dental Wings (Dental Wings Operating System - DWOS) system or CAD/CAM 3Shape system.

**WARNING**

Implants with diameters of 3.3 mm and 3.5 mm in combination with customized components with angling greater than 17° are recommended only for replacement of incisors. Implants with diameters of 3.75 mm in combination with customized components with angling greater than 17° are recommended only for replacement of incisors, canines, and pre-molars.

**CONTRAINDICATIONS**

The Preface abutment is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: Titanium alloy Ti6Al4V-ELI.

**HANDLING**

**PLANNING / ABUTMENT DESIGN:** Using the CAD/CAM Dental Wings (Dental Wings Operating System - DWOS) system, select the prosthetic interface and proceed the Abutment design according the plaster model anatomy. For the Abutment drawing some rules must be followed: 1 - The

gingival emergence profile must be prepared using the Healing or Provisional Abutments before receiving the final Abutment from the Preface. 2 - Position the emergence profile Abutment edge on vestibular region, slightly subgingival and follow the gingival contour on palatine/lingual region. 3 - The Abutment thickness should not be less than 0.4mm. 4 - The emergence profile diameter must be at least 3.3mm, to CM interface. To HE interface for 4.1 / 4.3 the minimum emergence profile diameter must be 4.8mm. 5 - The maximum limit for angulation is 30° and the maximum gingival height (transmucosal) until 3.5mm from the Implant platform. 6 - In case of Angled Abutment, the cementable height from the emergence profile (prosthetic height) should not exceed 6mm for Morse Taper interface and 7.7mm for External Hex interface. 7 - The total height of the abutment (above the implant-abutment connection platform) must be no greater than 14.2 mm and no shorter than 4 mm. 8- The range for the abutment body diameter (above the implant-abutment connection platform) is 3.3 mm to 8.5 mm.

**MANUFACTURE ABUTMENT:** After finish the design, fix the Preface on a compatible device with Preface and milling machine, following the manufacturer’s guidelines. To produce the designed Abutment you should use a CAM software to determine the wear strategy according manufacturer’s guidelines.

**PROCESSING:** After Preface milling it becomes the Abutment previously designed. This should be removed from the rest of the block with the aid of a cutting disc. Be careful not to change the dimensions of the component. The Preface interface must remain intact, cannot be modified. If necessary, the finishing may be done by specific drill to metal wear.

**INSTALLATION/CEMENTATION:** A clinical proof of the structure and complete final prosthesis is indicated. Perform the necessary adjustments. For cementation in the mouth, the instructions of the cement’s manufacturer must be followed. In this case, it is worthwhile to remember that is it necessary to know the technique for avoiding excess cement and chemical burns to the tissues adjacent to the Implant. Protect access to the screw throughout the cementation process. For the External Hex Abutment fixation you should use the Hex Screw according prosthetic platform (see the chart below).

Prosthetic Platform	Neotorque Screw	Titanium Screw
4.1 e 4.3	Neotorque Screw 4.1 and 4.3	Screw 4.1 and 4.3

Note: The CM Preface comes along with a Screw. The recommended torque for Abutment installation must follow the chart below:

Prosthetic Component	Torque Connection (mm)	Recommended Torque (N.cm)
External Hex Preface	1.2	32
CM Preface	0.9	15

For more information about the preparation of the crown/coping, please check:  
 • In Brazil: the guidebook “IMPLANT PREPARATION GUIDELINES”, available on the Neodent website [www.neodent.com.br/neodent-digital](http://www.neodent.com.br/neodent-digital))

• other countries: the preparation guidance from the specialized laboratory.

**TRACEABILITY LABEL**

This product is accompanied by three labels that allow its traceability and should be attached to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

**SANITATION**

Before installing in mouth, you must sanitize and sterilize the final restoration. 1<sup>st</sup> step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer) 2<sup>nd</sup> step: Wash in ultrasonic cleaner for approximately 10 to 15 minutes. 3<sup>rd</sup> step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is recommended. 4<sup>th</sup> step: Dry with a clean and dry cloth or with compressed air. 5<sup>th</sup> step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent - 1<sup>st</sup> step - and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence.

**PRESENTATION AND STERILIZATION**

The External Hex Preface is unitarily provided. The CM Preface comes along with a Screw. They are conditioned in blister (surgical grade film and paper type) packing. They are delivered non-sterile and they are of single use. The final restoration must be sterilized before installing in mouth. The following sterilization methods are recommended: moist heat (steam) autoclave, gravity-displacement or dynamic-air-removal (fractionated vacuum) cycle, unwrapped, 3 minute exposure at 132 °C (270 °F). The product must be unwrapped on an appropriate tray. Use the sterilized restoration immediately after sterilization, do not store.

**PRECAUTIONS**

- The Preface prosthetic interface does not allow customization.
- This product is to be used only once, cannot be reprocessed and reused.
- Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Use material and a technique suitable for the impression of the implant.
- Make sure that the transfer used is compatible with the analog.
- Correct selection of the analog must be observed, according to the implant.
- Make sure that the diameter of the block is compatible with the size of the abutment.
- Beware with cement excess during the Preface cementation in mouth.
- Be sure to select the Preface compatible with the implant.
- Selection of the prosthesis structure should consider local aspects of the patient.

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- During the installation, make sure to align the abutment to the insertion axis of the implant.
- Make sure that it is perfectly seated on the implant. To do so, we recommend using periapical X-rays with the parallel technique.
- Be sure to use a prosthetic screw compatible with prosthetic interface and prosthetic component.
- If the design, processing, and installation of the abutment, from the Preface is not done by the same dental professional, this Instruction for Use must be forwarded to the dentist responsible for the installation of the abutment.
- Inappropriate planning can jeopardize the performance of the implant/prosthesis unit leading to failures of the system, as loss or fracturing of the implant, loosening or fracturing of the prosthetic screws.
- Do not use the product if the packaging is damaged.
- Do not use the product after expiry date.
- For immediate load application, check the indication of the implant installed.
- See the torque to be given on the prosthetic component to be used. The excess or lack of torque can lead to undesirable result.
- Insertion torque greater than that recommended can make the system inoperative.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Make the occlusal and interproximal adjustment after the installation of the prosthetic.
- The lack of adjustment or inappropriate occlusal adjustment can endanger the implant/prosthesis unit.
- Some of the possible causes of unsuccessful

bone integration and loss of the prosthesis during the treatment are the following: unsuitable osteotomy, infections, deficient oral hygiene, occlusal trauma, systemic problems or diseases, low remaining bone quantity or quality, lack or failure of irrigation, use of instruments which are not specific and/or without cutting power and lack of specific training.

- Before each procedure check the conditions of the **Neodent** surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the **Neodent** product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the **Neodent** Implant System and exempts any product warranty.
- It is the professional's responsibility to use the **Neodent** products according to the instructions for use, and to determine whether it suits the individual situation of each patient.

#### ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

#### MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The **Neodent** Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the **Neodent** Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

#### POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

#### STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

#### DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

#### DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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This device, in its initial configuration, is intended for specialized lab procedure, which should only be executed by qualified dental surgeons with specific training in CAD/CAM technique. The use of a CAD/CAM Dental Wings (Dental Wings Operating System - DWOS) system is recommended. To avoid unsatisfactory results for the patient, do not use the product without mastering the appropriate techniques. Always perform the correct procedures under conditions and with appropriate equipment.

## DESCRIPTION

Preface is made of titanium alloy, according to ASTM F136, supplied in block, has a cylindrical shape and connection (interface) for crown/coping. It is available for Cone Morse interfaces, in diameters 11.5mm and 15.8mm. After planning and customization, performed in a CAD/CAM software, it becomes an abutment to be used into Neodent implant and support the prosthesis (crown / coping).

The CM Preface has an attached Screw. Preface with overall (milling) diameter of 15.8 mm is only used to mill an angled abutment.

## INDICATIONS FOR USE

Preface Abutment is a titanium abutment to be used in fabricating a full custom abutment and placed onto Neodent dental implants to provide support for customized prosthetic restorations. Preface Abutments are indicated for screw-retained or cement-retained single restorations.

## APPLICATIONS

This product is indicated according to the available interocclusal space, existing transmucosal height and three-dimensional Implant position. It is indicated for screw-retained or cement-retained single restorations. Preface must be designed and milled using a CAD/CAM Dental Wings (Dental Wings Operating System - DWOS) system.

## WARNING

Implants with diameter 3.5 mm in combination with components with angulation greater than 17° are recommended only for incisor teeth replacement.

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## HANDLING

**PLANNING / ABUTMENT DESIGN:** Using the CAD/CAM Dental Wings (Dental Wings Operating System - DWOS) system, select the prosthetic interface and proceed with the Abutment design according to the plaster model anatomy. For the Abutment design the following rules must be followed: 1 - The gingival emergence profile must be prepared using the Healing or Provisional Abutments before receiving the final Abutment from the Preface. 2 - Position the emergence profile Abutment edge on vestibular region, slightly subgingival and follow the gingival contour on palatine/lingual region. 3 - The Abutment thickness should not be less than 0.4mm. 4 - The emergence profile diameter must be at least 3.3mm. 5 - The maximum limit for angulation is 30° and the maximum gingival height (transmucosal) until 3.5mm from the Implant platform. 6 - In case of Angled Abutment, the

cementable height from the emergence profile (prosthetic height) should not exceed 6mm for Morse Taper interface. 7 - The total height of the abutment (above the implant-abutment connection platform) must be no greater than 14.2 mm and no shorter than 4 mm. 8 - The range for the abutment body diameter (above the implant-abutment connection platform) is 3.3 mm to 8.5 mm.

**MANUFACTURE ABUTMENT:** After finish the design, fix the Preface on a compatible device with Preface and milling machine, following the manufacturer's guidelines. To produce the designed abutment you should use a CAD/CAM Dental Wings (Dental Wings Operating System - DWOS) system to determine the wear strategy according manufacturer's guidelines.

**PROCESSING:** After Preface milling it becomes the Abutment previously designed. This should be removed from the rest of the block with the aid of a cutting disc. Be careful not to change the design of the component. The Preface interface must remain intact, cannot be modified. If necessary, the finishing may be done by specific drill to metal wear.

**INSTALLATION/CEMENTATION:** A clinical proof of the structure and complete final prosthesis is indicated. Perform the necessary adjustments. For cementation in the mouth, the instructions of the manufacturer must be followed. In this case, it is worthwhile to remember that it is necessary to know the technique for avoiding excess cement and chemical burns to the tissues adjacent to the Implant. Protect access to the screw throughout the cementation process. The Morse taper Preface is provided with an abutment screw. The recommended torque for Abutment installation must follow the chart below:

Prosthetic Component	Torque Connection (mm)	Recommended Torque (N.cm)
CM Preface	0.9	15

For more information on the preparation of the crown/coping follow the preparation guidelines of the specialized laboratory.

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## PRESENTATION AND STERILIZATION

The CM Preface comes along with a Screw. It is conditioned in blister (surgical grade film and paper type) packing. It is delivered non-sterile and they are of single use. The final restoration must be sterilized before installing in mouth. The following sterilization methods are recommended: moist heat (steam) autoclave, gravity-displacement or dynamic-air-removal (fractionated vacuum) cycle, unwrapped, 3 minute exposure at 132 °C (270 °F). The product must be unwrapped on an appropriate tray. Use the sterilized restoration immediately after sterilization, do not store.

## PRECAUTIONS

- The Preface prosthetic interface does not allow customization.
- This product is to be used only once, cannot be reprocessed and reused.
- Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Use material and a technique suitable for the impression of the implant.
- Make sure that the transfer used is compatible with the analog.
- Correct selection of the analog must be observed, according to the implant.
- Make sure that the diameter of the block is compatible with the size of the abutment.
- Beware with cement excess during the Preface cementation in mouth.
- Selection of the prosthesis structure should consider local aspects of the patient.
- During the installation, make sure to align the abutment to the insertion axis of the implant.
- Make sure that it is perfectly seated on the implant. To do so, we recommend using periapical X-rays with the parallel technique.
- Be sure to use a prosthetic screw compatible with prosthetic interface and prosthetic component.
- If the design, processing, and installation of the abutment, from the Preface is not done by the same dental professional, this Instruction for Use must be forwarded to the dentist responsible for the installation of the abutment.
- Inappropriate planning can jeopardize the performance of the implant/prosthesis unit leading to failures of the system, as loss or fracturing of the implant, loosening or fracturing of the prosthetic screws.
- Do not use the product if the packaging is damage.
- Do not use the product after expiry date.
- For immediate load application, check the indication of the implant installed.
- See the torque to be given on the prosthetic component to be used. The excess or lack of torque can lead to undesirable result.
- Insertion torque greater than that recommended can make the system inoperative.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.

- Make the occlusal and interproximal adjustment after the installation of the prosthetic.
- The lack of adjustment or inappropriate occlusal adjustment can endanger the implant/prosthesis unit.
- Some of the possible causes of unsuccessful bone integration and loss of the prosthesis during the treatment are the following: unsuitable osteotomy, infections, deficient oral hygiene, occlusal trauma, systemic problems or diseases, low remaining bone quantity or quality, lack or failure of irrigation, use of instruments which are not specific and/or without cutting power and lack of specific training.
- Before each procedure check the conditions of the **Neodent** surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the **Neodent** product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the **Neodent** Implant System and exempts any product warranty.
- It is the professional's responsibility to use the **Neodent** products according to the instructions

for use, and to determine whether it suits the individual situation of each patient.

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#### **MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION**

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#### **POST-OPERATIVE PRECAUTIONS AND MAINTENANCE**

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