ST.330.295

Português  Transfer CM e CM Exact

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

This product is a Prosthetic Component with the specific function of transfer (molding). The Transfer must be used on Morse Cone Implants. It is manufactured in titanium alloy according to standard ASTM F136 and comes in one or two pieces. The transfer is recommended for open tray (dragging) or closed tray (replacement) according to the table below:

Transfer Model	Rotational or anti-rotational	Recommended molding technique
CM Implant Impression Coping	Rotational	Closed tray
CM Implant Impression Coping (Short)	Rotational	Closed tray
CM Implant Open-Tray Impression Coping	Rotational	Open tray
CM Implant Exact Impression Coping	Anti-rotational	Closed tray
CM Exact Implant Open-Tray Impression Coping	Anti-rotational	Open tray

APPLICATIONS

This product is a prosthetic component used to transfer the position and shape of the implant Cone Morse in accordance with each planned case, respecting diameters and tolerance indications of the components, aimed at aesthetics and function. The transfer is recommended for unit restorations or multiple restorations according to the table below:

Model of the Transfer	Type of restoration
CM Implant Impression Coping	Unit restorations
CM Implant Impression Coping (Short)	Unit restorations
CM Implant Open-Tray Impression Coping	Multiple restorations
CM Implant Exact Impression Coping	Unit restorations
CM Exact Implant Open-Tray Impression Coping	Unit restorations

CONTRAINDICATIONS

This product has no contraindications, providing their intended purposes are correctly followed. The Open Tray technique is not recommended for patients with little buccal opening. The Transfer of the Exact CM Implant - Open Tray cannot be used for multiple restorations.

HANDLING

The Transfer must be selected according to the advance prosthetic planning, depending on the clinical situation of the patient and the interface of the implant. Fit the Neodent Transfer according to the planned molding technique. Open Tray: Fit the body of the Transfer together with the selected implant and turn only the screw with the help of 1.2 Torque Connection until the whole assembly is quite firm (do not exceed the maximum torque of 15N). For multiple restoration with CM Transfer, the transfers must be joined with acrylic resin. Make an opening in the tray to allow for the fitting of 1.2 Torque Connection after the molding. Prepare the molding material according to the recommendations of the manufacturer. Load the tray and introduce it into the patient's mouth. Wait the recommended time for setting and, using the 1.2 Torque Connection, completely loosen the screw of the transfer. Remove the molding and verify that the transfer went with the mold. Carefully fit the analog on the transfer. Make sure not to move the transfer during the fitting of the analog. Cast rubbery material to simulate the patient's soft tissue. Prepare the plaster model according to the appropriate techniques. Remove the transfer from the plaster model.

Closed Tray: screw the Transfer onto the Implant using the Transfer wrench until the whole assembly is very firm. Make sure that it is well-seated. Prepare the molding material according to the manufacturer's recommendations. Load the tray and introduce it into the patient's mouth. Wait the recommended setting time and remove the tray. With the Transfer wrench, remove the Transfer from the Implant, and then screw the Transfer onto the Analog. Carefully fit the Transfer-Analog assembly onto the mold. For the Exact CM transfer, make sure that the faceted surfaces of the Transfer (brim) are aligned with the faceted surfaces of the mold. Cast rubbery material to simulate the patient's soft tissue. Prepare the plaster mold according to the appropriate techniques. Remove the transfer from the plaster model.

PRESENTATION AND STERILIZATION

The Transfer and the screw are supplied individually, packaged in a blister pack (surgical grade paper and film). It is delivered sterile and It of single use. Sterilization method: ethylene oxide.

PRECAUTIONS

• As to the systemic aspect, consider the patient's general health condition according to the applicable literature. As to the local aspect, please observe the intra-oral tissues conditions. • This product is of single use and it may not be re-sterilized. • Reuse of this

product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products. • Do not use the product if the package is violated. • This product must be used immediately after opening the packaging, at the moment of surgery. If it is not used, discard it. • An insertion torque higher than that recommended and improper Drivers may cause damages to the material and make the system unusable. • Better results are achieved with the utilization of Neodent products. The utilization of instruments and/or Implants and/or Prosthetic Components of other systems do not assure the perfect function and voids any product guarantee. • Make sure that the parts are not swallowed or inhaled by the person who is receiving the product (patient). • The improper planning may endanger the Implant/prosthesis set performance resulting in system faults, such as Implant loss or fracturing, loosening or fracturing of the Prosthetic Screws. • During the installation, make sure to align it to the Implant insertion axle, avoiding jamming and damage to the threads. • Make sure that it is perfectly seated in the Implant. For that periapical radiographs are recommended with the parallelism technique.

ADVERSE EFFECTS

No adverse effects will occur, provided that the indications of this product are followed.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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