ST.330.296

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Español  Macho de Rosca Alvim

English  Alvim Bone Tap

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Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
	Não estéril / Non-sterile / No estéril / Non sterile / Nicht steril / Non stérile

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Consérvelo seco / Conservare in luogo asciutto / Trocken lagern / Garder au sec
	Manter afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conserver à l'abri de la lumière solaire
	Não utilizar se a embalagem estiver danificada / Do not use if package is damaged / No lo utilice si el envoltorio está dañado / Non utilizzare se la confezione originale è danneggiata / Bei beschädigter Verpackung Produkt nicht verwenden / Ne pas utiliser si l'emballage a été abimé

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Marcação CE para comercialização na Comunidade Europeia / CE Mark for European Community market / Marca CE para comercialización en la Comunidad Europea / Marchio CE per la commercializzazione in Europa / CE Kennzeichen zum Vertrieb in der Europäischen Gemeinschaft / Marquage CE pour une commercialisation en Europe
	Notificação exigida pelo FDA para comercialização nos Estados Unidos / Notification required by FDA for United States market / Notificación exigida por el FDA para comercialización en los Estados Unidos / Notifica necessaria per la commercializzazione negli Stati Uniti / Vom FDA geforderte Kennzeichnung für den Vertrieb in den USA / Notification exigée par FDA pour une commercialisation aux États-Unis
	Representante na Comunidade Europeia / Representative in the European Community / Representante en la Comunidad Europea / Rappresentante nella Comunità europea / Vertreter in der Europäischen Gemeinschaft / Représentant dans la Communauté Européenne
	Fabricante / Manufacturer / Fabricante / Produttore / Hersteller / Fabricant

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

The Alvim Bone Tap is produced in titanium in compliance with standard ASTM F136. At one end it has a hexagon which allows the adaptation of the Torque Wrench Connection (105.001 Hex Connection - Torque Wrench (Short) and 105.018 Hex Connection - Torque Wrench (Long)). The active part of the Alvim Bone Tap has cutting edges for forming threads and laser markings referring to the length of the Implants Alvim.

APPLICATIONS


The Alvim Bone Tap is indicated for forming threads in the surgical alveolus before installing the Implants Alvim bone bed type I or II, according to the bone classification of Lekholm and Zarb (1985). Lekholm & Zarb's jawbone quality classification: • Bone Type I = Homogeneous compact bone. • Bone Type II = Compact bone with core of trabecular bone. • Bone Type III = Thin cortical bone layer surrounding dense trabecular bone of favorable strength. • Bone Type IV = Thin cortical bone layer surrounding low density trabecular bone.

CONTRAINDICATIONS

• The Alvim Bone Tap must not be used in bone types III and IV. • The Alvim Bone Tap must not be used during the installation of Implants which are not Implants Alvim. • Do not use the Alvim Bone Tap for installation of implants Drive. • The handling as well as all the material to be used during the procedure must be sterile.

HANDLING

For using the Alvim Bone Tap employ, first, the sequence of burs as per planning and implant chosen. For use with the Torque Wrench, take the Alvim Bone Tap with the Hex Connection. Insert the Alvim Bone Tap in the surgical alveolus, maintaining the piercing axis. Make insertion movements, in the clockwise direction, slightly pressing the Torque Wrench, with maximum torque of 60 Ncm. Insert the Alvim Bone Tap until the marking related to the implant chosen (see image below). It may be necessary to alternate the direction of insertion and removal a few times to form the thread.

To remove the Alvim Bone Tap from the surgical

alveolus invert the direction of the Torque Wrench to the counterclockwise direction and proceed with the complete removal of the Male Thread, as otherwise it may impair the shape of the thread.

SANITATION

This product should be correctly sanitized after each use. Proceed as follows: 1st step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer) 2nd step: Wash in ultrasonic cleaner for approximately 10 to 15 minutes. 3rd step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is recommended. 4th step: Dry with a clean and dry cloth or with compressed air. 5th step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent - 1st step - and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence. 6th step: Select the packing in accordance to the sterilization process. Use preferably surgical grade paper self-sealing packing with laminated film. ATTENTION: Do not use anticrust products, avoid using enzymatic liquid in concentrations above 10% and dry parts that still have cleaning solution wastes, because these procedures cause oxidation. The use of steel brushes is also inadvisable.

PRESENTATION AND STERILIZATION

It is conditioned in blister (surgical grade film and paper type) packing and provided unitarily. It is reusable, delivered non-sterile and should, therefore, be sterilized before use. Sterilize the product on the day before or on the day of the procedure. It is not recommended to autoclave the product in its original packaging. We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard BS EN ISO 17665-1:2006 (Sterilization of health care products. Moist heat Requirements for the development, validation and routine control of a sterilization process for medical devices). Suitable parameters: 121°C, 30 min, at 1 atm. Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Sterilization validity: 7 to 15 days, provided the product is conditioned in clean, dry and protected from sunlight environment.

PRECAUTIONS

• To perform the preparation of surgical alveolus according to Alvim implant instrumental sequence chosen and planning surgery. • Execute the preparation of the surgical alveolus in accordance with the instrument sequence of the implant Alvim chosen and planning of the surgery. • For use with the Contra-Angle, adjust the surgical motor used in the procedure to a maximum speed of 30 RPM. • Do not use the product if the package is violated. • Do not sterilize the product in its original package. • This product should be sterile. • Make sure the parts are properly coupled before use. • This product should be handled firmly, avoiding contact with

soft tissues near beyond the region of interest. • Ensure that the patient does not breathe in the parts. • Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear. • Some of the possible causes of unsuccessful bone integration and loss of the prosthesis during the treatment are the following: unsuitable osteotomy, infections, deficient oral hygiene, occlusal trauma, systemic problems or diseases, low remaining bone quantity or quality, lack or failure of irrigation, use of instruments which are not specific and/or without cutting power and lack of specific training. • Inappropriate planning can jeopardize the performance of the Implant/prosthesis unit leading to failures of the system, as loss or fracturing of the Implant, loosening or fracturing of the Prosthetic Screws. • It is the dentist's responsibility to use the Neodent products in compliance with the instruction for use, as well as establish if it is suited to the individual situation of each patient. • Better results are achieved with the utilization of Neodent products. The utilization of instruments and/or Implants and/or Prosthetic Abutments of other systems does not assure the perfect function and voids any product guarantee.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and the use of instrument are improper.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

LIFE CYCLE

This product should be disposed of when it loses its functionality.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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