ST.330.298

Português 

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
1	103.170	Initial Drill
2	103.171	Alvim Twist Drill
3	103.172	Alvim Drill
4	103.173	Alvim Drill
5	103.174	Alvim Drill
6	103.216	CM Pilot Drill
7	103.219	CM Pilot Drill
8	103.220	CM Pilot Drill
9	103.091	Drill Extension
10	105.075	CM Implant Driver - Contra-Angle
11	105.073	CM Implant Driver - Torque Wrench (Short)

12	105.074	Cm Implant Driver - Torque Wrench (Long)
13	105.115	Narrow Smart Implant Driver - Contra-Angle
14	105.116	Regular Smart Implant Driver - Contra-Angle
15	105.117	Wide Smart Implant Driver - Contra-Angle
16	105.001	Smart/Ws Implant Driver - Torque Wrench (Short)
17	105.018	Hex Connection - Torque Wrench (Long)
18	111.036	Alvim Bone Tap
19	111.037	Alvim Bone Tap
20	111.038	Alvim Bone Tap

21	103.366	Countersink Drill
22	103.211	Countersink Drill
23	103.242	Sf Countersink Drill
24	128.015	CM Height Measurer
25	104.041	Manual Screwdriver (Medium)
26	104.012	Manual Screwdriver (Medium)
27	128.019	Direction Indicator
28	128.022	Direction Indicator
29	128.023	Direction Indicator
30	104.050	Torque Wrench
	129.001	Titanium Tweezers

The instruments contained in the case are sold separately.

This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

These are boxes for the storage of surgical and/or prosthetic instruments. They are made in polymer that can be sterilized in autoclave. They are provided with silicon supports for conditioning and safe fastening of each tool. They have markings guiding the tools use during the procedures. The Case may be equipped according to the work procedure.

APPLICATIONS

It is suitable for safe conditioning of surgical and/or prosthetic instruments, orientation of the instruments use sequence, as well as a support for sterilization of such instruments.

CONTRAINDICATIONS

It is contraindicated for sterilization under dry heat (ovens).

HANDLING

- Remove the instrumentation access cover pressing lightly the locking system between the base and the cover. To close the Case, position the cover on the base and press on it.
- For sterilization, the disinfected, clean and dry instruments should be conditioned in the Case, which should be packaged in surgical grade paper for autoclave.
- During the sterilization, the Case should not lean on the autoclave walls.
- During surgery the case should be placed on a sterile surface and handled with sterile gloves.
- The colors of the placeholders correspond to the indication of the instruments to be used according to the diameter of the implant to be installed (see table below):

Color	Implant diameter
Purple	3.5
Yellow	4.3
Orange	5.0

CASE ASSEMBLY

For the assembly of this case, each placeholder is related to a number from the instrumental table, check the picture above. The titanium tweezers stay under the tray. The instruments contained in the case are sold separately and must be used according to their own Instructions.

SANITATION

The carrying cases should be properly sanitized after its utilization. Do not let the product contact the moisture beyond the time required for the sanitation. For such, you should proceed as follows: Step 1: Disassemble the Case. Step 2: Submerge the part fully into enzymatic detergent (diluted according to the manufacturer). Step 3: Wash in ultrasonic washer for approximately 10 to 15 minutes. Step 4: Rinse with plenty of distilled water until the solution wastes are fully removed. The utilization of nylon brushes is recommended. Step 5: Dry with a clean and dry cloth or with compressed air. Step 6: Perform a visual inspection, watching if there are faults in the cleaning process. If any wastes are still left, the part should be submerged again in detergent - step 1 - and if necessary, the cleaning should be

made with the help of a nylon brush. Repeat the rinsing and drying sequence. Step 7: Select the package according to the sterilization process. Use preferably surgical grade paper self-sealing package with laminated film. ATTENTION: Do not use descaling products and avoid drying parts still containing cleaning solution wastes because such procedures favor oxidation of the instruments in the case. The use of steel brushes is not recommended.

PRESENTATION AND STERILIZATION

This product is supplied unitarily, conditioned in cardboard packaging. It is reusable, delivered non-sterile. It must be sterilized after insertion of instruments and before use. Sterilize the product on the eve or in the procedure day. Preferably, it is recommended to follow the steam autoclave sterilization method and also, the parameters and procedures established in the standard BS EN ISO 17665-1:2006 (Sterilization of health care products. Moist heat. Requirements for the development, validation and routine control of a sterilization process for medical devices). Recommended parameters: 121 °C, 30 min, at 1 atm. Alternatively, you may proceed according to the autoclave manufacturer's instruction.

PRECAUTIONS

- This product should be used in the sterile condition.
- The products contained in the case must be cleaned and/or sterilized according to their own instructions.
- Do not sterilize under dry heat (oven).
- The Case should not lean on the Autoclave walls, to prevent it from melting.
- Place the water measure for

autoclave recommended by the manufacturer. Lack of water during the sterilization cycle may compromise the effectiveness of the autoclave, lead to melting of the case, as well as damage to the instruments.

- The contact of the Case with acrylic resin liquid results in damages of the Case polymer structure.
- The chemical sterilization is not recommended.
- Inappropriate planning can jeopardize the performance of the implant/prosthesis unit leading to failures of the system, as loss or fracturing of the implant, loosening or fracturing of the prosthetic screws.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Before each procedure check the conditions of the **Neodent** surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the **Neodent** product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the **Neodent** Implant System and exempts any product warranty.
- It is the dentist's responsibility to use the **Neodent** products in compliance with the instruction for use, as well as establish if it is suited to the individual situation of each patient.

ADVERSE EFFECTS

Adverse effects will only occur if the selection and use of the product are improper.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION


Written on the label.

LIFE CYCLE

This product should be disposed of when it loses its functionality.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

1	103.170	Fresa Lanza Plus
2	103.171	Fresa Helicoidal Alvim Plus
3	103.172	Fresa Alvim Plus
4	103.173	Fresa Alvim Plus
5	103.174	Fresa Alvim Plus
6	103.216	Fresa Piloto CM Plus
7	103.219	Fresa Piloto CM Plus
8	103.220	Fresa Piloto CM Plus
9	103.091	Prolongador de Fresas
10	105.075	Conexion CM para Contra-Angulo
11	105.073	Conexion CM Corta para Carraca

12	105.074	Conexion CM Larga para Carraca
13	105.115	Conexion para Contra-Angulo Smart Estrecho
14	105.116	Conexion para Contra-Angulo Smart Regular
15	105.117	Conexion para Contra-Angulo Smart Ancho
16	105.001	Conexion Hexagonal Corta para Carraca
17	105.018	Conexion Hexagonal Larga para Carraca
18	111.036	Macho de Rosca Alvim
19	111.037	Macho de Rosca Alvim
20	111.038	Macho de Rosca Alvim

21	103.366	Broca Countersink Plus
22	103.211	Fresa Countersink Plus
23	103.242	Fresa Countersink Sf Plus
24	128.015	Medidor de Altura CM
25	104.041	Llave Digital Media
26	104.012	Llave Digital Media
27	128.019	Paralelizador
28	128.022	Paralelizador
29	128.023	Paralelizador
30	104.050	Llave Carraca Torquimetro
	129.001	Pinza de Titanio

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