ST.330.299

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This device is intended for a specialized procedure, which should be performed by professionals qualified with specific training on the CAD/CAM technology, on the CEREC® Sirona system for digitalization and production of custom prostheses. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions.

DESCRIPTION

The Titanium Base for CEREC® is an intermediary prosthetic component made of titanium alloy according to ASTM standard F136, to be installed between the implant and the prosthesis, which is made using the CEREC® Sirona system for digitalization and production of custom prostheses.

The Titanium Base for CEREC® has a cylindrical shape, with indexing guide for fitting the Coping, height of cementable area of 4.7 mm, and has variations of transmucosal heights (see table). The geometry of the cementable area of the Titanium Base for CEREC® is compatible with the L blocks of zirconia or lithium disilicate available for CEREC®.

The Titanium Base for CEREC® is available for adaptation on Neodent Implants with Morse Taper, External Hex (prosthetic platform 3.3 and 4.1/4.3) and II Plus interfaces. On prosthetic interfaces for External Hex and II Plus Implants, they come with a screw. Titanium Bases for CEREC® with interface for Morse Taper Implants come with a screw already coupled.

| Prosthetic Interface | Transmucosal Height (mm) |
|-------------------------------|---------------------------------------|
| Morse Taper | 0.8; 1.5; 2.5; 3.5 ; 4.5; 5.5; 6.5 |
| External Hex and Internal Hex | 0.5; 2 ; 3 ;4 |

INDICATIONS FOR USE

The Titanium Base is a titanium component that is placed over Neodent Implants to provide support for custom prosthetic restorations, such as copings or crowns. It is indicated for single-tooth screw-mounted restorations.

APPLICATIONS

The Titanium Base for CEREC® is indicated according to the interocclusal space available, the existing transmucosal height, and the three-dimension position of the Implant. It is indicated for single-tooth restorations on screw-mounted prosthesis over dental implants. It is recommended to use the CAD/CAM technique through the CEREC® Sirona system.

WARNING

The Titanium Base for CEREC® is provided without angling and should not be modified. Implants in combination with Titanium Bases for CEREC® with angled structures are recommended according to the table below:

| Ø Implant (mm) | Region of the Mouth for installation | Zirconia | IPS e.max |
|----------------|--------------------------------------|----------|-----------|
| 3.3 | Mandibula 1-2, Maxilla 2 | ● | ● |
| | Maxilla 1 | ● | |
| ≥3.5 | 1-2 | ● | ● |
| ≥3.75 | 1-5 | ● | ● |
| ≥4.0 | 1-8 | ● | |

Implants with a diameter equal to or less than 3.5 mm in combination with angled structures are recommended for the incisor region, with the exception of HE implants with 3.3 mm diameter in combination with angled crowns/copings made of IPS e.max® CAD, which are not recommended for the upper central incisor region.

Implants with 3.75 mm diameters, in combination with angled crowns/copings, are recommended only for the regions of incisors, canines, and pre-molars.

Implants with diameters equal to or greater than 4.0 mm in combination with angled crowns/copings are recommended for all regions of the mouth, except when in combination with angled crowns/copings made of IPS e.max® CAD, which are not recommended for the molars regions.

CONTRAINDICATIONS

The Titanium Base for CEREC® does not allow customization, is contraindicated for cases of multiple teeth applications, and is not compatible with the Sirona Galileos® Software. This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: Titanium alloy Ti6Al4V.

HANDLING

FOR INTRAORAL SCANNING: Select and install the Titanium Base for CEREC® on the implant according to the preplanning. Insert the selected scanbody* over the Titanium Base for CEREC® compatible with the equipment of the CAC/CAM system used and perform the scan according to the manufacturer's instructions. Ensure the correct fit of the transfer on the Titanium Base for CEREC®.

FOR PLASTIC MODEL SCANNING: Perform the molding of the implant in accordance with the preplanning, the clinical situation of the patient and the implant interface. Use artificial gingiva on the plaster model to simulate the soft tissue of the patient. Create the plaster model according to the proper techniques. Select and install the Titanium Base for CEREC® on the analog. Insert the selected scanbody* compatible with the equipment of the CAC/CAM system used and perform the scan according to the manufacturer's instructions. Ensure the correct fit of the transfer on the Titanium Base for CEREC®.

DESIGN OF THE STRUCTURE OF THE PROsthESIS: CEREC® Titanium Bases for Neodent implants (CM, HE, and II Plus) can be used with any of the specific CEREC® software libraries in conjunction with Sirona products, according to the table below:

| Library | Sirona's Products | | | |
|---------------|-------------------|----------------------|------------------------------|-------------------|
| | Scanbody Sirona | REF Scanbody Omnicam | REF Scanbody Bluecam / Ineos | Grinding Block |
| Tibase CEREC® | | | | |
| NBB 3.4 L | L | 6431329 | 6431303 | inCoris ZI meso L |
| NB A 4.5 L | L | 6431329 | 6431303 | inCoris ZI meso L |
| SSO 3.5 L | L | 6431329 | 6431303 | inCoris ZI meso L |
| S BL 3.3 L | L | 6431329 | 6431303 | inCoris ZI meso L |
| S BL 4.1 L | L | 6431329 | 6431303 | inCoris ZI meso L |
| BO 3.4 L | L | 6431329 | 6431303 | inCoris ZI meso L |

Design the structure of the prosthesis on the software according to the prosthetic planning and the material to be used. • Design the external form of the meso-structure according to the preparation guidelines for the superstructure required. • Keep a minimum thickness of 0.5 mm circularly around the wall of the screw channel. Be sure not to exceed an angle of 20° between the axis of the implant and the axis of the restoration. On restorations designed for the HE 3.3 prosthetic interface, do not exceed 17°. • If the meso-structure is designed to receive aesthetic ceramic, make sure that this will not narrow the channel for the screw. The cavity for fitting the meso-structure onto the Titanium Base for CEREC® must not be coated. • Be sure that there are not sharp edges or corners in the design of the meso-structure. Before creating the structure of the prosthesis, make sure that the size of the machining block is compatible with the design being undertaken. Create the piece as indicated.

CLINICAL PROOF OF THE STRUCTURE OF THE PROSTHESIS: Perform a clinical proof of the created structure on the Titanium Base for CEREC®. Make sure of the proper fit between the pieces, clinically and using X-ray. Make sure there is occlusal and interproximal space and that the structure of the prosthesis still meets the aesthetic and functional requirements. Apply the aesthetic ceramic according to the prosthetic planning and the appropriate techniques.

CEMENTATION OF THE TITANIUM BASE FOR CEREC®: It is recommended to screw the Titanium Base for CEREC® onto an analog of the implant throughout the process of finishing and finalization of the prosthesis. Position the Titanium Base for CEREC® over the analog and put in the mounting screw over it, tightening it slightly with the Digital Driver indicated (see table below). The surfaces of the Titanium Base for CEREC® intended for cementation of the structure of the prosthesis must be blasted (aluminum oxide, 50 µm, maximum pressure of 2 bar). Protect the access of the screw (with teflon and resin compound) during the cementing process. NOTE: It is recommended to use chemically activated resin cements for bonding on metals (e.g. Panavia - Kuraray). Handling of the cement should follow the instructions of the manufacturer. Apply the cement to the external portion of the Titanium Base for CEREC® and press the Coping, fitting it according to the indexing guide. The cemented Coping should be pressed over the Titanium Base for CEREC® and the excess cement which overflows through its hole should be immediately removed. After the cement sets, the infrastructure is unscrewed from the Analog and the excess cement remaining on the edge of the Titanium Base for CEREC® is removed. When Lithium Disilicate is applied, the use of high-opacity cement is mandatory.

INSTALLATION OF THE CUSTOM PROSTHESIS: Sanitize and sterilize the assembly, Titanium Base for CEREC® and custom prosthesis, before putting it in the mouth. After these procedures, position the assembly in the mouth using the Connection and the torque indicated according to the table below. Make the necessary occlusal and interproximal adjustments.

| Model of Prosthetic Screw that comes with the Titanium Base for CEREC® | Insertion torque (N.cm) | Screwdriver / Connection |
|--|-------------------------|--------------------------|
| Hex Screw | External Hex: 32 | Hex 1.2 |
| | II Plus: 20 | |
| CM Screw | CM: 15 | Hex 0.9 |

*The scanbody compatible with the Titanium Base for CEREC® is an accessory that is not provided by Neodent.

TRACEABILITY LABEL

This product is accompanied by three labels that allow its traceability and should be attached to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

SANITATION

Before installing in mouth, you must sanitize and sterilize the final restoration. For Titanium bases for CEREC® in the External Hexagon and II Plus interfaces, disengage the screw before the sterilization.

1st step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer) 2nd step: Wash in ultrasonic cleaner for approximately 10 to 15 minutes. 3rd step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is recommended. 4th step: Dry with a clean and dry cloth or with compressed air. 5th step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent - 1st step - and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence.

PRESENTATION AND STERILIZATION

This product is indicated for single use and provided sterilized by the ethylene oxide method, packed individually into a package. The final prosthesis must be sterilized before installation in the mouth. The following sterilization methods are recommended: moist heat (steam) autoclave, gravity-displacement or dynamic-air-removal (fractionated vacuum) cycle, unwrapped, 3 minute exposure at 132 C (270 F). The product must be unwrapped on an appropriate tray. Use the sterilized restoration immediately after sterilization, do not store.

PRECAUTIONS

- Inform the laboratory that any modification or damage to the prosthetic interface requires that the piece be discarded and remade.
- Ensure use of components that have the original Neodent prosthetic interface, which must not be machined in the laboratory.
- Make sure that the manufacturer’s recommendations are shared between the dentist and the laboratory during the fabrication process of the prosthesis. The blocks are for single use and should not be reprocessed.
- When choosing the Titanium Base for CEREC®, make sure that it is compatible with the diameter of the implant interface installed.
- The Titanium Base for CEREC® is not indicated for multiple restorations, nor for cementation in

the mouth.

- Special care must be taken during design of the prosthesis so that the interocclusal space will not be insufficient or excessive.
- Surgical planning and/or inadequate prosthesis can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of components and/or prosthetic screws.
- The selection of the material of the structure of the prosthesis must take local aspects of the patient into account.
- This product is of single use and cannot be reprocessed.
- Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if the packaging is broken.
- The sterilization is guaranteed only if the packaging is not damaged.
- Do not use the product after expiry date.
- For immediate load application, check the torque indication of the implant installed.
- See the torque to be given on the prosthetic component to be used. The excess or lack of torque can lead to undesirable result.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- During installation, make sure to align it with the insertion axis of the implant. Make sure it is perfectly seated on the implant. To do so, it is recommended to use periapical X-rays with the parallel technique.
- Make sure to use a prosthetic screw that is compatible with the prosthetic interface and with the prosthetic component.
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/prosthesis assembly.
- Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the Neodent product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the professional’s responsibility to use the Neodent products according to the instructions for use.

ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

**MAGNETIC RESONANCE IMAGING (MRI)
- SAFETY INFORMATION**

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

**POST-OPERATIVE PRECAUTIONS AND
MAINTENANCE**

Instruct the patient as to the need of a professional medical monitoring after o

procedimento and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored, in its original packaging, in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

After use, all the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

The libraries cited are protected by copyright and belong to their respective authors. All rights are reserved by these authors. The libraries in question are available in the CEREC software from Sirona Dental Systems GmbH. All rights reserved. CEREC® is a registered trademark of Sirona Dental Systems GmbH. Galileos is a registered trademark of Sirona Dental Systems GmbH, Alemanha.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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