# ST.330.300.01

Português 

Macho de Rosca Neodent Guided Surgery Alvim e Guia de Macho de Rosca Neodent Guided Surgery Alvim

Italiano 

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants, with knowledge of the Neodent Guided Surgery technique. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

### DESCRIPTION

The Alvim Neodent Guided Surgery Bone Tap is made of titanium in compliance with standard ASTM F136. On one end, it has a hexagon that allows for fitting the Hex Connection for Ratchet and Contra-angle Handpiece. The active portion of the Alvim Neodent Guided Surgery Bone Tap has cutting edges for formation of threads and laser markings for the length of the Alvim Implants, using the Neodent Guided Surgery technique. The Alvim Neodent Guided Surgery Bone Tap Guide has a surgical stainless steel handle and titanium alloy tap guide, in accordance with standard ASTM F136, and is available in diameters of Ø 3.5 and Ø 4.3 mm. It has a cylindrical geometry and a Stop for reference use on the grommet of the Surgical Guide.

### INDICATIONS FOR USE

The Neodent Guided Surgery technique is indicated for surgical procedures on maxilla or mandible bones with a sequence of guided instruments, allowing for flapless surgery for installation of implants in edentulous areas.

### APPLICATIONS

The Alvim Neodent Guided Surgery Bone Tap is indicated for formation of threads in the surgical socket, using the Neodent Guided Surgery technique, before installation of Alvim Implants in type I and II bone bed, according to the Lekholm and Zarb bone classification (1985). Lekholm and Zarb bone classification (1985):

- Bone Type I = Homogeneous compact bone.
  - Bone Type II = Compact bone with trabecular bone core.
  - Bone Type III = Thin layer of cortical bone around dense trabecular bone with good resistance.
  - Bone Type IV = Thin layer of cortical bone around dense trabecular bone with low density.
- The Alvim Neodent Guided Surgery Bone Tap Guide is an instrument that serves to direct the Alvim Neodent Guided Surgery Bone Tap on the Surgical Guide during the formation of threads in the surgical socket. It should be used according to the table below:

Color	Diameter of the Neodent Guided Surgery Bone Tap (in mm)
Purple	Ø 3.5
Yellow	Ø 4.3

### CONTRAINDICATIONS


The Alvim Neodent Guided Surgery Bone Tap is contraindicated for bone types III and IV. The Alvim Neodent Guided Surgery Bone Tap is contraindicated for installation of implants other than CM Alvim Implants. Do not use the Neodent Guided Surgery Bone Tap for installation of CM Drive implants.

### HANDLING

To use the Alvim Neodent Guided Surgery Bone Tap, first use the sequence of drills according to

the planning and the implant selected. To use the Alvim Neodent Guided Surgery Bone Tap Ø 3.5 and Ø 4.3, select the corresponding Alvim Neodent Guided Surgery Bone Tap Guide. Make sure to use the set correctly. Fit the end of the Neodent Alvim Neodent Guided Surgery Tap Guide inside the grommet of the Surgical Guide, holding it firmly by the handling during instrumentation. The Alvim Neodent Guided Surgery Bone Tap Ø 3.5 Straight does not require use of the Bone Tap Guide.

For use with contra-angle handpieces, attach the Alvim Neodent Guided Surgery Bone Tap to the Hex Connection for Contra-angle Handpiece. Set the surgical motor to a rotation speed of no more than 30 rpm and maximum torque of 35 N.cm. Bring the set to the hole of the Alvim Neodent Guided Surgery Bone Tap Guide. Start the motor of and insert the Alvim Neodent Guided Surgery Bone Tap into the hole of the Alvim Neodent Guided Surgery Bone Tap Guide, respecting the drilling axis, until stability is attained. Proceed with the Torque Meter Ratchet Wrench, using the Hex Connection for Ratchet. Make clockwise insertion motions with light pressure on the Ratchet Wrench, with maximum torque of 60 N.cm. Insert the Alvim Neodent Guided Surgery Bone Tap until you reach the marking for the implant chosen (see image below). It may be necessary to alternate the direction between insertion and removal a few times to form the thread. To remove the Alvim Neodent Guided Surgery Bone Tap from the surgical socket, invert the direction of the Torque Meter Ratchet Wrench to counterclockwise and proceed with full removal of the Alvim Neodent Guided Surgery Bone Tap, otherwise it can damage the shape of the thread.

Alvim Neodent Guided Surgery Bone Tap

ATTENTION: The figure is merely illustrative. It does not represent the real dimension and characteristics of the product. The Alvim Neodent Guided Surgery Bone Tap has an additional length, compensating for the thickness of the soft tissue, the height of the grommet, and the height of the Neodent Alvim Neodent Guided Surgery Tap Guide, when applicable.

Nominal size (mm)	Real size (mm)	
	Neodent Alvim Neodent Guided Surgery Bone Tap - Regular	Neodent Alvim Neodent Guided Surgery Bone Tap - Straight
8	19	17
10	21	19

Nominal size (mm)	Real size (mm)	
	Neodent Alvim Neodent Guided Surgery Bone Tap - Regular	Neodent Alvim Neodent Guided Surgery Bone Tap - Straight
11.5	22.5	20.5
13	24	22
16	27	25

### SANITIZATION

These products must be correctly sanitized after each use. Do not leave it in contact with dampness for prolonged periods. Follow these steps for cleaning: 1st step: Remove drill from contra-angle (if applicable). 2nd step: Immerse drill completely in enzymatic detergent (10% solution). 3rd step: Wash it in an ultrasonic washer for approximately 10 minutes. 4th step: Rinse it with plenty of distilled water, until all solution residues have been completely removed. The use of nylon brushes is permitted. 5th step: Dry it with a paper towel or compressed air. 6th step: Perform a visual inspection, observing whether there are any faults in the cleaning process. If you still find any residues, reimmerse the drill in detergent - 2nd step - and clean it using a nylon brush, if necessary. Repeat rinsing and drying sequence. 7th step: If applicable, place the drill in an appropriate, properly cleaned and dried case. 8th step: Select packaging according to sterilization process. Preferably use self-sealing package of medical grade paper with laminated film. IMPORTANT: Never store drills without drying them thoroughly. Drying is very important to prevent oxidation. CAUTION: Do not use incrustation removal products, as they tend to darken drills and cause oxidation. Avoid using enzymatic liquid in concentrations above 10%, and drying drills that still contain residues from cleaning solutions, as these procedures favor oxidation. The use of steel brushes is also not recommended. The drying of the parts is of utmost importance before storage and sterilization, because the accumulation of moisture on the products is harmful and may cause oxidation.

### PRESENTATION AND STERILIZATION

These products are reusable and supplied non-sterile, packed individually in a package that provides double protection. These products must be correctly sanitized and sterilized before use. Sterilize the products the day before or on the day of the procedure. ATTENTION: It is not recommended to autoclave this product in its original packaging. Sterilize the product on the day before or on the day of the procedure. We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard BS EN ISO 17665-1:2006 (Sterilization of health care products. Moist heat Requirements for the development, validation and routine control of a sterilization process for medical devices). Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Suitable parameters: 121°C, 30 min, at 1 atm. Sterilization validity: 7 to 15 days, provided that the product is conditioned in clean, dry and protected from the sunlight environment.

## PRECAUTIONS

- Inappropriate planning can jeopardize the performance of the implant unit leading to failures of the system, as loss or fracturing of the implant.
- Some of the possible cause of potential failure in osseointegration and loss of the prosthesis during the treatment are: inadequate osteotomy, infections, deficient oral hygiene, occlusal trauma, diseases or systemic problems, low quantity or quality of remaining bone, lack or failure in irrigation, use of non-specific instruments and/or instruments without cutting power, and lack of specific training.
- The **Neodent** Guided Surgery technique does not allow for installation of implants with diameters equal to or greater than 5.0 mm.
- Be aware in cases patients that show signs of allergy or hypersensitivity to the chemical components of the material: titanium alloy and titanium Ti6Al4V-ELI.
- Do not use the product if the packaging is broken.
- Do not use the product with the validity expired.
- Before each procedure, make sure the pieces are properly seated.
- Be sure you have all the instruments needed to perform the procedure according to the surgical planning.
- Before each procedure check the conditions of the **Neodent** surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- During the surgical procedure, make sure you

don't exert lateral movements with the surgical instruments, avoiding the loss of the mounting of the Grommet.

- Perform the preparation of the surgical alveolus according to the instrument sequence of the Alvim implant chosen and surgical planning.
- This product must be handled securely, avoiding contact with nearby soft tissues beyond the area of interest.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Always use the product **Neodent** sequence for the **Neodent** Guided Surgery technique. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the **Neodent** Implant System and exempts any product warranty.
- It is the professional's responsibility to use the **Neodent** products according to the instructions for use.

## ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

## POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

## STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

## DISPOSAL OF MATERIAL

After use, every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

## DATE OF EXPIRATION

Written on the label.

## LIFE CYCLE

The **Neodent** Guided Surgery Bone Tap Guide must be discarded when it loses its functionality. The Alvim **Neodent** Guided Surgery Bone Tap has a useful lifespan of up to 30 instrumentations for bone qualities I and II, provided that its conditions for use are respected, including adequate irrigation and rotation. Regardless of the lifespan mentioned, the professional should always check the sharpness condition of the instrumental. Products with improper conditions for cutting should be discarded. Note: Classification of the bone quality according to Lekholm and Zarb (1985).

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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