ST.330.305

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants, with knowledge of the **Neodent Guided Surgery** technique. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

The **Neodent Guided Surgery Drill 1.3** is made of surgical stainless steel. On one end, it has a coupling for a Contra-Angle handpiece according to ISO 1797-1 - *Dental rotary instruments - Shank*. On its active end, it has a helical format, to perform bone drilling. The **Neodent Guided Surgery Clamp** is a cylindrical pin that has one end with a sharp active tip, for bone anchorage, and the other end is spherical, for handling. It is made of titanium alloy, according to ASTM standard F136, and it has a *Stop* for reference for use on the *Sleeve*.

INDICATIONS FOR USE

The **Neodent Guided Surgery** technique is indicated for surgical procedures on maxilla or mandible bones with a sequence of guided instruments, allowing for flapless surgery for installation of implants in edentulous areas.

APPLICATIONS

The **Neodent Guided Surgery Drill 1.3** is an instrument used for bone drilling, allowing for the insertion of the *Clamp* on the prototyped *Surgical Guide*. The *Clamp* serves as a mounting support for the *Surgical Guide* during the implant installation procedure.

CONTRAINDICATIONS

There are no contraindications, provided that it is used correctly for the purposes indicated.

HANDLING

Bring the *Surgical Guide* to the mouth and make sure it is seated in the correct position. Fit the **Neodent Guided Surgery Drill 1.3** on the contra-angle handpiece. Set the surgical motor with rotation speed of 500-800 rpm. Bring the 1.3 Drill to the *Sleeve* of the mounting pin. Activate the motor and perform the drilling with continuous movements of insertion and removal, with plenty of irrigation, until the *Stop* of the 1.3 Drill touches the *Sleeve* of the mounting pin. The drilling direction must follow the direction of the *Sleeve*, without lateral movements. The irrigation can be manual or combined with the irrigation from the motor. After drilling, manually insert the *Clamp* until the *Stop* touches the *Sleeve*. After installation of all the implants, manually remove the *Clamps* and the *Surgical Guide*. NOTE: Do not interrupt the rotation of the motor during drilling.

SANITATION

This products should be correctly cleaned after each use. Do not leave it in contact with dampness for prolonged periods. Follow these steps for cleaning: 1st step: Remove drill from contra-

angle (if applicable). 2nd step: Immerse drill completely in enzymatic detergent (10% solution). 3rd step: Wash it in an ultrasonic washer for approximately 10 minutes. 4th step: Rinse it with plenty of distilled water, until all solution residues have been completely removed. The use of nylon brushes is permitted. 5th step: Dry it with a paper towel or compressed air. 6th step: Perform a visual inspection, observing whether there are any faults in the cleaning process. If you still find any residues, reimmerse the drill in detergent - 2nd step - and clean it using a nylon brush, if necessary. Repeat rinsing and drying sequence. 7th step: If applicable, place the drill in an appropriate, properly cleaned and dried case. 8th step: Select packaging according to sterilization process. Preferably use self-sealing package of medical grade paper with laminated film. IMPORTANT: Never store drills without drying them thoroughly. Drying is very important to prevent oxidation. CAUTION: Do not use incrustation removal products, as they tend to darken drills and cause oxidation. Avoid using enzymatic liquid in concentrations above 10%, and drying drills that still contain residues from cleaning solutions, as these procedures favor oxidation. The use of steel brushes is also not recommended.

The drying of the parts is of utmost importance before storage and sterilization, because the accumulation of moisture on the products is harmful and may cause oxidation.

PRESENTATION AND STERILIZATION

These products are reusable and supplied non-sterile, packed individually in a package that provides double protection. These products must be correctly sanitized and sterilized before use. Sterilize the products the day before or on the day of the procedure. ATTENTION: It is not recommended to autoclave these products in their original packaging.

We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard BS EN ISO 17665-1:2006 (Sterilization of health care products. Moist heat Requirements for the development, validation and routine control of a sterilization process for medical devices). Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Suitable parameters: 121°C, 30 min, at 1 atm. Sterilization validity: 7 to 15 days, provided that the product is conditioned in clean, dry and protected from the sunlight environment.

PRECAUTIONS

- Be sure you have all the instruments needed to perform the surgery according to the surgical planning.
- Inappropriate planning can jeopardize the performance of the implant unit leading to failures of the system, as loss or fracturing of the implant.
- Be aware in cases patients that show signs of allergy or hypersensitivity to chemical components of the material. *Clamp*: titanium alloy Ti6Al4V-ELI. *Drill 1.3*: surgical stainless steel.

- The material to be used during the procedure must be sterile.
- Do not use the product if the packaging is broken.
- Do not use the product with the validity expired.
- Before each procedure, make sure the pieces are properly seated.
- Before each procedure check the conditions of the **Neodent** surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- During the surgical procedure, make sure you don't exert lateral movements with the surgical instruments, avoiding the loss of the mounting of the *Sleeve*.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Always use the sequence of **Neodent** products for the **Neodent Guided Surgery** technique. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the **Neodent Implant System** and exempts any product warranty.
- The **Neodent Guided Surgery 1.3** drill should not be resharpened.
- It is the professional's responsibility to use the **Neodent** products according to the instructions for use.

ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

After use, all the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

LIFE CYCLE

This product should be disposed of when it loses its functionality.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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