# ST.330.306

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This device is intended for a specialized procedure, which should be performed by qualified professionals with specific training on the CAD/CAM technique. The CAD/CAM Dental Wings (*Dental Wings IntraOral • DWIO*) is recommended. To avoid unsatisfactory results for the patient, do not use the product without knowledge of the proper techniques. Always follow the correct procedures, under proper conditions and with proper equipment.

## DESCRIPTION

The Intraoral Scanbody is made of PEEK (high-performance polymer, specific for medical use). On External Hex, II Plus (HI), Morse Taper, and WS prosthetic interfaces, they come in a set with a specific screw made of titanium alloy, according to standard ASTM F136. The Intraoral Scanbody has a marking to identify the prosthetic interface, is available in the rotational or anti-rotational forms in relation to the coupling of the prosthetic component or implant, and is available according to the table below:

Prosthetic interface	Identification on the component	Form of presentation	Type of Restoration (CAD design)
HE 3.3	1	Anti-rotational	Multiple or single-tooth
HE 4.1 4.3	2	Anti-rotational	Multiple or single-tooth
HE 5.0	3	Anti-rotational	Multiple or single-tooth
Universal Abutment height 4	4	Anti-rotational	Single-unit
II Plus (HI)	5	Anti-rotational	Multiple or single-tooth
Universal Abutment height 6	6	Anti-rotational	Single-unit
CM Exact	8	Anti-rotational	Single-unit
Micro Pillar	9	Rotational	Multiple
Mini Conical Abutment 4.1 4.3	10	Rotational	Multiple
Mini Conical Abutment 5.0	11	Rotational	Multiple
Conical Abutment 4.1 4.3	12	Anti-rotational	Multiple or single-tooth
Conical Abutment 5.0	13	Anti-rotational	Multiple or single-tooth
CM   GT Abutment and WS Abutment	14	Anti-rotational	Multiple or single-tooth
WS	15	Rotational	Single-unit
Facility	16	Anti-rotational	Single-unit

## INDICATIONS FOR USE

This product is a prosthetic component used to transfer the position and orientation of the implant or prosthetic component, according to each case planned, respecting dimensions and indications of tolerance of the implants and prosthetic components, aiming for aesthetics and functionality.

## APPLICATIONS

This product is used over an implant and/or prosthetic components to transfer their positions during the intraoral scan using the CAD/CAM technique.

## CONTRAINDICATIONS

This product is not indicated for extraoral scanning.

This product is not indicated for customization.

This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: PEEK.

## HANDLING

Identify the Intraoral Scanbody compatible with the interface of the implant/prosthetic component installed. Screw or fit the Intraoral Scanbody onto the implant/prosthetic component.

In transfer of External Hex, II Plus (HI), Morse Taper, and WS prosthetic interfaces, the fit onto the implant is done using a 1.2 mm screwdriver (maximum torque of 10 N.cm).

With the Intraoral Scanbody over an Abutment, the fit is done with a “click”.

With the Intraoral Scanbody over a Facility implant, the fit is done with a slight manual pressure.

Check the fit between the pieces, and make sure there is no rotational or vertical looseness. Looseness indicates that the Intraoral Scanbody is not properly coupled and should be repositioned. Make sure to align the Intraoral Scanbody to the axis of insertion of the implant/prosthetic component, preventing lock-ups and damage to the system.

Follow the instructions provided with the scanner used for the scanning procedure; if necessary, spray the region for scanning, avoiding the Intraoral Scanbody, with the spray or liquid for scanning, according to the recommendations of the scanner manufacturer.

After scanning, check the alignment between the Intraoral Scanbody and the CAD/CAM library.

## PRESENTATION AND STERILIZATION

This product is indicated for single use and provided in non-sterile form, and must therefore be sterilized before use. It is packaged individually.

It may come with a titanium screw, according to the table below.

Interface type	Components in the package
External hex (HE)	Scanbody and screw
Universal Abutment	Scanbody
II Plus / HI	Scanbody and screw
CM Exact	Scanbody with attached screw
Micro Pillar	Scanbody and screw
Mini Conical Pillar	Scanbody and screw
Conical Pillar	Scanbody and screw
CM   GT Abutment and WS Abutment	Scanbody and screw
WS	Scanbody and screw
Facility	Scanbody

For Intraoral Scanbodies with HE, II Plus (HI), Micro Abutment, Conical Mini Abutment, Conical Abutment, CM | GM Abutment, WS Abutment, and WS interfaces, the screw must be removed for sterilization.

It is recommended, preferably, to follow the method of sterilization by steam autoclave (moist heat), gravitational or vacuum autoclave, with exposure time of 3 minutes at 132 ° C (270 ° F). The product must be unpacked and placed on an appropriate tray. Use it immediately after sterilization, do not store.

## PRECAUTIONS

- Use of the Intraoral Scanbody on Facility, WS, and Morse Taper for scanning and fabrication of multiple prostheses is not recommended.
- Use of the Intraoral Scanbody on II Plus (HI) implants with angles between them for scanning and fabrication of multiple prostheses is not recommended.
- The Intraoral Scanbody, when used with II Plus (HI) angled implants, and Facility, WS, and Morse Taper implants has a limitation for fabrication of multiple prostheses on Titanium Bases or directly on the head of the implant.
- For multiple restorations on Facility, Morse Taper, and WS implants, an intermediary must be used over the implant and its respective scanbody during scanning.
- Make sure not to use the spray/liquid for scanning outside the oral cavity.
- Before using the product, make sure that it is clean and that it has the same prosthetic interface as the implant/prosthetic component to be scanned.
- When using the Intraoral Scanbody on a Universal Abutment, make sure that the heights match.
- The use of the Intraoral Scanbody with a damaged platform may lead to unsatisfactory scanning results.
- Make sure that the stability of the implant is sufficient to support the procedures of insertion and removal of the Intraoral Scanbody.

- Surgical planning and/or inadequate prosthesis can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of components and/or prosthetic screws.
- This product is of single use and cannot be re-sterilized.
- Reuse of this product may cause:
  - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
  - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if the packaging is broken.
- Do not use the product after expiry date.
- For immediate load application, check the torque indication of the implant installed.
- See the torque to be given on the prosthetic component to be used. The excess or lack of torque can lead to undesirable result.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Before each procedure check the conditions

of the **Neodent** surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.

- Always use the **Neodent** product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the **Neodent** Implant System and exempts any product warranty.
- It is the professional's responsibility to use the **Neodent** products according to the instructions for use.

#### ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

#### MAGNETIC RESONANCE IMAGING (MRI)

##### • SAFETY INFORMATION

The **Neodent** Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the **Neodent** Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

#### POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the procedure and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

#### STORAGE CONDITIONS

This product should be stored, in its original packaging, in a clean and dry location, in a maximum temperature of 40 °C and protected from direct sunlight.

#### DISPOSAL OF MATERIAL

After use, all the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

#### DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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