ST.330.308

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions.

Neodent titanium Prosthetic Components are used as intermediaries between the implant placed intraorally and the prosthesis, planned according to each case, respecting the prosthetic interfaces, aiming for aesthetics and function.

DESCRIPTION

It is an intermediary prosthetic component to be installed between the Implant and the prosthesis (crown), made of titanium alloy according to standard ASTM F136.

It has a conical shape and smooth surface. Its cementable region can be customized according to the mesio-distal and interocclusal space available. The Preparable Abutment is available for External Hex (4.1 / 5.0), Internal Hex, and Morse Taper Implants. It is available as shown in the table below:

Description / Prosthetic Interface	Diameter	Trans-mucosal height (mm)	Cementable area	Form of presentation in relation to the implant
Preparable Abutment (Ext. Hex)	5.0 for platform 4.1	2.3	8.7	Anti-rotational
	6.0 for platform 4.1	2.5	8.5	
	6.0 mm for platform 5.0	2.5	8.5	
Preparable Abutment II	5.0 for platform 4.3	2.3	8.7	
CM Preparable Abutment	3.3 and 4.5	0.2	4 and 6	Not applicable

INDICATIONS FOR USE

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading. Multiple teeth applications can be rigidly splinted.

APPLICATIONS

It is indicated for single- or multiple-tooth prostheses cemented onto implants, according to the available interocclusal space, transmucosal height, and three-dimensional position of the implant. It can be used in immediate or conventional rehabilitation procedures on the maxilla or mandible. The CM Preparable Abutment is indicated for cases where the Implant is exposed to the oral cavity.

CONTRAINDICATIONS

For customized components on CM en Ext. Hex prosthetic interfaces with angling greater than 17° in combination with implants of diameter equal to or less than 3.75 mm, they are not recommended for single-tooth final applications in the molar region.

This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: Titanium alloy Ti6Al4V.

Contraindicated for insufficient interocclusal space and unsatisfactory three-dimensional position of the implant.

HANDLING

For use of the Preparable Abutment in two-stage procedures, prior preparation of the soft tissues can be done with use of a healing cap.

For use in the mouth (applicable for the II and Ext. Hex interface: select the Preparable Abutment according to the prosthetic planning and install it on the implant, using the torque and connection recommended (table below). Make sure of the correct fit between the Preparable Abutment and the implant. Protect the access for the screw with teflon and resin compound). Perform the customization of the Preparable Abutment according to the need for interocclusal and interproximal space; it is clinically recommended to maintain a minimum height of 4 mm above the transmucous. The customization of the Preparable Abutment must be done with cooling and appropriate techniques, following the guidance of the information below. It is indicated to make a provisional prosthesis for gingival conditioning. Take the impression of the Preparable Abutment according to the appropriate technique. Tests should be carried out on passivity and the fit of the structure of the prosthesis. Cement the definitive prosthesis onto the Preparable Abutment. For prior use in the laboratory (applicable to the CM, II and External Hex interfaces): take the impression of the implant according to the interface of the implant and the prosthetic planning, using the appropriate technique. Prepare the plaster model and select the Preparable Abutment according to the prosthetic planning and fit it onto the plaster model, not exceeding the torque recommended (table below). Make sure of the correct fit between the Preparable Abutment and the analog. Perform the customization of the Preparable Abutment according to the need for interocclusal and interproximal space; it is clinically recommended to maintain a minimum height of 4 mm above the transmucous. The customization of the Abutment must be done in accordance with the appropriate techniques, following the guidelines according to the information below. It is indicated to make a provisional prosthesis for gingival conditioning. Tests should be carried out on passivity and the fit of the structure of the prosthesis. During the installation, protect the access of the screw (with teflon and resin compound). Cement the definitive prosthesis on the Preparable Abutment. For selection of the CM analog, ask the sales consultant.

PARAMETERS FOR CUSTOMIZATION:
 For the CM and External Hex interface, the minimum thickness of the abutment in the region of the seating base around the channel for the screw must be 0.43 mm; for the II interface it must be 0.5 mm. The minimum diameter of the

emergence profile for the CM interface must be 3.3 mm, for the Ext. Hex interface it must be 4 mm, and for the II interface, it must be 4.3 mm. For the CM interface, the minimum height of the emergence profile is 1.5 mm, for the Ext. Hex interface, the minimum height of the emergency profile must be 2 mm, and for the II interface it must be 3 mm. The maximum limit for angling is 30°.

For cementation in the mouth, follow the instructions of the manufacturer of the cement. In this case, it is worth remembering that it is necessary to take note of the technique to prevent excess cement and chemical burns on the tissues adjacent to the Implant. Protect the access of the screw (with teflon and resin compound) during the cementing process.

Prosthetic interface	Torque connection (mm)	Recommended torque (N.cm)
Morse Taper	0.9	15
II	1.2	20
HE		32

TRACEABILITY LABEL

This product is accompanied by three labels that allow its traceability and should be attached to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

SANITIZATION

This product must be correctly sanitized if it is not installed directly in mouth. For sanitization, proceed as follows: 1st step: Immerse completely in enzymatic detergent (diluted according to the manufacturer). 2nd step: Wash in an ultrasound washer for approximately 10 to 15 minutes. 3rd step: Rinse with plenty of distilled water, until the solution residues are completely removed. It is recommended to use nylon brushes. 4th step: Dry with a clean, dry cloth or with compressed air. 5th step: Conduct a visual inspection, observing whether there are failures in the cleaning process. If there are still residues, the piece must be immersed in detergent again - 1st step - and, if necessary, the cleaning should be done with the aid of a nylon brush. Repeat the sequence of rinsing and drying. ATTENTION: Do not use descaling agents, and avoid drying pieces that still contain residues of cleaning solution, because these procedures promote oxidation. The use of steel brushes is not recommended.

PRESENTATION AND STERILIZATION

This product is indicated for single use and provided sterilized by the ethylene oxide method, packed individually into a package.

PRECAUTIONS

- If the Prepable Abutment is used for a single-tooth application, it is clinically recommended to maintain a minimum height of 4 mm, in addition to the transmucosal height, in the creation of the prosthesis.
- Excessive wear may compromise the mechanical properties and/or render the product unusable.
- The CM Prepable Abutment is indicated only for customization in the laboratory.
- Surgical planning and/or inadequate prosthesis can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of components and/or prosthetic screws.
- The selection of the material of the structure of the prosthesis must take local aspects of the patient into account.
- This product is of single use.
- Reuse of this product may cause:
 - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
 - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if the packaging is broken.
- The sterilization is guaranteed only if the packaging is not damaged.
- Do not use the product after expiry date.
- For immediate load application, check the torque indication of the implant installed.
- See the torque to be given on the prosthetic

component to be used. The excess or lack of torque can lead to undesirable result.

- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- During installation, make sure to align it with the insertion axis of the implant. Make sure it is perfectly seated on the Implant. To do so, it is recommended to use periapical X-rays with the parallel technique.
- Make sure to use a prosthetic screw that is compatible with the prosthetic interface and with the prosthetic component.
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/prosthesis assembly.
- Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the Neodent product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the professional's responsibility to use the Neodent products according to the instructions for use.

ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the

MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the procedure and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored, in its original packaging, in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

After use, all the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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