ST.330.309

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English  Angled Abutment

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Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Tamanho do produto / Product size / Tamaño del producto / Dimensioni del prodotto / Produktmaße / Taille du produit
	Código do produto / Product code / Código del producto / Codice del prodotto / Produktcode / Code du produit
	Número do lote / Batch Number / Número de Partida / Numero di lotto / Chargennummer / Numéro de lot
	Data de fabricação / Date of manufacture / Fecha de fabricación / Data di fabbricazione / Herstellungsdatum / Date de fabrication
	Prazo de validade / Shelf life / Plazo de validad / Durata di conservazione / Mindesthaltbarkeit / Date de péremption
	Quantidade / Quantity / Cantidad / Quantità / Qualität / Quantité
	Material de fabricação do produto / Material used / Material de fabricación del producto / Materiale di fabbricazione del prodotto / Herstellungsmaterial des Produkts / Matériau de base du produit
	Requisito Brasileiro / Brazilian Requirement / Requisito Brasileño / Requisito Brasiliano / Brasilianische Anforderungen / Exigences brésiliennes
	Esterilizado por óxido de etileno / Sterile - ethylene oxide / Esterilizado con óxido de etileno / Sterile - ossido di etilene / Sterilisiert mit Ethylenoxid / Stérilisé à l'oxyde d'éthylène

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Não reutilizar / Do not reuse / No lo reutilice / Non riutilizzare / Nicht wiederverwendbar / Ne pas réutiliser / Nicht wiederverwendbar / Ne pas réutiliser
	Consulte as instruções de utilização / Refer to instructions for use / Consulte las instrucciones de utilización / Consultare le istruzioni per l'uso / Siehe Bedienungsanleitung / Consulter les instructions d'utilisation
	Limite superior de temperatura / Upper limit of temperature / Limite superior de temperatura / Limite di temperatura massima / Maximaltemperatur / Limite de température supérieure
	Conservar seco / Keep dry / Consérvelo seco / Conservare in luogo asciutto / Trocken lagern / Garder au sec
	Manter afastado da luz solar / Keep protected from sunlight / Manténgalo lejos de la luz solar / Tenere al riparo dalla luce solare / Vor Sonneneinstrahlung schützen / Conservar à l'abri de la lumière solaire
	Não utilizar se a embalagem estiver danificada / Do not use if package is damaged / No lo utilice si el envoltorio está dañado / Non utilizzare se la confezione originale è danneggiata / Bei beschädigter Verpackung Produkt nicht verwenden / Ne pas utiliser si l'emballage a été abîmé

Simbologia Symbology Simbología Simbologia Symbolik Symboles	Descrição / Description / Descripción / Descrizione / Beschreibung / Description
	Marcação CE para comercialização na Comunidade Europeia / CE Mark for European Community market / Marca CE para comercialización en la Comunidad Europea / Marchio CE per la commercializzazione in Europa / CE Kennzeichen zum Vertrieb in der Europäischen Gemeinschaft / Marquage CE pour une commercialisation en Europe
	Notificação exigida pelo FDA para comercialização nos Estados Unidos / Notification required by FDA for United States market / Notificación exigida por el FDA para comercialización en los Estados Unidos / Notifica necessaria per la commercializzazione negli Stati Uniti / Vom FDA geforderte Kennzeichnung für den Vertrieb in den USA / Notification exigée par FDA pour une commercialisation aux États-Unis
	Representante na Comunidade Europeia / Representative in the European Community / Representante en la Comunidad Europea / Rappresentante nella Comunità europea / Vertreter in der Europäischen Gemeinschaft / Représentant dans la Communauté Européenne
	Fabricante / Manufacturer / Fabricante / Produttore / Hersteller / Fabricant

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions.

Neodent titanium Prosthetic Components are used as intermediaries between the implant placed intraorally and the prosthesis, planned according to each case, respecting the prosthetic interfaces, aiming for aesthetics and function.

DESCRIPTION

It is an intermediary prosthetic component to be installed between the Implant and the prosthesis (crown), made of titanium alloy according to standard ASTM F136. They are provided in a form that is anti-rotational in relation to the implant. It has a tapered shape and smooth surface, and is available in angulations of 17° and 30°. It comes in two parts (Body and Screw). The Angled Abutment for External Hex allows for 12 positions for coupling with the implant and the model for Internal Hex allows for 6 positions. It is available as shown in the table below:

Description / Prosthetic Interface	Prosthetic platform	Trans-mucosal height (mm)	Cementable Area (mm)
Angled Abutment 17° (HE)	3.3	2.1, 3, and 4	7.7
	4.1	2.2, 3, and 4	
	4.3	2.2, 3, and 4	
	5.0	2.5, 3, and 4	
Angled Abutment 30° (HE)	4.1	3, 4, and 5	
Angled Abutment II Plus 17°	4.3 and 5.0	2, 3, and 4	
Angled Abutment II Plus 30°	4.3 and 5.0	3, 4, and 5	

INDICATIONS FOR USE

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading. Multiple teeth applications can be rigidly splinted.

APPLICATIONS

It is indicated for single- or multiple-tooth prostheses cemented onto implants, according to the available interocclusal space, transmucosal height, and three-dimensional position of the implant. It can be used in immediate or conventional rehabilitation procedures on the maxilla or mandible.

CONTRAINDICATIONS

30° Angled Abutments in combination with External Hex Implants with diameter of 3.75 mm are contraindicated for final rehabilitations in the molar region.

This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: Titanium alloy Ti6Al4V.

Contraindicated for insufficient interocclusal space and unsatisfactory three-dimensional position of the implant.

HANDLING

For use of the Angled Abutment in two-stage procedures, prior preparation of the soft tissues can be done with use of a healing cap.

For use in the mouth: select the Angled Abutment according to the prosthetic planning and install it on the implant, using the torque and connection recommended (table below). Make sure of the correct fit between the Angled Abutment and the implant. Protect the access for the screw with teflon and resin compound). It is indicated to make a provisional prosthesis for gingival conditioning. Take the impression of the Angled Abutment using the appropriate technique. Tests should be carried out on passivity and the fit of the structure of the prosthesis. Cement the definitive prosthesis on the Angled Abutment.

For prior use in the laboratory: take the impression of the implant according to the interface of the implant and the prosthetic planning, using the appropriate technique. Prepare the plaster model and select the Angled Abutment according to the prosthetic planning and fit it onto the plaster model, not exceeding the torque recommended (table below). Make sure of the correct fit between the Angled Abutment and the analog. It is indicated to make a provisional prosthesis for gingival conditioning. Tests should be carried out on passivity and the fit of the structure of the prosthesis. During the installation, protect the access of the screw (with teflon and resin compound). Cement the definitive prosthesis on the Angled Abutment. Cementation is done in the mouth following the instructions from the cement manufacturer. In this case, it is worth remembering that it is necessary to take note of the technique to prevent excess cement and chemical burns on the tissues adjacent to the Implant. Protect the access of the screw (with teflon and resin compound) during the cementing process.

Prosthetic component	Torque connection (mm)	Recommended torque (N.cm)
Angled Abutment (External Hex)	1.2	32
Angled Abutment II Plus		20

TRACEABILITY LABEL

This product is accompanied by three labels that allow its traceability and should be attached to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

SANITATION

This product must be correctly sanitized if it is not installed directly in mouth.

For sanitization, proceed as follows: 1st step: Immerse completely in enzymatic detergent (diluted according to the manufacturer). 2nd step: Wash in an ultrasound washer for approximately 10 to 15 minutes. 3rd step: Rinse with plenty of distilled water, until the solution residues are completely removed. It is recommended to use nylon brushes. 4th step: Dry with a clean, dry cloth or with compressed air. 5th step: Conduct a visual inspection, observing whether there are failures in the cleaning process. If there are still residues, the piece must be immersed in detergent again - 1st step - and, if necessary, the cleaning should be done with the aid of a nylon brush. Repeat the sequence of rinsing and drying. **ATTENTION:** Do not use descaling agents, and avoid drying pieces that still contain residues of cleaning solution, because these procedures promote oxidation. The use of steel brushes in not recommended.

PRESENTATION AND STERILIZATION

This product is indicated for single use and provided sterilized by the ethylene oxide method, packed individually into a package.

PRECAUTIONS

- Care should be taken in relation to passivity in cases of impression of multiple angled implants.
- It is clinically advisable to maintain the minimum height of 4 mm above the transmucosal height.
- The Angled Abutment should be used without additional angling during any customization.
- Excessive wear may compromised the mechanical properties and/or render the product unusable.
- Surgical planning and/or inadequate prosthesis can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of components and/or prosthetic screws.
- The selection of the material of the structure of the prosthesis must take local aspects of the patient into account.
- This product is of single use.
- Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if the packaging is broken.
- The sterilization is guaranteed only if the packaging is not damaged.
- Do not use the product after expiry date.
- For immediate load application, check the torque indication of the implant installed.
- See the torque to be given on the prosthetic component to be used. The excess or lack of torque can lead to undesirable result.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.

- During installation, make sure to align it with the insertion axis of the implant. Make sure it is perfectly seated on the Implant. To do so, it is recommended to use periapical X-rays with the parallel technique.
- Make sure to use a prosthetic screw that is compatible with the prosthetic interface and with the prosthetic component.
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/prosthesis assembly.
- Before each procedure check the conditions of the **Neodent** surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the **Neodent** product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the **Neodent** Implant System and exempts any product warranty.
- It is the professional's responsibility to use the **Neodent** products according to the instructions for use.

ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

MAGNETIC RESONANCE IMAGING (MRI) - SAFETY INFORMATION

The **Neodent** Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the **Neodent** Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the procedure and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored, in its original packaging, in a clean and dry location, in a maximum temperature of 40° C and protected from direct sunlight.

DISPOSAL OF MATERIAL

After use, all the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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