# ST.330.310

**Português**  Macho de Rosca Neodent Guided Surgery Facility

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants, with knowledge of the **Neodent Guided Surgery** technique. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

## DESCRIPTION

The Facility **Neodent Guided Surgery Bone Tap** is made of surgical stainless steel. At one end, it has a coupling which allows it to fit with the Facility Bone Tap Ratchet Connection. The active portion of the Facility **Neodent Guided Surgery Bone Tap** has cutting edges for formation of threads and laser markings for the length of the Facility Implants, using the **Neodent Guided Surgery** technique.

## INDICATIONS FOR USE

The **Neodent Guided Surgery** technique is indicated for surgical procedures on maxilla or mandible bones with a sequence of guided instruments, allowing for flapless surgery for installation of implants in edentulous areas.

## APPLICATIONS

The Facility **Neodent Guided Surgery Bone Tap** is indicated for formation of threads in the surgical socket, using the **Neodent Guided Surgery** technique, before installation of Facility Implants in type I and II bone bed, according to the Lekholm and Zarb bone classification (1985). Lekholm and Zarb bone classification (1985):

- Bone Type I = Homogeneous compact bone.
- Bone Type II = Compact bone with trabecular bone core.
- Bone Type III = Thin layer of cortical bone around dense trabecular bone with good resistance.
- Bone Type IV = Thin layer of cortical bone around dense trabecular bone with low density.

## CONTRAINDICATIONS

The Facility **Neodent Guided Surgery Bone Tap** is contraindicated for bone types III and IV.


## HANDLING

To use the Facility **Neodent Guided Surgery Bone Tap**, first use the drills according to the planning and the implant selected.

Fit the Facility **Neodent Guided Surgery Bone Tap** on the Contra-Angle handpiece. Set the surgical motor to a rotation speed of no more than 30 rpm and maximum torque of 35 N.cm. Start the motor and insert the Facility **Neodent Guided Surgery Bone Tap** into the surgical socket, maintaining the drilling axis, until stability is attained. Proceed with the Torque Meter Ratchet Wrench, using the Hex Connection for Ratchet. Make clockwise insertion motions with light pressure on the Ratchet Wrench, with maximum torque of 45 N.cm. Insert the Facility **Neodent Guided Surgery Bone Tap** until you reach the marking for the implant chosen (see image below). It may be necessary to alternate the direction between insertion and removal a few times to form the thread.

To remove the Facility **Neodent Guided Surgery Bone Tap** from the surgical socket, invert the direction of the Torque Meter Ratchet Wrench to counterclockwise and proceed with full removal of the Facility **Neodent Guided Surgery Bone Tap**,

otherwise it can damage the shape of the thread.

Facility Neodent Guided Surgery Bone Tap

ATTENTION: The figure is merely illustrative. It does not represent the real dimension and characteristics of the product. The **Neodent Facility Neodent Guided Surgery** has an additional length, compensating for the thickness of the soft tissue, the height of the Sleeve.

Nominal size (mm)	Real size (mm)
10	16
12	18
14	20

## SANITIZATION

This product must be correctly sanitized after each use.

Do not leave it in contact with dampness for prolonged periods. Follow these steps for cleaning: 1st step: Remove drill from contra-angle (if applicable). 2nd step: Immerse drill completely in enzymatic detergent (10% solution). 3rd step: Wash it in an ultrasonic washer for approximately 10 minutes. 4th step: Rinse it with plenty of distilled water, until all solution residues have been completely removed. The use of nylon brushes is permitted. 5th step: Dry it with a paper towel or compressed air. 6th step: Perform a visual inspection, observing whether there are any faults in the cleaning process. If you still find any residues, reimmerse the drill in detergent - 2<sup>nd</sup> step - and clean it using a nylon brush, if necessary. Repeat rinsing and drying sequence. 7th step: If applicable, place the drill in an appropriate, properly cleaned and dried case. 8th step: Select packaging according to sterilization process. Preferably use self-sealing package of medical grade paper with laminated film. IMPORTANT: Never store drills without drying them thoroughly. Drying is very important to prevent oxidation. CAUTION: Do not use incrustation removal products, as they tend to darken drills and cause oxidation. Avoid using enzymatic liquid in concentrations above 10%, and drying drills that still contain residues from cleaning solutions, as these procedures favor oxidation. The use of steel brushes is also not recommended.

The drying of the parts is of utmost importance before storage and sterilization, because the accumulation of moisture on the products is harmful and may cause oxidation.

## PRESENTATION AND STERILIZATION

This product is reusable and supplied non-sterile, packed individually in a package that provides double protection. This product must be correctly sanitized and sterilized before use. Sterilize the products the day before or on the day of the procedure. ATTENTION: It is not recommended to autoclave this product in its original packaging. Sterilize the product on the day before or on the day of the procedure. We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard BS EN ISO 17665-1:2006 (Sterilization of health care products. Moist heat Requirements for the development, validation and routine control of a sterilization process for medical devices). Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Suitable parameters: 121°C, 30 min, at 1 atm. Sterilization validity: 7 to 15 days, provided that the product is conditioned in clean, dry and protected from the sunlight environment.

## PRECAUTIONS

- Inappropriate planning can jeopardize the performance of the implant unit leading to failures of the system, as loss or fracturing of the implant.
- Some of the possible cause of potential failure in osseointegration and loss of the prosthesis during the treatment are: inadequate osteotomy, infections, deficient oral hygiene, occlusal trauma, diseases or systemic problems, low quantity or quality of remaining bone, lack or failure in irrigation, use of non-specific instruments and/or instruments without cutting power, and lack of specific training.
- Be aware in cases patients that show signs of allergy or hypersensitivity to the chemical components of the material: surgical stainless steel.
- Do not use the product if the packaging is broken.
- Do not use the product with the validity expired.
- Before each procedure, make sure the pieces are properly seated.
- Be sure you have all the instruments needed to perform the procedure according to the surgical planning.
- Before each procedure check the conditions of the **Neodent** surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- During the surgical procedure, make sure you do not exert lateral movements with the surgical instruments, avoiding the loss of the mounting of the Sleeve.
- Perform the preparation of the surgical alveolus according to the instrument sequence of the Alvim implant chosen and surgical planning.
- This product must be handled securely, avoiding contact with nearby soft tissues beyond the area of interest.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Always use the product **Neodent** sequence for the **Neodent Guided Surgery** technique. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the **Neodent Implant System** and exempts any product warranty.

- It is the professional's responsibility to use the Neodent products according to the instructions for use.

#### **ADVERSE EFFECTS**

No adverse effects will occur, provided that the indications are complied with.

#### **POST-OPERATIVE PRECAUTIONS AND MAINTENANCE**

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

#### **STORAGE CONDITIONS**

This product should be stored in a clean and dry

location, in a maximum temperature of 40°C and protected from direct sunlight.

#### **DISPOSAL OF MATERIAL**

After use, every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

#### **DATE OF EXPIRATION**

Written on the label.

#### **LIFE CYCLE**

The Facility Neodent Guided Surgery Bone Tap has a useful lifespan of up to 30 instrumentations for bone qualities I and II, provided that its conditions for use are respected, including adequate irrigation and rotation. Regardless of the lifespan mentioned, the professional should always check the sharpness condition of the instrumental. Products with improper conditions for cutting should be discarded. Note: Classification of the bone quality according to Lekholm and Zarb (1985).

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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
Macho de Rosca **Neodent Guided Surgery Facility**

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