# ST.330.314

Português 

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
1	103.395	Guided Surgery Drill 1.3
2	125.100	Guided Surgery Guide Clamp
3	103.381	Alvim Guided Surgery Drill 2.0
4	103.382	Alvim Guided Surgery Drill 3.5
5	103.383	Alvim Guided Surgery Drill 4.3
6	103.396	Guided Surgery Pilot Drill 3.5
7	103.388	Guided Surgery Pilot Drill 4.3
8	103.390	Narrow Guided Surgery Pilot Drill

9	111.039	Alvim Guided Surgery Bone Tap 3.5
10	111.040	Alvim Guided Surgery Bone Tap 4.3
11	111.041	Alvim Guided Surgery Bone Tap 3.5
12	125.101	Guided Surgery Guide Stabilizer
13	125.102	Narrow Guided Surgery Guide Stabilizer
14	105.127	Guided Surgery CM Connection - Contra-Angle
15	105.126	Narrow Guided Surgery CM Connection - Contra-Angle
16	105.001	Smart/Ws Implant Driver - Torque Wrench (Short)

17	105.002	Smart/Ws Implant Driver - Contra-Angle
18	104.012	Manual Screwdriver (Medium)
19	125.095	Alvim Guided Surgery Drill Guide 4.3
20	125.097	Alvim Guided Surgery Drill Guide 2.0/3.5
21	125.098	Alvim Narrow Guided Surgery Bone Tap Guide 3.5/4.3
22	125.097	Alvim Narrow Guided Surgery Drill Guide 2.0/3.5

This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants, with knowledge of the Neodent Guided Surgery technique. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

**DESCRIPTION**

Case for storage of instruments for conducting surgical procedures using the Neodent Guided Surgery technique. It is made of autoclavable polymer. It has silicone holders for accommodation and secure attachment of each instrument. It has markings which guide the sequence of use of the instruments and can be equipped according to the procedure.

**INDICATIONS FOR USE**

The Neodent Guided Surgery technique is indicated for surgical procedures on maxilla or mandible bones with a sequence of guided instruments, allowing for flapless surgery for installation of implants in edentulous areas.

**APPLICATIONS**

The Conical Implants Neodent Guided Surgery Surgical Kit Case is indicated for use during installation of CM Alvim and CM Drive implants of diameters Ø 3.5 and Ø 4.3 mm, using the Neodent Guided Surgery technique. It offers secure accommodation of the instruments, guidance on the sequence of use of the surgical instruments, as well as support for their sterilization.

**CONTRAINDICATIONS**

There are no contraindications, provided that it is used correctly for the purposes indicated.

**HANDLING**

- Remove the lid for access to the instruments by pressing gently on the latch system between the base and the lid. To close the Case, place the lid on the base and press down.
- During the surgical procedure, the case should be placed on a sterile surface and handled with a sterile glove.
- The colors of the reserved spaces correspond to the indication of the instruments to be used according to the diameter of the Implant to be installed (see table below):

Color	Diameter of the Implant (in mm)
Purple	3.5
Yellow	4.3

**SETTING UP THE CASE**

To set up this case, each reserved space is related to a number from the instrument table; see the image shown. The instruments that go in the case are sold separately and should be used according to their own Instructions for Use.

**SANITATION**

This product should be correctly cleaned after each use. Do not leave it in contact with dampness for prolonged periods. Follow these steps for cleaning: 1st step: Remove drill from contra-angle (if applicable). 2nd step: Immerse drill completely in enzymatic detergent (10% solution). 3rd step: Wash it in an ultrasonic washer for approximately 10 minutes. 4th step: Rinse it with plenty of distilled water, until all solution residues have been completely removed. The use of nylon brushes is permitted. 5th step: Dry it with a paper towel or compressed air. 6th step: Perform a visual inspection, observing whether

there are any faults in the cleaning process. If you still find any residues, reimmerse the drill in detergent - 2nd step - and clean it using a nylon brush, if necessary. Repeat rinsing and drying sequence. 7th step: If applicable, place the drill in an appropriate, properly cleaned and dried case. 8th step: Select packaging according to sterilization process. Preferably use self-sealing package of medical grade paper with laminated film. **IMPORTANT:** Never store drills without drying them thoroughly. Drying is very important to prevent oxidation. **CAUTION:** Do not use incrustation removal products, as they tend to darken drills and cause oxidation. Avoid using enzymatic liquid in concentrations above 10%, and drying drills that still contain residues from cleaning solutions, as these procedures favor oxidation. The use of steel brushes is also not recommended.

The drying of the parts is of utmost importance before storage and sterilization, because the accumulation of moisture on the products is harmful and may cause oxidation.

### PRESENTATION AND STERILIZATION

This product is reusable and supplied non-sterile, packed individually in a package. This product must be correctly sanitized and sterilized before use. Sterilize the products the day before or on the day of the procedure. **ATTENTION:** It is not recommended to autoclave this product in its original packaging.

We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard BS EN ISO 17665-1:2006 (Sterilization of health care products. Moist heat Requirements for the development, validation and routine control of a sterilization process for medical devices). Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Suitable parameters: 121 °C, 30 min, at 1 atm. Sterilization validity: 7 to 15 days, provided that the product is conditioned in clean, dry and protected from the sunlight environment.

### PRECAUTIONS

- Inappropriate planning can jeopardize the performance of the implant unit leading to failures of the system, as loss or fracturing of the

implant.

- The **Neodent Guided Surgery** technique, using the **Conical Implants Neodent Guided Surgery Surgical Kit Case**, allows for installation of **CM Alvim** and **CM Drive Implants**.
- The **Neodent Guided Surgery** technique does not allow for installation of implants with diameter of 5.0 mm or with length of 18 mm.
- For performance of the procedure, the patient must have sufficient interocclusal space in the area for installation of the implant, for handling of the instruments of the **Neodent Guided Surgery** technique.
  - Do not use the product if the packaging is broken.
  - Do not use the product with the validity expired.
  - Before each procedure, make sure the pieces are properly seated.
  - Be sure you have all the instruments needed to perform the procedure according to the surgical planning.
  - This product must be used sterile. Do not sterilize under dry heat (oven).
  - During sterilization, the case must not touch the walls of the autoclave, to avoid melting.
  - Place the amount of water recommended for the autoclave by the manufacturer. Lack of water during the sterilization cycle may compromise the effectiveness of the autoclaving and lead to melting of the case, as well as damage to the instruments.
  - Contact of the case with acrylic resin liquid results in damage to the structure of the polymer of the case.
  - Chemical sterilization is not recommended.
  - For sterilization, the instruments must be sanitized individually, according to their own instructions, and accommodated in the Case, which must be wrapped in surgical grade paper for autoclaving.
  - Before each procedure check the conditions of the **Neodent** surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
  - Ensure that the parts are not swallowed or aspirated by the patient.
  - Always use the sequence of **Neodent** products

for the **Neodent Guided Surgery** technique. The use of surgical instruments and Sleeves from other manufacturer does not guarantee the function of the **Neodent Implant System** and is exempt from the warranty of the related products.

- It is the professional's responsibility to use the **Neodent** products according to the instructions for use.

### ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

### POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

### STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

### DISPOSAL OF MATERIAL

After use, all the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

### DATE OF EXPIRATION

Written on the label.


### LIFE CYCLE

This product should be disposed of when it loses its functionality.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

1	103.395	Broca Guided Surgery 1.3
2	125.100	Fijador de Guia Guided Surgery
3	103.381	Broca Guided Surgery Alvim 2.0
4	103.382	Broca Guided Surgery Alvim 3.5
5	103.383	Broca Guided Surgery Alvim 4.3
6	103.396	Broca Piloto Guided Surgery 3.5
7	103.388	Broca Piloto Guided Surgery 4.3
8	103.390	Broca Piloto Guided Surgery Estrecho

9	111.039	Macho de Rosca Guided Surgery Alvim 3.5
10	111.040	Macho de Rosca Guided Surgery Alvim 4.3
11	111.041	Macho de Rosca Guided Surgery Alvim 3.5
12	125.101	Estabilizador de Guia Guided Surgery
13	125.102	Estabilizador de Guia Guided Surgery Estrecho
14	105.127	Conexion CM Guided Surgery Para Contra-Angulo
15	105.126	Conexion CM Guided Surgery Para Contra-Ángulo Estrecho
16	105.001	Conexion Hexagonal Corta Para Carraca

17	105.002	Conexion Hexagonal Para Contra-Angulo
18	104.012	Llave Digital Media
19	125.095	Guia de Broca Guided Surgery Alvim 4.3
20	125.097	Guia de Broca Guided Surgery Alvim 2.0/3.5
21	125.098	Guia de Macho Guided Surgery Estrecho Alvim 3.5/4.3
22	125.097	Guia de Broca Guided Surgery Estrecho Alvim 2.0/3.5

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